



YSP- SCIENTIFIC SOLUTION

A to Z Pricelist




Your One-Stop Shop for **Quality and Value.**

BOOK NOW

Contact Us:

 yspscientific@gmail.com

 +91 73387 21875

+91 96262 97366

Why Choose YSP?

- ✓ Product Information & Consultation
- ✓ Custom Pricing Solution
- ✓ Flexible Shipping & Delivery
- ✓ Technical Support and Installation

24/7 Customer Support



OURS BRANDS






1. ANEMOMETER LUTRON / HTC

An anemometer is a crucial instrument used to measure wind speed and direction. Ideal for meteorology, environmental monitoring, and various industrial applications, it provides accurate and reliable data to help in weather forecasting, aviation, and more.

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
A1	AM4201	LUTRON	Model AM 4201 Portable Anemometer Lutron, Taiwan made with fan type vane, range 0.4 to 30m/s knots, feet and kms	₹12,250
A2	AM4202	LUTRON	Model AM 4202 Without RS-232 AP with temperature ,Taiwan made, range 0.4 to 30m/s	₹17,500
A3	AM4204	LUTRON	Model AM 4204 Hot Wire Anemometer Range: 0.2 to 20 m/s (knots, feet, and km) Temperature: 0°C to 50°C	₹26,100
A4	AM4213	LUTRON	HTC AVM-06 Anemometer made with fan type vane, range 0.4 to 30 m/s (knots, feet, and km)	₹34,800
A5	AM4221	LUTRON	Model AM 4221 Lutron, Anemometer Cup counter Taiwan made	₹16,675
A6	AM4222	LUTRON	Model AM 4222 Lutron, Taiwan handy with vane	₹8,750
A7	HTCAVM-06	HTC	HTC AVM-07 Anemometer range 0.4 to 45m/s (knots, feet and kms)	₹7,613
A8	HTC AVM-07	HTC	HTC AVM-07 Anemometer range 0.4 to 45m/s (knots, feet and kms)	₹12,250
A9	HTC AVM-08	HTC	HTC AVM-08 Hot Wire Anemometer 0.1~25.0 m/s (knots, feet and kms) Temp - 0-50'C	₹17,500

ANEMOMETER COMBO - LUTRON



The Lutron Anemometer Combo is a versatile instrument designed to measure wind speed and other related environmental parameters. HVAC systems, and various field measurements where accurate wind data is crucial.

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
A10	LM8000	LUTRON	Model LM 8000 Taiwan 4 in 1 for Velocity, Temperature, Humidity and Lux. Handy built-in vane	₹12,600
A11	LM8010	LUTRON	Model LM 8010 Anemo / Temp/ Air Flow / RH / LUX 5 in 1	₹14,438
A12	LM8102	LUTRON	Model LM 8102 Anemo / Sound Level / Temp / RH / LUX 5 in 1	₹16965

ANEMOMETER POCKETLINE RANGE - TESTO



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
A13	TESTO4101	TESTO	Testo-410-1 Velocity Meter(Vane Anemometer), Two variants with velocity, humidity & temperature, Timed mean calculation, Variant with humidity also calculates psychometric values	₹16,100
A14	TESTO4102	TESTO	Testo-410-2 Velocity Meter (3 in 1 Anemometer), Two variants with velocity, humidity & temperature, Timed mean calculation, Variant with humidity also calculates psychometric values	₹20,554



ANEMOMETER SMART - TESTO

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
A15	TESTO405I	TESTO405I	Testo-405i Smart Thermal Anemometer (Thermal Anemometer) Measurement of air flow velocity , temperature and volume flow, range 0 - 30 m/s and -20°C to +60°C, Telescopic shaft, extendable up to 400 mm, wireless connection to a smartphone or tablet via Bluetooth free app, All measurement data in graph or tabular format, can be directly sent as PDF or excel files, testament 415	₹25,121
A16	TESTO410I	TESTO	Testo-410i Smart Vane Anemometer Measurement of air flow velocity, volume flow and temperature at air outlets, range of 0.4 to 30 m/s and -20°C to +60°C, Telescopic shaft, extendable upto 400 mm, wireless connection to a smartphone or tablet via Bluetooth free app, All measurement data in graph or tabular format, can be directly sent as PDF or excel files	₹21,351



Enhance your wind speed measurements with the **ANEMOMETER SMART - TESTO**. Get accurate, dependable readings for both professional and industrial use.





AUTOCLAVE - YSP

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
A17	YSPBUCKET20	YSP	YSP-BUCKET-20 Bucket type Aluminum-30mm (12"x15") with Pressure Gauge , Heater 20 lt	₹10,938
A18	YSPCOOKER20	YSP	YSP-COOKER-20 Cooker Type S.Steel-30mm(12"x12") / Bucket Model Vertical Type with Pressure Gauge, Heater 20 lt	₹18,850
A19	YSP-COOKER-25	YSP	YSP-COOKER-25 Cooker Type S.Steel-30mm(12"x14") / Bucket Model Vertical Type with Pressure Gauge, Heater 25 lt	₹19,938

AUTOCLAVE VERTICAL ECONOMY (WING NUT TYPE) - YSP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
A20	YSP203	YSP	FT-203 Size 12" X 20" Verical Economy Autoclave, Double Chamber construction, Wing Nut Type locking system, Inner & Outer made of SS 202 grade, support ring is mad of cast aluminium dully polished,Lid Mad of 202 Grade SS , Pressure adjustable from 5 psi to 15 psi	POR
A21	YSP204	YSP	FT-204 Size 12" X 20" Verical Economy Plus Autoclave,Double Chamber construction, Wing Nut Type locking system, Inner & Outer made of SS 202 grade, support ring is mad of cast aluminium dully polished,Lid Mad of 202 Grade SS , Pressure adjustable from 5 psi to 15 psi	₹34,075

AUTOCLAVE VERTICAL LABORATORY (RADIAL LOCKING) - YSP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
		YSP	Laboratory Autoclaves are made of Double chamber, Stainless steel 304 grade. Supplied with standard accessories such as: Radial locking system, paddle lifting arrangement, pressure control switch, low water cut off, manual timer. Dual safety valve, Pressure guage, steam release valve, water outlet valve.	POR
		YSP	SEMI AUTOMATIC: Incorporated with low water cut off device, microprocessor based controller (Programmable by user) digital display of temp and timer. After achieving 121 degree timer starts automatically which is also preset at 30 Minutes (Adjustable). After cycle end unit gets switched off and steam can be released manually. Standard working temperature 121 degree. However as per requirement it can be upgraded upto 134 degree (On Extra Cost)	POR
A22	YSPS301	YSP	FTS-301 , Semi automatic Size 12X 20(300 X 300 mm) Capacity 40 Ltrs, 3 kw load	₹ 97,875
A23	YSPS302	YSP	FTS-302 , Semi automatic Size 12X 20(300 X 300 mm) Capacity 50 Ltrs, 3 kw load	₹110,128
A24	YSPS303	YSP	FTS-303 , Semi automatic Size 12X 20(300 X 300 mm) Capacity 78 Ltrs, 4.5 kw load	₹147,900
A25	YSPS303	YSP	FTS-303 , Semi automatic Size 12X 20(300 X 300 mm) Capacity 78 Ltrs, 4.5 kw load	₹190,856



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
			FULLY AUTOMATIC: Incorporated with low water cut off device, microprocessor based controller (Programmable by user) digital display of temp and timer. Enable reproducibility of results as full cycle included purging of stale air, Sterilization hold time, cycle ends with automatic steam release. Requires no intervention. Standard working temperature 121 degree, However as per requirement it can be upgraded upto 134 degree (On Extra Cost)	POR
A26	YSPF301	YSP	FTF-301 , Fully automatic Size 12X 20(300 X 300 mm) Capacity 40 Ltrs, 3 kw load	₹ 117,450
A27	YSPF304	YSP	FTF-304 , Fully automatic Size 18X 24(450 X 600 mm) Capacity 90 Ltrs, 6 kw load	₹210,431



Sterilize with confidence using the
**AUTOCLAVE VERTICAL LABORATORY
 RADIAL LOCKING - YSP.** Easy to use and
secure, perfect for any lab



AGATE MORTAR AND PASTEL (STONE)

Agate mortar and pestle sets are highly valued in various scientific and laboratory settings for their exceptional durability and chemical resistance. They are commonly used to grind, mix, and homogenize materials, ensuring precise and contamination-free results.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
A28		GENERAL	Agate Black Stone 2"	₹ 963
A29		GENERAL	Agate Black Stone 2.5"	₹ 963
A30		GENERAL	Agate Black Stone 3"	₹1,575
A31		GENERAL	Agate Black Stone 4"	₹2,625
A32		GENERAL	Agate Black Stone 5"	₹4,638
A33		GENERAL	Agate Black Stone 6"	₹6,563
A34		GENERAL	Agate Gray Stone 2"	₹4,375
A35		GENERAL	Agate Gray Stone 3"	₹8,750
A36		GENERAL	Agate Gray Stone 4"	₹14000
A37		GENERAL	Agate Gray Stone 5"	₹25,375
A38		GENERAL	Agate Gray Stone 6"	₹43,500

**We specialize in the sales of all
laboratory products.**

For inquiries or orders,



**96262 97366
73387 21875**



We're here to help!

ALUMINA CRUCIBLE - LARK

Alumina crucibles, especially those manufactured by LARK, are essential laboratory and industrial tools used for high-temperature applications. Known for their superior thermal and chemical properties, these crucibles are ideal for processes that require durability and purity.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
A39	AC15	LARK	Alumina Crucible 15 ml	₹ 613
A40	AC25	LARK	Alumina Crucible 25 ml	₹788
A41	AC30	LARK	Alumina Crucible 30 ml	₹ 875
A42	AC50	LARK	Alumina Crucible 50 ml	₹ 1,050
A43	AC100	LARK	Alumina Crucible 100 ml	₹ 2,013
A44	AC200	LARK	Alumina Crucible 200 ml	₹ 2,975
A45	AC250	LARK	Alumina Crucible 250 ml	₹ 3,281
A46	AC500	LARK	Alumina Crucible 500 ml	₹ 7000
A47	AC1000	LARK	Alumina Crucible 1000 ml	₹ 15,750



**We're your go-to source for
all laboratory products.**

For inquiries or orders,



**96262 97366
73387 21875**

We're here to help!

ALUMINA LID - LARK

Alumina lids, particularly those produced by LARK, are essential components in various high-temperature applications, particularly within laboratory and industrial settings. Known for their excellent thermal and chemical properties, these lids provide reliable performance under extreme conditions.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
A48	AL15	LARK	Alumina Lid 15 ml	₹306
A49	AL25	LARK	Alumina Lid 25 ml	₹394
A50	AL30	LARK	Alumina Lid 30 ml	₹420
A51	AL50	LARK	Alumina Lid 50 ml	₹525
A52	AL100	LARK	Alumina Lid 100 ml	₹1,050
A53	AL200	LARK	Alumina Lid 200 ml	₹1313
A54	AL250	LARK	Alumina Lid 250 ml	₹1,575
A55	AL500	LARK	Alumina Lid 500 ml	POR
A56	AL1000	LARK	Alumina Lid 1000 ml	POR

BAROMETER ANALOG

Analog barometers are traditional instruments used to measure atmospheric pressure, providing valuable data for weather forecasting and various scientific applications.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B1	BARIGO 183	BARIGO	Baraometer BARIGO Cat no.183 Dial size 85 mm German Make	₹ 16,188

YSP- SCIENTIFIC SOLUTION

BUNSEN BURNER ELECTRIC & BURETTE STAND - YSP

The Electric Bunsen Burner by YSP is a modern, safe, and efficient alternative to the traditional gas-powered Bunsen burner. It is designed for use in laboratories where open flames and gas lines pose safety concerns.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B2	YSPBB	YSP	YSP-BB Electric Bunsen with energy regulator up to 350°C Rating 230V	₹ 4,813

YSP- SCIENTIFIC SOLUTION

ROCKER LAB BURNER

The Rocker Lab Burner is a versatile and modern laboratory burner designed for precise and safe heating applications in scientific research, educational institutions, and industrial laboratories. It combines innovative technology with user-friendly features to enhance safety and efficiency in the lab.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B3	DRAGON200	ROCKER	Dragon 200 Lab Burner. Max Temp 1300 deg Celsius. Adjustable flame and Refillable butane tank capacity 35g. Safety Switch is also provided. 	₹ 9,800
B4	DRAGON100	ROCKER	Dragon 100 Lab Burner. Max Temp 1300 deg Celsius. Adjustable flame and Refillable butane tank capacity 40g. Safety Switch is also provided. 	₹ 9,100

BOD INCUBATOR - YSP

The BOD (Biochemical Oxygen Demand) Incubator by YSP is a specialized laboratory instrument designed to provide precise temperature control for the incubation of biological samples.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B5	YSPBOD4	YSP	YSP- BOD-4 Size 400x450x600 mm 4CFT CAPACITY 105LTRS INNER SS OUTER MS Powder coated with digital temp cont 5 to 100c	₹89,900
B6	YSPBOD4	YSP	YSP- BOD-6 Size 410x500x850 mm 6CFT CAPACITY 175LTRS INNER SS OUTER MS Powder coated with digital temp cont 5 to 100c	₹105,850

BOD MEASUREMENT (Installation and DEMO charges applicable)

BOD Measurement (Biochemical Oxygen Demand) is a crucial laboratory procedure used to assess the amount of oxygen required by aerobic microorganisms to decompose organic matter in a sample.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B7	2444460	LOVIBOND	BOD Measurement System LOVIBOND BD 600 6 measuring sensors code no 2444460	₹585,293
B8	2444465	LOVIBOND	BOD Measurement System LOVIBOND BD 606 12 measuring sensors code no 2444465	₹1,428,420

BOMB CALORIMETER - ADVANCE / RICO

Bomb Calorimeters are sophisticated instruments used for measuring the calorific value of substances. They are crucial in various fields, including chemistry, materials science, and energy research. The Advance and RICO models represent advanced solutions for accurate and efficient calorimetric analysis.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B9	RSBT5	RICO	RSBT-5 Digital Bomb Calorimeter With Water Jacket With Puff Insulation, Safety Device, Sensor Probe, High Pressure Flexible Tubing Set & Certified Bomb. Temp. Resolution 0.01°C.	₹ 147,900
B10	RSBT6	RICO	RSBT-6 Automatic Bomb Calorimeter Microprocessor Based 20 X 4 Lcd Display With Keypad On Front Panel, Automatic Calculation & Memory Storage Via Usb Interface In Pc. Supplied With Water Jacket With Puff Insulation, Safety Device, Sensor Probe, High Pressure Flexible Tubing Set & Certified Bomb. Temp. Resolution 0.001°C.	₹ 273,325
B11	RSBT7	RICO	RICO RSBT 7 Bomb Calorimeter With Temp Accuracy 0.001°C Auto Calculation of calorific value with Printer	₹ 317,188
B12	BCMA	ARICO	Advance Make Bomb Calorimeter Model BCM-A With Printer with Auto Gas Filling System, with Printer complete with 0.01°C accuracy Digital Temp. Indicator Push Fit Type Gas Filling System latest Combustion Bomb Design that does not use a Schrader Valve	₹ 297,975



Shop our high-precision weights E1, E2, F1, F2, & M1

Perfect for calibration and lab accuracy. Ensure reliable measurements every time



96262 97366 - 73387 21875

BOMB CALORIMETER :OPTIONAL ACCESSORIES

Bomb Calorimeters are sophisticated instruments used for measuring the calorific value of substances. They are crucial in various fields, including chemistry, materials science, and energy research.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B13	RSBTS2	RICO	RSBT-S-2 Complete Bomb Certified By Testing Lab	₹50,750
B14	RSBTS8	RICO	RSBT-S-8 Stirrer Motor	₹21,750
B15	RSBTS9	RICO	RSBT-S-9 Pellet Press	₹18,125
B16	RSBTS13	RICO	RSBT-S-13 High Pressure Flexible Tube (Set Of 3 Tubes)	₹15,575
B17	RSBTS14	RICO	RSBT-S-14 Cotton Reel (Approx 35 Meters)	₹263
B18	RSBTS15	RICO	RSBT-S-15 Stainless Steel Crucible 8 Ml. Capacity	₹3,500
B19	RSBTS17	RICO	RSBT-S-17 Schrader Valve (For Three Pieces)	₹1,050
B20	RSBTS18	RICO	RSBT-S-18 Ignition Wire (Nichrome Wire) (Minimum 50 Mtrs)	₹4,813
B21	RSBTS22	RICO	RSBT-S-22 Bomb 'O' ring set of 3	₹1,050
B22	RSBTS23	RICO	RSBT-S-23 Valve body 'O' ring Set Of 2	₹1,050
B23	RSBTS27	RICO	RSBT-S-27 Sensor for digital beckmann	₹13,125
B24	RSBTS32	RICO	RSBT-S-32 Connecting wire	₹1,575



BOMB CALORIMETER AUTOMATIC - LABOHOLIC

Experience precise and automated calorimetric analysis with the LABOHOLIC Bomb Calorimeter, designed for accurate energy measurement in various samples, ensuring reliable and repeatable results.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B25	LH1160	LABOHOLIC	LH-1160 Bomb Calorimeter Automatic Bomb Calorimeter Microprocessor Based (Graphical Display) 128x64 Graphical LCD, Temperature Resolution 0.010 °C, Internal Memory 1000 logs with date and time, USB interface to PC, Inbuilt thermal printer, Time Clock, Realtime clock with battery backup. Measurement & calculation of calorific value & water equivalent, USB interface for sending data to PC, full test report printout by inbuilt thermal printer, Weight of sample water equivalent & sensitivity can be set through keypad on single page, Fusewire open indication error adjustable, No paper or printer door open detection, Audio indication on firing as well as after test, Automatic detection of rature rise & fall, Power supply 230 + 10% AC, Complete with all accessories.	₹176,900
B26	LH160	LABOHOLIC	LH-160 Digital Bomb Calorimeter, manual measurement and calculation, Initial temperature rise and fall, 3 Digit LED Display, Temp. accuracy 0.1°C, Temp. resolution 0.1°C, RTD Temperature Detector, Temp. range RT to 99.9°C NEW ADDITION	₹110,563



+91 96262 97366
+91 73387 21875



yspscientific@gmail.com
sales.ypscientific@gmail.com



B: 25 G1, Ground Floor, #7/10
Annai Avanthika Vengaivasal,
Chennai - 600100

BOMB CALORIMETER :OPTIONAL ACCESSORIES

Bomb Calorimeters are sophisticated instruments used for measuring the calorific value of substances. They are crucial in various fields, including chemistry, materials science, and energy research.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B27	RBMT2	RICO	RBMT-1 Ball mill 1 kg. (wth aluminium jar)	₹ 40,600
B28	RBMT3	RICO	RBMT-2 Ball mill 1 kg. (with s.s. jar)	₹46,183



**We're your go-to source for
all laboratory products.**

For inquiries or orders,



96262 97366

73387 21875

We're here to help!

 **YSP- SCIENTIFIC SOLUTION**
BALANCE AP-AD SERIES - SHIMADZU

The BALANCE AP-AD Series by Shimadzu is a line of high-precision analytical balances designed for accurate and reliable measurements in various laboratory and industrial applications.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B29	AP135WAD	SHIMADZU	Model AP135W-AD Super Range(sm) built-in calibration 0.01mg capacity 135g Pan size Ø 91	POR
B30	AP225WAD	SHIMADZU	Model AP225W-AD Super Range(sm) built-in calibration 0.01mg capacity 220g Pan Size Ø 91	POR
B31	AP125WDA D	SHIMADZU	Model AP125WD-AD Dual Range(sm) built-in calibration 0.1mg/0.01mg Cap 120g/52g Pan Size Ø 91	POR
B32	AP225WDA D	SHIMADZU	Model AP225WD-AD Dual Range(sm) built-in calibration 0.1mg/0.01mg Cap 220g/102g Pan Size Ø 91	POR
B33	AP224WAD	SHIMADZU	Model AP224W-AD Analytical built-in calibration 0.1mg Cap 220g Pan Size Ø 91	POR
B34	AP324WAD	SHIMADZU	Model AP324W-AD Analytical built-in calibration 0.1mg Cap 320g Pan Size Ø 91	POR



Get precise measurements effortlessly with the BALANCE AP-AD SERIES - SHIMADZU. Reliable and user-friendly for all your lab needs.



BALANCE AP SERIES - SHIMADZU

The Shimadzu AP Series balances are advanced precision instruments designed for highly accurate weighing applications in laboratories, industrial settings, and research institutions. Known for their reliability, precision, and user-friendly operation, these balances are suitable for a wide range of weighing tasks.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
ADVANCE PERFORMANCE UNIBLOC BALANCES				
B35	AP225W	SHIMADZU	Model AP225W Super Range(sm) Internal calibration 0.01mg capacity 220g Pan size Ø 91	₹717,750
B36	AP135W	SHIMADZU	Model AP135W Super Range(sm) Internal calibration 0.01mg capacity 135g Pan size Ø 91	₹449,500
B37	AP225WD	SHIMADZU	Model AP225WD Semi micro Internal calibration 0.01mg/0.1mg Cap 102g/220g Pan Size 91Dia	₹39,1500
B38	AP125WD	SHIMADZU	Model AP125WD Semi micro Internal calibration 0.01mg/0.1mg Cap 52g/120g Pan Size 91Dia	₹362,500
B39	AP324W	SHIMADZU	Model AP324W Analytical Internal calibration 0.1 mg Cap 320g Pan Size 91Dia	₹377,000
B40	AP224W	SHIMADZU	Model AP224W Analytical Internal calibration 0.1 mg Cap 220g Pan Size 91Dia	₹326,250
B41	A124W	SHIMADZU	Model AP124W Analytical Internal calibration 0.1 mg Cap 120g Pan Size 91Dia	₹282,750
B42	A324X	SHIMADZU	Model AP324X Analytical Internal calibration 0.1 mg Cap 320g Pan Size 91Dia	₹275,500
B43	A224X	SHIMADZU	Model AP224X Analytical Internal calibration 0.1 mg Cap 220g Pan Size 91Dia	₹261,000
B44	A124X	SHIMADZU	Model AP124X Analytical Internal calibration 0.1 mg Cap 120g Pan Size 91Dia	₹253,750
B47	A124Y	SHIMADZU	Model AP124Y Analytical External calibration 0.1 mg Cap 120g Pan Size 91Dia	₹246,500

BALANCE AUW-D SERIES - SHIMADZU

The Shimadzu AUW-D Series balances are high-precision analytical balances designed for accurate and reliable weighing in various laboratory and industrial applications.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B48	AUW220D	SHIMADZU	Model AUW220D Semi micro Internal calibration 0.01mg/0.1mg Cap 82g/220g Pan Size 80 Dia	₹ 288,550

BALANCE AT-R SERIES - NEW (Responsiveness & stability) - SHIMADZU

The Shimadzu AT-R Series balances are advanced precision instruments designed for high responsiveness and stability in various laboratory and industrial applications.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B49	ATX324R	SHIMADZU	Model ATX 324R Capacity 320gm Resolution 0.1mg PAN SIZE 90 DIA	₹188,500
B50	ATX224R	SHIMADZU	Model ATX 224R Capacity 220gm Resolution 0.1mg PAN SIZE 90 DIA	₹174,000
B51	ATX124R	SHIMADZU	Model ATX 124R Capacity 120gm Resolution 0.1mg PAN SIZE 90 DIA	₹172,550
B52	ATX84R	SHIMADZU	Model ATX 84R Capacity 82gm Resolution 0.1mg PAN SIZE 90 DIA	₹172,100
B53	ATY324R	SHIMADZU	Model ATY324R Capacity 320gm Resolution 0.1mg PAN SIZE 90 DIA	₹166,750
B54	ATY224R	SHIMADZU	Model ATY224R Capacity 220gm Resolution 0.1mg PAN SIZE 90 DIA	₹152,250

B55	ATY 124R	SHIMADZU	Model ATY 124R Capacity 120gm Resolution 0.1mg PAN SIZE 90 DIA	₹150,800
B56	ATY 64R	SHIMADZU	Model ATY 64R Capacity 62gm Resolution 0.1mg PAN SIZE 90 DIA	₹150,075

YSP- SCIENTIFIC SOLUTION

BALANCE UP SERIES UPX/UPY - SHIMADZU

The Shimadzu UP Series balances, including the UPX and UPY models, are designed for high-precision weighing applications in various laboratory and industrial settings.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
X SERIES WITH INTERNAL CALIBRATION WEIGHTS				
B57	UP223X	SHIMADZU	Model UP223X Internal Calibration Capacity 220gm Resolution 1mg/0.001g PAN SIZE 108x105 mm	₹143,550
B58	UP423X	SHIMADZU	Model UP423X Internal Calibration Capacity 420gm Resolution 1mg/0.001g PAN SIZE 108x105 mm	₹155,150
B59	UP623X	SHIMADZU	Model UP623X Internal Calibration Capacity 620gm Resolution 1mg/0.001g PAN SIZE 108x105 mm	₹172,550
B60	UP823X	SHIMADZU	Model UP823X Internal Calibration Capacity 820gm Resolution 1mg/0.001g PAN SIZE 108x105 mm	₹200,100
B61	UP1023X	SHIMADZU	Model UP1023X Internal Calibration Capacity 1020gm Resolution 1mg/0.001g PAN SIZE 108x105 mm	₹224,750
B62	UP422X	SHIMADZU	Model UP422X Internal Calibration Capacity 420gm Resolution 10mg/0.01g PAN SIZE 108x105 mm	₹134,850

B63	UP822X	SHIMADZU	Model UP822X Internal Calibration Capacity 820gm Resolution 10mg/0.01g PAN SIZE 108x105 mm	₹156,600
B64	UP2202X	SHIMADZU	Model UP2202X Internal Calibration Capacity 2200gm Resolution 10mg/0.01g PAN SIZE 170x180 mm	₹149,350
B65	UP4202X	SHIMADZU	Model UP4202X Internal Calibration Capacity 4200gm Resolution 10mg/0.01g PAN SIZE 170x1180 mm	₹160,950
B66	UP6202X	SHIMADZU	Model UP6202X Internal Calibration Capacity 6200gm Resolution 10mg/0.01g PAN SIZE 170x180 mm	₹174,000
B68	UP8201X	SHIMADZU	Model UP8201X Internal Calibration Capacity 8200gm Resolution 100mg/0.1g PAN SIZE 170x180 mm	₹169,650
B67	UP4201X	SHIMADZU	Model UP4201X Internal Calibration Capacity 4200gm Resolution 100mg/0.1g PAN SIZE 170x180 mm	₹158,050

Y Series with External Calibration Weights

B69	UP223Y	SHIMADZU	Model UP223Y External Calibration, Capacity 220gm, Resolution 1mg/0.001g, PAN SIZE 108x105 mm	₹12,950
B70	UP423Y	SHIMADZU	Model UP423Y External Calibration, Capacity 420gm, Resolution 1mg/0.001g, PAN SIZE 108x105 mm	₹140,650
B71	UP623Y	SHIMADZU	Model UP623Y External Calibration, Capacity 620gm, Resolution 1mg/0.001g, PAN SIZE 108x105 mm	₹158,050
B72	UP823Y	SHIMADZU	Model UP823Y External Calibration, Capacity 820gm, Resolution 1mg/0.001g, PAN SIZE 108x105 mm	₹185,600
B73	UP1023Y	SHIMADZU	Model UP1023Y External Calibration Capacity 1020gm Resolution 1mg/0.001g PAN SIZE 108x105 mm	₹210,250

B74	UP422Y	SHIMADZU	Model UP422Y External Calibration Capacity 420gm Resolution 10mg/0.01g PAN SIZE 108x105 mm	₹120350
B75	UP822Y	SHIMADZU	Model UP822Y External Calibration Capacity 820gm Resolution 10mg/0.01g PAN SIZE 108x105 mm	₹150,800
B76	UP2202Y	SHIMADZU	Model UP2202Y External Calibration Capacity 2200gm Resolution 10mg/0.01g PAN SIZE 170x180 mm	₹134850
B77	UP4202Y	SHIMADZU	Model UP4202Y External Calibration Capacity 4200gm Resolution 10mg/0.01g PAN SIZE 170x180 mm	₹146,450
B78	UP6202Y	SHIMADZU	Model UP6202Y External Calibration Capacity 6200gm Resolution 10mg/0.01g PAN SIZE 170x180 mm	₹159500
B79	UP4201Y	SHIMADZU	Model UP4201Y External Calibration Capacity 4200gm Resolution 100mg/0.1g PAN SIZE 170x180 mm	₹143,550
B80	UP8201Y	SHIMADZU	Model UP8201Y External Calibration Capacity 8200gm Resolution 100mg/0.1g PAN SIZE 170x180 mm	₹155,150

We're your go-to source for all laboratory products.



Probe Sonicator



ice maker



ph meter

Upgrade your lab with top-quality equipment! Get the Probe Sonicator for efficient sample preparation, an Ice Maker for reliable cooling, and a PH Meter for accurate readings. Everything you need for precise lab work!

**BALANCE TW/TX/TXB SERIES - SHIMADZU**

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
TOP-LOADING BALANCES				
B81	TW223L	SHIMADZU	Model TW 223L Capacity 220gm Resolution 1mg/0.001g PAN SIZE Ø 110	₹130,500
B82	TW323L	SHIMADZU	Model TW 323L Capacity 320gm Resolution 1mg/0.001g PAN SIZE Ø 110	₹137,750
B83	TW423L	SHIMADZU	Model TW 423L Capacity 420gm Resolution 1mg/0.001g PAN SIZE Ø 110	₹143,550
B84	TX223L	SHIMADZU	Model TX 223L Capacity 220gm Resolution 1mg/0.001g PAN SIZE & Ø 110	₹108,750
B85	TX323L	SHIMADZU	Model TX 323L Capacity 320gm Resolution 1mg/0.001g PAN SIZE & Ø 110	₹110,200
B86	TX423L	SHIMADZU	Model TX 423L Capacity 420gm Resolution 1mg/0.001g PAN SIZE & Ø 110	₹111,650
B87	TX2202L	SHIMADZU	Model TX 2202L Capacity 2.2kg Resolution 1mg/0.001g PAN SIZE Ø 110	₹94425
B88	TX3202L	SHIMADZU	Model TX 3202L Capacity 3.2kg Resolution 1mg/0.001g PAN SIZE 110	₹97150
B89	TX4202L	SHIMADZU	Model TX 4202L Capacity 4.2kg Resolution 1mg/0.001g PAN SIZE Ø 110	₹101500
B90	TXB622L	SHIMADZU	Model TXB622L Capacity 620gm Resolution 10mg/0.01g PAN SIZE Ø 110	₹77,575
B91	TXB6201L	SHIMADZU	Model TXB6201L Capacity 6.2Kg Resolution 100mg/0.1g PAN SIZE Ø 160	₹77,575

B92	TWC323L	SHIMADZU	Model TWC323L Carat Scale internal calibration Capacity 320ct Resolution 0.001ct PAN SIZE Ø 80	POR
B93	TXC323L	SHIMADZU	Model TXC323L Carat Scale internal calibration Capacity 320ct Resolution 0.001ct PAN SIZE Ø 80	POR
B94	TXC623L	SHIMADZU	Model TXC623L Carat Scale Capacity 620ct Resolution 0.001ct PAN SIZE Ø 80	POR

YSP- SCIENTIFIC SOLUTION BALANCE ELB/BL SERIES - SHIMADZU

The BALANCE ELB/BL series by Shimadzu is designed for precision weighing in a variety of settings, including laboratories and industrial environments.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
PORTABLE ELECTRONIC BALANCES				
B95	ELB300	SHIMADZU	Model ELB300 Capacity 300 gm Resolution 0.01g (10mg) PAN SIZE 110 DIA	₹52,200
B96	ELB3000	SHIMADZU	Model ELB3000 Capacity 3KG Resolution 0.1g (100mg) PAN SIZE 170 x 130 DIA	₹52,200
B97	ELB12KH	SHIMADZU	Model ELB12KH Capacity 12KG Resolution 1g (1000mg) PAN SIZE 170 x 130 DIA	₹52,200
B98	BL220H	SHIMADZU	Model BL220H Capacity 2200 gm Resolution 0.01g(10mg) PANSIZE 160X124	₹62,350
B99	BL2200H	SHIMADZU	Model BL2200H Capacity 2200 gm Resolution 0.01g(10mg) PANSIZE 160X124	₹62,350



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Electronic Printer				
B100	EP100	SHIMADZU	Model EP-100 Printer with AC Adaptor	POR
B101	EP110	SHIMADZU	Model EP-110 Printer with AC Adaptor	₹105,850
B102	STABLOAP	SHIMADZU	Model STABLO-AP Ionizer-Speed AC Discharge for Static Electricity Removal	₹159,500
B103	ANMUSP100	SHIMADZU	Model ANM USP100 Printer (Equivalent to Shimadzu EP-110 printer)	₹50,750
B104	UDK	SHIMADZU	Universal Density kit	₹50,750

Shimadzu Balances Warranty - 1 year * IQ OQ PQ Doumnets charges extra



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
EMFC Based fully Auto Internal Calibration Analytical Balance				
High grade chemicals resistance Metallic housing				
B105	SES224	SAFFRON	Model SES224 Capacity 220g, Resolution 0.0001g, Pan Size 90mmØ, LCD Display, Auto internal calibration	₹63,438
B106	SES324	SAFFRON	Model SES324C Capacity 320g, Resolution 0.0001g, Pan Size 90mmØ, LCD Display, Auto internal calibration	₹88,450

Advanced Cutting edge Brilliant Technology

Competent products with High Quality Brand

EMFC- MONOLOG SINGLE BLOCK TECHNOLOGY (Analytical, Laboratory,Pharma,Chemicals,Precision, Institution, Gems-Jewelry)

B107	SES323S	SAFFRON	Model SES323S Capacity 320g, Resolution 0.001g, Pan Size 120mm∅;,LCD Display,Auto External calibration	₹36250
B108	SES2202	SAFFRON	Model SES2202 Capacity 2200g, Resolution 0.01g, Pan Size 150mmØ;,LCD Display,Auto External calibration	₹36,975



YSP- SCIENTIFIC SOLUTION

**BALANCE LOAD CELL BASE TECHNOLOGY
ANALYTICAL LAB - SAFFRON**



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B109	SES200LCD	SAFFRON	Model SES200 (LCD) Capacity 200g, Resolution 0.001g, Pan Size 80mmØ,LCD Display, External calibration	₹17,063
B110	SES300LCD	SAFFRON	Model SES300 (LCD) Capacity 300g, Resolution 0.001g, Pan Size 80mmØ,LCD Display, External calibration	₹14,718
B111	SES400LCD	SAFFRON	Model SES400 (LCD) Capacity 400g, Resolution 0.001g, Pan Size 80mmØ,LCD Display, External calibration	₹19,285
B112	SES500LCD	SAFFRON	Model SES500 (LCD) Capacity 500g, Resolution 0.001g, Pan Size 80mmØ,LCD Display, External calibration	₹21,750

B113	SES605LCD	SAFFRON	Model SES605 (LCD) Capacity 600g, Resolution 0.005g, Pan Size 120mmØ,LCD Display, External calibration	₹15,050
B114	SES1005LCD	SAFFRON	Model SES1005 (LCD) Capacity 1000g, Resolution 0.005g, Pan Size 120mmØ,LCD Display, External calibration	₹14,718
B115	SES602LCD	SAFFRON	Model SES602 (LCD) Capacity 600g, Resolution 0.01g, Pan Size 120mmØ,LCD Display, External calibration	₹13,738
B116	SES1202LCD	SAFFRON	Model SES1202 (LCD) Capacity 1200g, Resolution 0.01g, Pan Size 120mmØ,LCD Display, External calibration	₹16,450
B117	SES2002LCD	SAFFRON	Model SES2002 (LCD) Capacity 2000g, Resolution 0.01g, Pan Size 120mmØ,LCD Display, External calibration	₹15,878
B118	SES3002LCD	SAFFRON	Model SES3002 (LCD) Capacity 3000g, Resolution 0.01g, Pan Size 170X170mm,LCD Display, External calibration	₹22,693
B119	SES5002LCD	SAFFRON	Model SES5002 (LCD) Capacity 5000g, Resolution 0.01g, Pan Size 170X170mm,LCD Display, External calibration	₹34,003
B120	SES6002LCD	SAFFRON	Model SES6002 (LCD) Capacity 6000g, Resolution 0.01g, Pan Size 170X170mm,LCD Display, External calibration	₹38,570
B121	SES6001LCD	SAFFRON	Model SES6001 (LCD) Capacity 6000g, Resolution 0.1g, Pan Size 170X170mm,LCD Display, External calibration	₹18,125
B122	SES14LCD	SAFFRON	Model SES14 (LCD) Capacity 10g, Resolution 0.0001g, Pan Size 45mmØ,LCD Display, External calibration	₹31,755
B123	SES24LCD	SAFFRON	Model SES24 (LCD) Capacity 20g, Resolution 0.0001g, Pan Size 45mmØ,LCD Display, External calibration	₹36,250

BALANCE LOAD CELL BASE PRECISION SCALE - SAFFRON



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
External Calibration				
B124	SES202SLED	SAFFRON	Model SES202S (LED) Capacity 200g, Resolution 0.01g, Pan Size 120mmØ, Red LED Display	₹7,875
B125	SES302SLED	SAFFRON	Model SES302S (LED) Capacity 300g, Resolution 0.01g, Pan Size 120mmØ, Red LED Display	₹9,625
B126	SES602SLED	SAFFRON	Model SES602S (LED) Capacity 600g, Resolution 0.01g, Pan Size 120mmØ, Red LED Display	₹9,625
B127	SES1202SLED	SAFFRON	Model SES1202S (LED) Capacity 1200g, Resolution 0.01g, Pan Size 120mmØ, Red LED Display	₹12,075
B128	SES302LED	SAFFRON	Model SES302 (LED) Capacity 300g, Resolution 0.01g, Pan Size 120mmØ, Green LED Display	₹10,938
B129	SES602LED	SAFFRON	Model SES602 (LED) Capacity 600g, Resolution 0.01g, Pan Size 120mmØ, Green LED Display	₹10,938
B130	SES1202LED	SAFFRON	Model SES1202 (LED) Capacity 1200g, Resolution 0.01g, Pan Size 120mmØ, Green LED Display	₹13,825
B131	SES605LED	SAFFRON	Model SES605 (LED) Capacity 600g, Resolution 0.005g, Pan Size 120mmØ, Green LED Display	₹13,825

B132	SES1005LED	SAFFRON	Model SES1005 (LED) Capacity 1000g, Resolution 0.005g, Pan Size 120mmØ, Green LED Display	₹15,050
B133	SES2002LED	SAFFRON	Model SES2002 (LED) Capacity 2000g, Resolution 0.01g, Pan Size 120mmØ, Green LED Display	₹15,950

**Complete wash down SS 304 Table Top, Bench scale for
Pharmaceutical & Chemistry Industry**

B134	SES6THSS	SAFFRON	Model SES6THSS Capacity 6Kg, Resolution 0.1g, Pan Size 25mm, LCD/LED Display	₹27,188
B135	SES10THSS	SAFFRON	Model SES10THSS Capacity 10Kg, Resolution 0.1g, Pan Size 25mm, LCD/LED Display	₹31,755
B136	SES20THSS	SAFFRON	Model SES20THSS Capacity 20Kg, Resolution 0.1g, Pan Size 25mm, LCD/LED Display	₹36,250
B137	SES30THSS	SAFFRON	Model SES30THSS Capacity 30Kg, Resolution 0.1g, Pan Size 25mm, LCD/LED Display	₹40,818
B138	SES50THSS	SAFFRON	Model SES50THSS Capacity 50Kg, Resolution 1g, Pan Size 30mm, LCD/LED Display	₹40,818



Durable and easy to clean, our SS 304 Table Top and Bench Scale are perfect for the pharmaceutical and chemistry industries. Reliable precision for your lab needs!



BALANCE EMFC 5 DIGIT SEMI MICRO - SAFFRON

(Analytical, Laboratory, Pharma, Chemicals, Precision, Insitution)



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B139	SES45	SAFFRON	Model SES45 Capacity 40g, Resolution 0.01mg, Pan Size 90mmØ,LCD Display, Auto internal calibration	₹208,438
B140	SES65	SAFFRON	Model SES65 Capacity 60g, Resolution 0.01mg, Pan Size 90mmØ,LCD Display, Auto internal calibration	₹213,005
B141	SES165	SAFFRON	Model SES165 Capacity 60/100g, Resolution 0.01/0.1mg, Pan Size 90mmØ,LCD Display	₹217,500
B142	SES265	SAFFRON	Model SES265 Capacity 60/220g, Resolution 0.01/0.1mg, Pan Size 90mmØ,LCD Display	₹213,005
B143	SES285	SAFFRON	Model SES285 Capacity 80/220g, Resolution 0.01/0.1mg, Pan Size 90mmØ,LCD Display,	₹237,945

EMFC- MONOLOG SINGLE BLOCK TECHNOLOGY

(Analytical, Laboratory, Pharma, Chemicals, Precision, Institution, Gems-Jewelry)

Auto External calibration

B144	SES220C	SAFFRON	Model SES220C Capacity 220g, Resolution 0.0001g, Pan Size 90mmØ,LCD Display	₹95,193
B145	SES623	SAFFRON	Model SES623 Capacity 620g, Resolution 0.001g, Pan Size 120mmØ,LCD Display	₹70,253
B146	SES3202	SAFFRON	Model SES3202 Capacity 3200g, Resolution 0.01g, Pan Size 220X220mm,LCD Display	₹47,705
B147	SES6202	SAFFRON	Model SES6202 Capacity 6200g, Resolution 0.01g, Pan Size 220X220mm,LCD Display	₹108,750



BALANCE TABLE TOP SCALE LOAD CELL TECHNOLOGY DP MOQ5 - SAFFRON

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B148	SES6TH	SAFFRON	Model SES6TH Capacity 6Kg, Resolution 0.1g, Pan Size 240x290, LED Display, External calibration	₹15,050
B149	SES3T	SAFFRON	Model SES3T Capacity 3Kg, Resolution 0.1g, Pan Size 170x220, LED Display, External calibration	₹6,038
B150	SES6T	SAFFRON	Model SES6T Capacity 6Kg, Resolution 0.5g, Pan Size 240x290, LED Display, External calibration	₹6,038
B151	SES10T	SAFFRON	Model SES10T Capacity 10Kg, Resolution 0.5g, Pan Size 240x290, LED Display, External calibration	₹6,038
B152	SES20T	SAFFRON	Model SES20T Capacity 20Kg, Resolution 1g, Pan Size 240x290, LED Display, External calibration	₹6,038
B153	SES30T	SAFFRON	Model SES30T Capacity 30Kg, Resolution 1g, Pan Size 240x290, LED Display, External calibration	₹6,038

We specialize in the sales of all laboratory products.

For inquiries or orders,



96262 97366
73387 21875



BALANCE OIML - STANDARD WEIGHT BOX & INDIVIDUAL WEIGHTS - SAFFRON

For Weight Box Please Refer "w"



BREATH ALCOHOL DETECTOR [Not serviceable] (Battery Operated)

B154		GENERAL	Model AT 6000 Breath Alcohol meter China.	₹2,625
------	--	---------	-------------------------------------------	--------

BLOOD CELL COUNTER - LABOHOLIC



B155	LH21	LABOHOLIC	LH-21 Differential Blood Cell Counter, a microcontroller based system for quick differential counts, has agility to count 8 types of cell and audible beep facility, 16 digit alphanumeric LCD Display shows percentage (%) collations, individual cell count with built-in stop watch and timer functions.	₹8,750
------	------	-----------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------



BLOOD ROLLER MIXER - LABOHOLIC

B156	LH43	LABOHOLIC	LH-43 Blood Roller Mixer to make blood homogeneous with 4 rollers 230mmlong.	₹6,694
------	------	-----------	------------------------------------------------------------------------------	--------

BLOOD CELL COUNTER



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B157		GENERAL	Blood cell counters with 5 keys and one totalizer with bell	₹1,488
B158		GENERAL	Blood cell counters with 8 keys and one totalizer with bell	₹2,100
B159		GENERAL	Blood cell counter Digital Economy Model with 8 keys & totalizer china make	₹4,375

BOTTLE TOP DISPENSER FULLY AUTOCLAVABLE - SCITUS Microlit (New Model)



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B160	SCI25	MICROLIT	SCI-2.5 Range 0.25 - 2.5 ml Microlit	₹15,341
B161	SCI5	MICROLIT	SCI-5 Range 0.5 - 5 ml Microlit	₹15,341
B162	SCI10	MICROLIT	SCI-10 Range 1-10 ml Microlit	₹15,341
B163	SCI30	MICROLIT	SCI-30 Range 2.5 - 30 ml Microlit	₹22,939
B164	SCI60	MICROLIT	SCI-60 Range 5- 60 ml Microlit	₹22,939
B165	SCI100	MICROLIT	SCI-100 Range 10- 100 ml Microlit	₹27,623

BOTTLE TOP DISPENSER WITH RE-CIRCULATION VALVE - BEATUS Microlit (New Model)



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B166	BEAT25	MICROLIT	BEAT-2.5 Range 0.25 - 2.5 ml Microlit	₹20,721
B167	BEAT5	MICROLIT	BEAT-5 Range 0.5 - 5 ml Microlit	₹20,721
B168	BEAT10	MICROLIT	BEAT-10 Range 1 -10 ml Microlit	₹20,721
B169	BEAT30	MICROLIT	BEAT-30 Range 2.5 - 30 ml Microlit	₹30,987
B170	BEAT60	MICROLIT	BEAT-60 Range 5 - 60 ml Microlit	₹30,987
B171	BEAT100	MICROLIT	BEAT-100 Range 10 - 100 ml Microlit	₹37,338

BOTTLE DISPENSER - TOP TECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B172	T15	TOP TECH	Bottle Top Dispenser 0-5 ml	₹7,875
B173	T16	TOP TECH	Bottle Top Dispenser 2-10 ml	₹7,875

BURSTING STRENGTH TESTER DIGITAL - RICO



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B174	RSBT 2	RICO	RICO RSBT- 2 bursting strength tester (digital model) capacity – 0 to 35kg / cm ² , to determine bursting strength of knitted fabrics, paper, corrugated paper etc. latest design with control panel with maximum capacity 35 kgs. glycerin used as test fluid	₹162,400

BUFFER PH SOLUTIONS - EUTECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Buffer Solutions 480 ml Bottle Calibrated with Certificate of analysis from EUTECH				
B175	PH168	EUTECH	pH Buffer Solution 1.68pH	₹3,063
B176	PH401	EUTECH	pH Buffer Solution 4.01pH	₹3,063
B177	PH686	EUTECH	pH Buffer Solution 6.86pH	₹3,063
B178	PH7	EUTECH	pH Buffer Solution 7pH	₹3,063
B179	PH9	EUTECH	pH Buffer Solution 9pH	₹3,063
B180	PH918	EUTECH	pH Buffer Solution 9.18pH	₹3,063
B181	PH1001	EUTECH	pH Buffer Solution 10.01pH	₹3,063
B182	PH1245	EUTECH	pH Buffer Solution 12.45pH	₹3,063

BUFFER TDS SOLUTIONS - EUTECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Buffer Solutions 480 ml Bottle Calibrated with Certificate of analysis from EUTECH				
B183	TDS50	EUTECH	TDS Buffer Solution 50ppm	₹3,238
B184	TDS300	EUTECH	TDS Buffer Solution 300ppm	₹3,238
B185	TDS1000	EUTECH	TDS Buffer Solution 1000ppm	₹3,238
B186	TDS3000	EUTECH	TDS Buffer Solution 3000ppm	₹3,238
B187	TDS5000	EUTECH	TDS Buffer Solution 5ppt (5000ppm)	₹4,725
B188	TDS25000	EUTECH	TDS Buffer Solution 25ppt (25000ppm)	₹4,725
B189	TDS45000	EUTECH	TDS Buffer Solution 45ppt (45000ppm)	₹4,725

BUFFER ORP SOLUTION - EUTECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Buffer Solutions 480 ml Bottle Calibrated with Certificate of analysis from EUTECH				
B190	ORP86MV	EUTECH	Buffer Solutions for ORP 86 mv	₹7,875
B191	ORP263MV	EUTECH	Buffer Solutions for ORP 263mv	₹7,875
B192	ORP475MV	EUTECH	Buffer Solutions for ORP 475mv	₹7,875

BUFFER CONDUCTIVITY SOLUTIONS - EUTECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Buffer Solutions 480 ml Bottle Calibrated with Certificate of analysis from EUTECH				
B193	COND84	EUTECH	Conductivity Buffer Solution 84 μ s	₹3,238
B194	COND100	EUTECH	Conductivity Buffer Solution 100 μ s	₹3,238
B195	COND500	EUTECH	Conductivity Buffer Solution 500 μ s	₹3,238
B196	COND1413	EUTECH	Conductivity Buffer Solution 1413 μ s	₹3,238
B197	COND2764	EUTECH	Conductivity Buffer Solution 2764 μ s	₹3,238
B198	COND5	EUTECH	Conductivity Buffer Solution 5 ms	₹3,238
B199	COND12880	EUTECH	Conductivity Buffer Solution 12880 μ s	₹3,456
B200	COND1118	EUTECH	Conductivity Buffer Solution 111.8 ms	₹3,456

BUFFER CONDUCTIVITY SOLUTIONS - HAMILTON



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B201	COND1.3US	HAMILTON	HAMILTON COND BUFFER SOLUTION 1.3 US WITH NIST CERTIFICATE	₹23,128
B202	COND5US	HAMILTON	HAMILTON COND BUFFER SOLUTION 5 US WITH NIST CERTIFICATE	₹21,678

BUFFER ELECTRODE STORAGE SOLUTION - EUTECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Buffer Solutions 480 ml Bottle Calibrated with Certificate of analysis from EUTECH				
B203	ECRE005	EUTECH	ELECTRODE STORAGE SOLUTION EUTECH 480ML ECRE005	₹14,175

BUFFER SOLUTION - AQUASOL



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Buffer Solutions 480 ml Bottle Calibrated with Certificate of analysis from EUTECH				
B204	AQUASOL PH BUFFER SOLUTION	AQUASO L	Aquasol Buffer pH 4, 7, 10, 6.86, 9.18, 12.45 each 500ml bottle	₹700
B205	AQUASOL PH 1.68 BUFFER SOLUTION	AQUASO L	Aquasol Buffer pH 1.68 each 500ml bottle	₹1,269
B206	AQUASOL TURBIDITY BUFFER SOLUTION	AQUASO L	Aquasol Turbidity NTU 0, 10, 20, 50, 100, 500, 800, 1000, 2000, 4000 each 100ml bottle	₹1,269

BIO CHEMISTRY ANALYSER - LABOHOLIC

Bulk Density Apparatus Refer "P"



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Bio Nano Spectrophotometer - Laboholic Refer "S"				
B207	LH1201	LABOHOLIC	LH-1201 Auto Bio Chemistry Analyser For End Point, Kinetic, Fix Point, Multi Standard and coagulation tests, Quartz flowcell system, Built in Thermal printer with graphics capability, Upto 250test items, 240 x 128 LCD display, Paltrier Controlled thermostat, open reagent system, user friendly. (Imported)	₹176,900
B208	LH1202	LABOHOLIC	LH-1202 Semi Automatic Bio Chemistry Analyser, Real time coloured graph and absorbance, Windows Embedded compact 7, 4.3 inch true color touchscreen, 300 programmable test, 1000 Result storage with editable patient information, universal power supply, End point with standard/factor, Two point kinetics, Multi standard upto 8 standards, biochromatic, dual wavelength, absorbance, photometric analysis with interference filter, aborbance range 0.0000-3.0000 (Linearity:+/-1%, Repeatability: 0.002, Resolution: 0.0001), Spectral range selection: Automatic via Filter wheel, Light source: 6V 10 W Long life tungsten Halogen lamp, Detector: Photodiode.	₹152,250

**BLOOD BANK INSTRUMENTS
(AGAINST ORDER 2 TO 3
WEEKS TIME) - GLOLAB**



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
B209	GLBTS42	GLOLAB	GLBTS-4.2 Blood Bag Tube Sealer is a compact equipment to seal the Blood Bag pilot tube by radio frequency, Sealed tube can be easily separated by pulling the sealing section by both sides, Scan Here For Detailed Specifications, Oscillation Frequency : 40.68 MHz, Maximum Sealing Time - 1.5 seconds, Radio Frequency Output :20W	₹108,750
B210	GLP1	GLOLAB	GLP-1 Plasma Expressor, Pressure plate: Acrylic, Manual control, size: 232 x 165 x 243 mm	₹17,325
B211	GLP2	GLOLAB	GLP-2 Electronic Plasma Expressor, Power supply 230v/110v AC, Infra Red sensor, Clamping: Motor activated, Acrylic Pressure plate, alarm: audio visual, automatic control, Size 252 x 165 x 243 mm	₹71,775
B212	GLBM14	GLOLAB	GLBM - 14 Blood Collection Monitor, Control System : Micro controller based, Alarm : Audio Visual Alarm if the flow is less than 20 ml per minute at the end of the Collection, Outer body : Mild steel powder coated, Speed of rocker : 10 RPM, Accuracy : ± 2 RPM, Settable volume : 1 to 999 ml Readability : 1 ml/1g, Pinching : Motorized pinching at the end of collection, Display : user friendly LCD display with messages, Blood bag Tray : ABS molded, Display parameters : Weight and volume of collection, duration of collection Previous collection records, Calibration : Auto calibration from the control panel	₹108,750
B213	GLAP48	GLOLAB	GLAP-48 Platelet Agitator is very well-known for its uniqueness in design. S.S racks with high ventilation capability allows maximum oxygen transfer to the platelets. Sliding stainless steel trays makes the cleaning purpose easier as well as enables easy removal or loading of bags. Capacity 48 bags, No. of tray: 8, Oscillation: 70/Minute, Accuracy: ± 5/ minute"	₹97,875



CONDUCTIVITY METER PEN TYPE -EUTECH/ HANNA / HM

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C1	ECOECTESTLOW	EUTECH	Model ECOEC TEST LOW Range 0-1999 μ S/cm, NO TEMP.Display, Resolution 10 μ S/cm, Sensor Pin	₹12,128
C2	ECOECTESTHIGH	EUTECH	Model ECOEC TEST HIGH Range 0-19.9MS, NO TEMP.Display, Resolution 0.1ms/cm, Sensor Pin	₹12,128
C3	ECTEST11	EUTECH	Model EC TEST 11 Range 0-2000 μ S/cm & 0-20MS/cm, TEMP.Display Yes, Resolution 10 μ S/cm, Sensor Pin	₹15,791
C4	ECTEST11PLS	EUTECH	Model EC TEST 11 PLS Range 200,2000 μ S/cm,20MS/cm, TEMP.Display YES, Resolution 0.1 μ S/cm, Sensor Cup	₹24,795
C5	HI98303	HANNA	Conductivity Meter Dist 3 Range 0-1990 μ S/cm, Resolution 1 μ S/cm, Temp.BOLD Display HI98303 New	₹12,031
C6	HI98304	HANNA	Conductivity Meter Dist 4 Range 0.01ms TO 19.99ms, Resolution .01ms, Temp.BOLD Display HI98304 New	₹12,031
C7	AQ98303P	ACETEQ	Aceteq Conductivity Meter Dist 3 Range 0-1990 μ S/cm, Resolution 1 μ S/cm, Temp.BOLD Display AQ98303P New	₹2,625
C8	AQ98304P	ACETEQ	Aceteq Conductivity Meter Dist 4 Range 0.01ms TO 19.99ms, Resolution .01ms, Temp.BOLD Display AQ98304P New	₹2,625
C9	AP2	HM	AP 2 Conductivity Meter Range 0- 9999 μ S/cm, Resolution 1 μ S/cm, TEMP.Small Display	₹1,313
C10	COM80	HM	COM80 Model COND /TDS HM Model Range 9990 μ S/cm/5000ppm, Resolution 10 μ S/cm/10ppm, TEMP.BOLD Display	₹2,800

CONDUCTIVITY METER PEN TYPE - EUTECH/ HANNA / HM




SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C11	CON603PLUSK	EUTECH	Model CON603PLUSK Parameter Conductivity Range 0 to /200MS/cm, Powered by battery IP54	₹57,420
C12	CONWP15003K	EUTECH	Model CONWP15003K Parameter Cond / TDS Range 0-200ms/200ppt, Memory 150 Powered by Main/Battery	₹71,775
C13	CONWP45003K	EUTECH	Model CONWP45003k Parameter Cond/TDS/Salinity, Range 200ms/200ppt/50ppt, Memory 500 Powered by Main/Battery	₹109,620
C14	CONWP61043K	EUTECH	Model CONWP61043K Parameter Cond/TDS/Res/Sal, Range 0 to 500MS/cm 500ppt, Memory 500 Powered by Main/Battery IP 67	POR

CONDUCTIVITY METER PORTABLE - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C15	LH17	LABOHOLIC	LH-17 Portable Conductivity Meter 3 1/2 digit LCD display, battery operated, ABS Plastic body.	₹7,525

CONDUCTIVITY METER BENCH MODEL - EUTECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C16	ECCON70043S	EUTECH	Model ECCON70043S Conductivity/ATC Probe CONSEN9501D ,K=1.0, Integral Electrode Holder & 100/240 VAC Adapter Cond/TDS, Range 0 to 200ms/cm:0 to 100 ppt, Powered by Main	₹95,265
C17	ECCON270043S	EUTECH	Model ECCON270043S With 4-Cell Conductivity Probe CONSEN9201D K=0.530, Integral Electrode Holder & 100/240 VAC Adapter, RS232 cableParameter Cond/TDS/Salinity/Restivity and Tepm.Range 0 to 500ms/cm 0 to 500ppt; 0 to 80 ppt,o to 200mohms	₹152,975
C18	CON2700AN270	EUTECH	ANTECH AN 270 Printer Suitable for ECCON270043S model .	₹52,200

EUTECH - PH & CONDUCTIVITY METER COMBO OFFER



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C19	CONPH700	EUTECH	PH 700 & CON 700	₹57,420

CONDUCTIVITY METER TABLE MODEL - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C20	LH5100	LABOHOLIC	LH-5100 Trutouch Conductivity Meter, 5 Inch touch screen with high resolution and sensitivity, Combined and Separate measurement facility cond. TDS. Sal. and temp., 6 auto ranges, Auto sensor check facility, 5 point calibration facility 1 for each range(EC-0-20 , 0-200S μ s, 0-2ms,0-20ms,0-200ms, TDS - 0-20, 0-200ppm, 0-20,0-200ppt, Salinity- 0-20,0-200ppm, 0-0,2-20,0-5ppt), Meets GLP norms, Password protected data and login facility, PC attachment through USB.,1000 data storing facility with file name. NEW ADDITION	₹48,648
C21		LABOHOLIC	EPSON LX-310 Printer Suitable for Trutouch Conductivity Meter(LH-5100) & Trutouch pH Meter(LH-5001)	₹22,475
C22	LH16	LABOHOLIC	LH-16 Auto Digital Conductivity Meter (AUTO RANGING) 3.5 digit LED display with facility for temperature compensation and digital cell constant adjustment along with conductivity cell & other accessories.	₹10,675
C23	LH26	LABOHOLIC	LH-26 Deluxe Conductivity Meter 3.5 digit LED display with facility for automatic temp. Compensation & digital cell constant adjustment along with cond. cell & temperature probe	₹11,743
C24	LH-601	LABOHOLIC	LH-601 Deluxe Conductivity Meter 3.5 Digit LED Display with Stirrer facility	₹11,813



MULTIPARAMETERS COND/PH/TDS COMBO METERS PEN/PORTABLE/BENCH MODEL

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Please Refer Page "PMultiparameter				
CONDUCTIVITY/TDS/SAL/PH/TEMP MULTIPARAMETER TABLE MODEL - LABOHOLIC				
C25	LH51	LABOHOLIC	<p>LH-51 Microprocessor Conductivity/TDS/Salinity/Temp. Meter (Graphical Display) a research grade instrument with 128 x 64 Large Graphical LCD display having automatic calibration, automatic range selection, automatic Temperature compensation, selectable reference temp, fluid resistant housing & rotational electrode stand along with accessories. Data Storage fealty up to 1000 samples, complete with cell. temperature probe, rotational electrode stand, dust cover& instructions manual, computer interlace through USB</p>	
C26	LH53	LABOHOLIC	<p>LH-53 Microprocessor pH/Conductivity/TDS/Salinity/Temp. Meter a research grade instrument with 16 x 2 Large LCD display having automatic calibration, automatic range selection, automatic Temperature compensation, selectable reference temp, fluid resistant housing & rotational electrode stand along with accessories. Data Storage fealty up to 1000 samples, complete with cell. temperature probe, rotational electrode stand, dust cover & instructions manual, computer interlace through USB NEW</p>	

C27	LH23	LABOH OLIC	LH-23 Digital Ph, Conductivity and Temperature Meter for measurement of pH. Conductivity & temp., 3.5 digit LED display supplied with pH combination electrode, conductivity cell and temp.probe	
C28	LH25	LABOH OLIC	LH-25 Digital TDS/Conductivity Meter Digital Conductivity Cum TDS meter 3.5 digit LED display with facility for temp. compensation & digital cell constant adjustment along with conductivity cell	


 **YSP- SCIENTIFIC SOLUTION**

CONDUCTIVITY / TDS CELL SPARE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C29	SWINGARM STAND	YSP	Laboholic swing arm stand with good base at best rate	₹1,313
C30	CONSEN91 W2SS	EUTECH	Cond.cell CONSEN 91W 2SS rings with ATC, Cell constant K=1.0, 6 pin connector, 1M cable for CON 510 / 150/ 450 / 11 / 110	₹39,150
C31	CONSEN91B	EUTECH	Cond Cell CONSEN91B Suitable FOR SALT 6+ / Con 6+ / TDS 6+	₹39,150
C32	CONSEN95 01D	EUTECH	Cond Cell CONSEN9501D COND / ATC Electrode K=1.0 Suitable for Con700 Models	₹39,150
C33	CONSEN92 01D	EUTECH	Cond Cell CONSEN9201D 4-Cell Conductivity/ATC Electrode K=0.530, Graphite Sensor, Epoxy-Body (2700 series)	₹78,300
C34	93X546101	EUTECH	93X546101 2-Cell Conductivity/ATC Electrode, 0.1 cell constant, SS Sensor, Epoxy-Body (700 & 2700 series)	₹49,590

COLONY COUNTER - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C35	LH36	LABOHOLIC	LH-36 Digital Colony Counter 3 digit LED display Range 0-999 with 125 memory, hold and count correction facility. Data restored even after powerfailure.	₹8,006
C36	LH37	LABOHOLIC	LH-37 Digital ColonyCounter4 digit LED display Range 0-9999. With 125 memory, hold and count correction facility. Data restored even after powerfailure	₹8,540
C37	LH-376	LABOHOLIC	LH-376 Digital Colony Counter 6 digit LED display Range 0-999 with 125 memory, hold and count correction facility. Data restored even after powerfailure.	₹13,125

COLORIMETER - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C38	LH116	LABOHOLIC	LH-116 Digital Fully Automatic Colorimeter 9 digital filters with facility for automatic zero setting with finger touch, digital wavelength selection. mode OD inbuilt 6V Ni-Cd Re-Chargeable Battery. 1 ml. solution measurement, Range400 to 700 nm, highly accurate & compact.	₹11,209
C39	LH111	LABOHOLIC	LH-116 Digital Fully Automatic Colorimeter 9 digital filters with facility for automatic zero setting with finger touch, digital wavelength selection. mode OD inbuilt 6V Ni-Cd Re-Chargeable Battery. 1 ml. solution measurement, Range400 to 700 nm, highly accurate & compact.	₹12,031

C40	LH114	LABOHOLIC	LH-114 Auto Photo Colorimeter 8 filters with facility for automatic zero setting with finger touch, 2.5 digit LED display, Range 400to 700nm, highly accurate & compact.	₹10,675
C41	LH112	LABOHOLIC	LH-112 Deluxe Photo Colorimeter 8 filters mains cum battery operated with inbuilt rechargeable battery. 2.5 digit LED display, Range 400 to 700nm, highly accurate & compact.	₹10,675
C42	LH12	LABOHOLIC	LH-12 Digital Photo Colorimeter 8 filters mains operated with optical density 2.5 digit LED display. Range 400 to 700 nm. lml solution measurement.	₹9,888

 **YSP- SCIENTIFIC SOLUTION**

COLORIMETER - LABOHOLIC



C43	LH790	LABOHOLIC	LH-790 Visual Colorimeter Measuring principle visual in terms of international colour scale units, measuring range 0.1-79.9 Red, 0.1 - 79.9 Yellow, 0.1- 49.9 Blue, 0.1 -3.9 Neutral, Resolution 0.1 colour scale unit, Light source Incandescent lamp, path length upto 180mm.	₹174,000
-----	-------	-----------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------

 **YSP- SCIENTIFIC SOLUTION**

COAGULATION ANALYZER - LABOHOLIC



C44	LH1232	LABOHOLIC	LH-1232 Coagulation Analyzer, 5 inch LCD Touch screen, 4 Testing channels, 24 Cuvette incubation positions, 6 Reagent incubation positions, Wavelength LED 470nm, Sample Type Plasma, 10,000 results storage memory, Sample volume 20-40ul, Reagent volume 20-40ul, Electronic pipette 20-200ul, Incubation Temp. 37.0± 1.0°C NEW ADDITION	₹229,100
-----	--------	-----------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------

CENTRIFUGE - YSP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C45	YSPDOCTOR	YSP	YSP-DOCTOR CENTRIFUGE YSP Doctor Centrifuge 8 x 15 ml Metal body Head ABS	₹5,250

CENTRIFUGE - DLAB



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C46	D2012 PLUS	DLAB	Cat.No.9032002121 D2012 Plus – 15,000 RPM High Speed Mini Centrifuge (CE marked), with A12-2P plastic Rotor	₹130,500
C47	D1008	DLAB	Cat.No.9031001012 D1008 - Mini Centrifuge with blue lid including two rotors A08-2 & A02-PCR8 and both adapter SA02P2 & SA05P2	₹20,300
C48	D1524R	DLAB	Cat.No.9013111121 D1524R High Speed Refrigerated Micro-Centrifuge, 15000 rpm with automatic chamber lock, customizable programs, adjustable acceleration/deceleration and Temp Range -20 to 40'C other upgrades with AS24-2 (2 ml / 1.5 ml x 24) aluminium alloy rotor kit 19400006	₹565,500

CENTRIFUGE MINI / MICRO - REMI



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C49	RM02PLUS	REMI	RM-02 PLUS Remi Mini Plus Centrifuge 6000 RPM With Speed Indicator, Capacity 8 x 1.5 ml	₹15,225
C50	RM03PLUS	REMI	RM-03 PLUS Remi Micro Centrifuge 15000 RPM LCD Display with 12x1.5ml /2 ml and PCR strip Rotor (2 x 8 x 0.2) ml	₹79,750

CENTRIFUGE MEDICO PLUS - REMI



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C51	MEDICOPLUS	REMI	Medico Plus Medico Centrifuge max.speed 4000 rpm with brushless induction motor, frequency drive, LED display of speed and time, safety lid lock, imbalance detector, dynamic brake supplied without rotors and accessories	₹37,700
C52	SM815	REMI	SM 8-15 Swing out rotor suitable for 8 tubes of 15ml (supplied with glass tubes)	₹8,400
C53	SM450	REMI	SM 4-50 Swing out rotor suitable for 4 tubes of 50ml (supplied with glass tubes)	₹7,350
C54	AM365	REMI	AM 36-5 Angle rotor suitable for 36 Vacutainer tubes of 5ml (12.5mm x 81mm) without tubes	₹8,400
C56	AM1810	REMI	AM 18-10 Angle rotor suitable for 18 Vacutainer tubes of 10ml (15.5mm x 106mm) without tubes	₹7,700

C57	AM1215	REMI	AM 12-15 Angle rotor suitable for 12 tubes of 15ml (supplied with glass tubes)	₹7,700
C58	AM650	REMI	AM 6-50 Angle rotor suitable for 6 tubes of 50ml (supplied with plastic tubes)	₹9,625

 **YSP- SCIENTIFIC SOLUTION**

CENTRIFUGE - REMI



C59	C852	REMI	C-852 Medical Centrifuge Buckets 4 Capacity 15ml RPM 3500 NO Display	₹12,075
C60	C854	REMI	C-854/4 Clinical Centrifuge Bucket 4 Capacity 15ml RPM 3500 NO Display	₹14,000
C61	C854 6	REMI	C-854/6 Clinical Centrifuge Bucket 6 Capacity 15ml RPM 3500 NO Display	₹15,400
C62	C8548	REMI	C-854/8 Clinical Centrifuge Bucket 8 Capacity 15ml RPM 3500 NO Display	₹16,240
C63	R303	REMI	R-303 Doctors Centrifuge Bucket 8 Capacity 15ml RPM 3800 NO Display	₹15,750
C64	R304	REMI	R-304 Doctors Centrifuge Bucket 24 Capacity 5ml RPM 7000 NO Display	₹20,010
C65	R4C	REMI	R-4C W/Timer Centrifuge Without Head RPM 4400 WITH Display	₹37,700
C66	R41	REMI	Head R-41 Spare Bucket 8 (SWING) Capacity 15ml RPM 3500, Suitable for R-4C	₹7,525
C67	R42	REMI	Head R-42 Spare Bucket 4 (SWING) Capacity 50ml RPM 3500 Suitable for R-4C	₹9,275
C68	R43	REMI	Head R-43 Spare Bucket 8 (ANGLE) Capacity 15ml RPM 4400 Suitable for R-4C	₹7,525
C69	R44	REMI	Head R-44 Spare Bucket 4 (ANGLE) Capacity 50ml RPM 4200 Suitable for R-4C	₹9,450

C70	R8C	REMI	R-8C Centrifuge Without Head RPM 5250 WITH Display	₹44,950
C71	R8CBL	REMI	R8C BL Centrifuge Without Head RPM 6000 WITH Display	₹104,400
C72	R81	REMI	Head R-81 Spare Bucket 16 (swing) Capacity 15ml RPM 4000 Suitable for R-8C	₹14,000
C73	R81A	REMI	Head R-81A Spare Bucket 12(swing) Capacity 15ml RPM 4000 Suitable for R-8C	₹12,250
C74	R81B	REMI	Head R-81B Spare Bucket 8(swing) Capacity 15ml RPM 4200 Suitable for R-8C	₹7,525
C75	R82	REMI	Head R-82 Spare Bucket 6 (swing) Capacity 50ml Suitable for R-8C	₹12,600
C76	R82A	REMI	Head R-82A Spare Bucket 4 (swing) Capacity 50ml RPM 4000 Suitable for R-8C	₹9,275
C77	R83	REMI	Head R-83 Spare Bucket 16 (Angle) Capacity 15ml Suitable for R-8C	₹12,250
C78	R83A	REMI	Head R-83A Spare Bucket 12 (Angle) Capacity 15ml RPM 5250 Suitable for R-8C	₹9,625
C79	R83B	REMI	Head R-83B Spare Bucket 8 (Angle) Capacity 15ml with polypropylene tubes suitable for R- 8C	₹7,700
C80	R84	REMI	Head R-84 Spare Bucket 6 (Angle) Capacity 50ml with polypropylene tubes Suitable for R- 8C	₹7,700
C81	R84A	REMI	Head R-84A Spare Bucket 4 (Angle) Capacity 50ml RPM 5250 Suitable for R-8C	₹7,700
C82	R88	REMI	Head R-88 Spare Angle Head 4 x 100 ml with Polypropylene tube For R-8C BL	₹17,500
C83	R89	REMI	Head R-89 Spare Angle Head 24 x 15 ml with Polypropylene tube For R-8C BL	₹15,225
C84	R90	REMI	Head R-90 Spare 4 x 100 ml swing out with glass tube & Biosafety cap for R-8CBL	₹30,450

C85	RM12C	REMI	RM-12C W/timer Micro Centrifuge w/o Head RPM 16000 with Display	₹39,150
C86	RM1210	REMI	Head RM-1210 Spare Bucket 24 Capillary Capacity 75mmx1mm RPM 12000 Suitable for RM - 12C	₹7,700
C87	RM1214	REMI	Head RM-1214 Spare Bucket 10 (Angle) Capacity 2.0ml RPM 16000 Suitable for RM - 12C	₹7,700
C88	RM1215	REMI	Head RM-1215 Spare Bucket 24 (Angle) Capacity 1.5ml RPM 11850	₹11,375
C89	RM1213	REMI	Head RM-1213 Spare (Reading Device for RM-1210)	₹8,400
C90	RM1216	REMI	Head RM 1216 Spare Bucket 8 Capacity 5ml	₹16,625
C91	R8K	REMI	R-8K Centrifuge For Milk & Oil W /O Rotor with Display	₹50,750
C92	R86	REMI	Head R-86 Spare Rotor for above suitable for R-8k Capacity 4 x 50 ml Head	₹12,600
C93	R8M PLUS	REMI	R-8M PLUS Centrifuge (Max 6000 RPM) With Brushless Induction Motor without Head	₹68,150
C94	R8CPLUS	REMI	R-8C PLUS Centrifuge (Max 6000 RPM) With Brushless Induction Motor without Head	₹56,550
C95	S2415	REMI	Head S 24 -15 Swing Out Rotor 24 x 15ml (Supplied with Glass Tubes)	₹19,575
C96	S815	REMI	Head S 8 -15 Swing Out Rotor 8 x15ml (Supplied with Glass Tubes)	₹11,900
C97	A2415	REMI	Head A24 - 15 Angle Rotor 24 x 15 ml (Supplied with Plastic Tubes)	₹17,690
C98	AR650	REMI	Head AR 6 - 50 Angle Rotor 6 x 50 ml (Supplied with Plastic Tubes)	₹10,500

CENTRIFUGE REMI NEYA SERIES - REMI



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C99	NEYA16	REMI	NEYA 16 Bench Top Centrifuge, (max. speed 16000 rpm) LCD display of speed and time, imbalance detector, dynamic brake safety lid lock, 10 acceleration / deceleration profiles, 10 programs and rotor identification, supplied without rotors and accessories	₹193,938
C100	NEYA16R	REMI	NEYA 16R Refrigerated Bench Top Centrifuge, (max. speed 16000 rpm) LCD display of speed and time, imbalance detector, dynamic brake, safety lid lock, 10 acceleration / deceleration profiles, 10 programs and rotor identification, supplied without rotors and accessories	₹328,933
C101	S4175	REMI	S 4-175 4 place Swing out Rotor	₹24,070
C102	S696MP	REMI	S 6-96 MP 2 place Swing out Rotor for 6 micro plates of 96 wells	₹51,185

CENTRIFUGE HIGH SPEED & COOLING & REFRIGERATED NORMAL& PLUS SERIES - REMI



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C103	R24	REMI	R-24 Revolutionary High Speed Centrifuge Without Rotor Head And Accessories	₹111,723
C104	R242	REMI	R-242 ANGLE Head Rotor With PP Tubes 6X50ML Suitable For R-24 And C-24	₹35,670

C105	CM12PLUS	REMI	CM-12 Plus REFRIGERATED Compact Bench Top Cooling Micro Centrifuges w/o rotor Heads and accessories (-20'c up to 5 ml)	₹271,513
C106	RM1214M	REMI	RM-1214M 10 x 2 ml Angle Head (with Polycarbonate Tubes)	₹11,638
C107	RM1215M	REMI	RM-1215M 24 x 1.5 ml Angle Head (with Polycarbonate Tubes)	₹15,750
C108	RM1216M	REMI	RM-1216M 8 x 5 ml Angle Head (with Polycarbonate Tubes)	₹16,313
C109	RM1217M	REMI	RM-1217M Angle Head suitable for PCR Strip Tubes - Capacity 4 strips of 8 x 0.2 ml	₹24,070
C110	VS02	REMI	VS-02 Servo controlled Voltage Stabilizer	₹20,300
C111	C23PLUS	REMI	C-23 Plus HIGH Volume Low Speed Table Top Cooling Centrifuge Max 7000 RPM / 7560 RCF,LED Display without rotor Heads	₹411,148
C112	CPR23PLUS	REMI	CPR-23 Plus HIGH Volume Low Speed Table Top Cooling Centrifuge Max 7000 RPM / 7400 RCF, LCD Display W/O rotor Heads	₹571,010
C113	R238M	REMI	R-238M Complete with Polypropylene Tubes (Capacity 12 x 50 ML)	₹46,545
C114	R239M	REMI	R-239M Complete with Polypropylene Tubes (Capacity 36 x 15 ML)	₹52,780
C115	VS03	REMI	VS-03 Servo controlled Voltage Stabilizer	₹23,345
C116	C24PLUS	REMI	C-24 Plus LOW VOLUME HIGH SPEED COOLING Centrifuge MAX 20000 RPM / 37570 RCF, LED Display W/O rotor Heads	₹403,390
C117	CPR24PLUS	REMI	CPR-24 Plus LOW VOLUME HIGH SPEED COOLING Centrifuge, MAX 20000 RPM / 37570 RCF, LCD Display W/O rotor Heads	₹564,775
C118	C30PLUS	REMI	C-30 Plus HIGH Volume High Speed FLOOR Model COOLING Centrifuge MAX 20000 RPM / 37570 RCF, LED Display, W/O rotor	₹496,480

C119	CPR30PLUS	REMI	CPR-30 Plus HIGH Volume High Speed FLOOR Model COOLING Centrifuge MAX 20000 RPM / 37570 RCF LCD Display w/o Rotor	₹605,085
C120	R243M	REMI	R-243M Complete with Polypropylene Tubes (Capacity 8 x 25 ML)	₹38,788
C121	R244M	REMI	R-244M Complete with Polypropylene Tubes (Capacity 12 x 15 ML)	₹38,788
C122	R247M	REMI	R-247M Complete with Polypropylene Tubes (Capacity 8 x 50 ML)	₹46,545
C123	R248M	REMI	R-248M Complete with Polypropylene Tubes (Capacity 24 x 1.5 ML)	₹29,508
C124	VS03	REMI	VS-03 Servo controlled Voltage Stabilizer	₹23,345

YSP- SCIENTIFIC SOLUTION
CYCLOMIXER/VORTEX
MIXER DLAB



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C125	MX-S	DLAB	Cat No 8031102000 MX-S Vortex Mixer (adjustable speed) 2500 rpm 4Dia	₹21,170
C126	MX-F	DLAB	Cat.N.8031101000 MX-F Vortex Mixer (Fixed speed)2500 rpm 4Dia	₹19,575

YSP- SCIENTIFIC SOLUTION
CYCLOMIXER / VORTEX MIXER -
REMI AND REMI PLUS SERIES



C127	CM101	REMI	REMI Model CM101 with single and Multi tube Head	₹11,900
C128	CM101PLUS	REMI	CM-101 Plus Vortex Mixer with brushless motor & suitable for Interchangeable mixing Heads	₹15,660



COLD CABINETS (LABORATORY REFRIGERATOR) - REMI

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C129	CC6PLUSLED	REMI	CC-6 Plus Capacity 200 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02) LED	₹271,585
C130	CC6PLUSLCD	REMI	CC-6 Plus Capacity 200 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02) LCD	₹310,300
C131	CC10PLUSLED	REMI	CC-10 Plus Capacity 280 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02) LED	₹285,505
C132	CC10PLUSLCD	REMI	CC-10 Plus Capacity 280 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02) LCD	₹322,770
C133	CC12PLUSLED	REMI	CC-12 Plus Capacity 370 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02) LED	₹301,020
C134	CC12PLUSLCD	REMI	CC-12 Plus Capacity 370 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02) LCD	₹338,285
C135	VS02	REMI	VS-02 Servo controlled Voltage Stabilizer	₹20,300

The high end REMI Items will be Supplied Against Order With Advance . Delivery 3 - 4 Weeks

CARBON MONOXIDE METER (CO) - LUTRON



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C136	GCO2008	LUTRON	Carbon Monoxide Range 0 - 1000 ppm with Temp - 0 - 50'C Model Lutron GCO 2008	₹22,475

CO2 METER - LUTRON



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C137	GCH2018	LUTRON	CO2 Meter Model GCH-2018 For Measuring Carbon Dioxide. Range 0-3000ppm, With Alarms And Warning. RH 10 to 95%,temp, 0 to 50c and dew point Lutron Taiwan Made	₹40,600
C138	GC2028	LUTRON	CO2 meter GC- 2028 With alarm Setting Lutron Taiwan Made	₹33,350

CO2 MONITOR - TESTO / HM



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C139	TESTO3153	TESTO	Testo 315-3 Ambient CO / CO2 Monitor Parallel and direct CO / CO2 measurement, TUV-tested according to EN 50543, Convenient, easy handling, optical and audible alarms, wireless data transfer	POR



COOKING OIL TESTER - TESTO

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C140	TESTO270	TESTO	Testo 270 Cooking Oil Tester Ergonomic and robust design, with temperature sensor along with TPM, clean alarms via backlit display, Rinsable under running water (IP65), HACCP certified (without calibration)	₹78,083
C141	TM-4000	HM	TM-4000 Cooking oil Thermometer	₹3,894

COOKING OIL TESTER - TESTO YSP- SCIENTIFIC SOLUTION



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C142	CAPILLARY TUBE	TOP TECH	BEO Melting Ponit Capillary 3"Long (Pkt 100)	₹70
C143	TT2	TOP TECH	BEO Melting Ponit Capillary 3"Long (Pkt 250)	₹158
C144	TT3	TOP TECH	BEO Melting Ponit Capillary 4"Long (Pkt 100)	₹105
C145	TT4	TOP TECH	BEO Melting Ponit Capillary 4"Long (Pkt 250)	₹236
C146	TT5	TOP TECH	OES Melting Ponit Capillary 3"Long (Pkt 100)	₹96
C147	TT6	TOP TECH	OES Melting Ponit Capillary 3"Long (Pkt 250)	₹201
C148	TT7	TOP TECH	OES Melting Ponit Capillary 4"Long (Pkt 100)	₹123
C149	TT8	TOP TECH	OES Melting Ponit Capillary 4"Long (Pkt 250)	₹263

CAPILLARY TLC - TOP TECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C150	T9	TOP TECH	TLC Capillary 100mm long (Pkt 250)	₹700

CAPILLARY HEPERNISED AND NON HEPERNISED - TOP TECH



C151	T10	TOP TECH	Hepernised Capillary 3" long (Pkt 100)	₹140
C152	T11	TOP TECH	Non Hepernised Capillary 3" long (Pkt 100)	₹105

CAPILLARY GRADUATED - TOP TECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C153	GRADUATED CAPILLARY TUBE	TOP TECH	Graduated Capillary 5 ul 1 Mark &	₹875
C154	T13	TOP TECH	Graduated Capillary 5 ul 5 Mark	₹875
C155	T14	TOP TECH	Graduated Capillary 10 ul 2 Mark	₹875

COMPARATOR - LOVIBOND



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C156	142020	LOVIBOND	142020 Comparator 2000+ (+2 cells)	₹28,348
C157	234060	LOVIBOND	234060 3/40E DISC CHLORINE DPD -TEST 0,02 - 0,3 MG/L - 0,3 MG/L	₹42,195
C158	511050BT	LOVIBOND	511050BT DPD NO.1 TABLETS BF IN PACK OF 100	₹4,725
C159	511080BT	LOVIBOND	511080BT DPD NO.3 TABLETS IN PACK OF 100	₹6,125
C160	511220BT	LOVIBOND	511220BT DPD NO.4 TABLETS IN PACK OF 100	₹6,125

COMPARATOR & DISCS INDIAN



C161		GENERAL	COMPARATOR with 1 pair of square 13.5 mm cells "INDIAN"	₹8,750
C162		GENERAL	NESSELIZER INDIAN WITH SAMPLE TUBES	₹12,250
C163		GENERAL	Spare colour disc price code A CHLORINE / pH / UNIVERSAL INDIAN each Disk	₹6,125
C164		GENERAL	Spare colour disc price code B - HAZEN INDIAN	₹6,125
C165		GENERAL	Spare colour disc price code C INDIAN	₹6,125
C166		GENERAL	Spare 13.5mm, 10ml moulded cells for comparator with stoppers (One Pair) INDIAN	₹2,188

CHECKIT COMPARATOR - SINGLE PARAMETER TEST KITS - LOVIBOND



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C167	147020	LOVIBOND	147020 TESTKIT CHECKIT COMPARATOR CHLORINE DPD 0 - 4 MG/L	₹15,313
C168	147275	LOVIBOND	147275 TESTKIT CHECKIT COMPARATOR OZONE DPD 0 - 1 MG/L	₹15,008
C169	147100	LOVIBOND	147100 TESTKIT CHECKIT COMPARATOR pH 6,5 - 8,4	₹15,313

COLOUR CARDS OR RAPID WATER TEST DROP TEST - LOVIBOND



C170	415020MK	LOVIBOND	IN415020MK pH 2- 13 pH (2, 4, 5,6 7, 8, 9 10,11, 12, 13) Pk. Sz. 100 Colour card ShuddhCheck	₹1,313
C171	415050MK	LOVIBOND	IN415050MK Turbidity 0-80 NTU (0, 0-10, 10-20, 20-40, 40-60, Pk. Sz. 100 Colour card ShuddhCheck	₹2,713
C172	415030MK	LOVIBOND	IN415030MK FRC 0-3 ppm (0, 0.2, 0.5, 1.0, 1.5, 2.0, 3.0) Pk. Sz. 100 Colour card ShuddhCheck	₹3,850
C173	415010MK	LOVIBOND	IN415010MK Fluoride 0- 3 ppm(0, 0.5, 1.0, 1.5, 2.0, 3.0) Pk. Sz. 100 Colour card ShuddhCheck	₹2,713
C174	415060MK	LOVIBOND	IN415060MK IRON HR 0- 4 ppm(0, 1.0, 1.5, 2.0, 3.0, 4.0) Pk. Sz. 100 Colour card ShuddhCheck	₹1,313

C175	415110MK	LOVIBOND	IN415110MK Silica 0-0.530 ppm as SiO ₂ (0, .021, 0.053,0.106, Pk. Sz. 100 Colour card ShuddhCheck	₹10,500
C176	415070MK	LOVIBOND	IN415070MK Iron LR 0-1.5 ppm(0, 0.2, 0.3, 0.4, 0.5, 1.0, 1.5) Pk. Sz. 100 Colour card ShuddhCheck	₹2,188
C177	415120MK	LOVIBOND	IN415120MK Iron ULR 0-0.5 ppm(0, 0.1, 0.2, 0.3, 0.4, 0.5) Pk. Sz. 100 Colour card ShuddhCheck	₹2,188
C178	415040MK	LOVIBOND	IN415040MK Nitrate Mini Kit 0-154 ppm NO ₃ (0, 22, 44, 66, 88, 110, 132, 154 ppm as NO ₃ 0,5,10,15,20,25,30,35 ppm as NO ₃ -N) Pk. Sz. 100ml Colour card ShuddhCheck	₹8,750
C179	416060MK	LOVIBOND	IN416060MK Nitrate 0-40 as NO ₃ -N & 0-175 ppm as NO ₃ () Pk. Sz. 100 Drop Test ShuddhCheck	₹4,375
C180	416010MK	LOVIBOND	IN416010MK Chloride 0-1000 ppm() Pk. Sz. 100 Drop Test ShuddhCheck	₹4,375
C181	416040MK	LOVIBOND	IN416040MK Total Hardness 0-3000 ppm() Pk. Sz. 100 Drop Test ShuddhCheck	₹2,625
C182	416030MK	LOVIBOND	IN416030MK Total Hardness 0-500 ppm() Pk. Sz. 100 Drop Test ShuddhCheck	₹2,188
C183	416020MK	LOVIBOND	IN416020MK Total Hardness 0-50 ppm() Pk. Sz. 100 Drop Test ShuddhCheck	₹2,625
C184	416050MK	LOVIBOND	IN416050MK Total Alkalinity 0-1000 ppm() Pk. Sz. 100 Drop Test ShuddhCheck	₹2,188

 **YSP- SCIENTIFIC SOLUTION**
PHOTOMETER - LOVIBOND



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C185	276000	LOVIBOND	276000 MD100 CHLORINE TABLETS 0,01-6,0 MG/L - PHOTOMETER IN CASE	₹128,470
C186	276160	LOVIBOND	276160 MD100 HAZEN PHOTOMETER IN CASE	₹144,420
C187	276150	LOVIBOND	276150 MD100 SUSPENDED SOLIDS PHOTOMETER IN CASE	₹118,900
C188	214025	LOVIBOND	214025 MD610 PHOTOMETER BLUETOOTH IN CARRYING CASE WITH STANDARD ACCESSORIES	₹271,875

 **YSP- SCIENTIFIC SOLUTION**
COATING THICKNESS GAUGE -
WORKZONE



C189	DFT-222	WORKZONE	DFT-222 Professional high accuracy Coating Thickness gauge, measuring range 0 to 2000um, Resolution 0-99.9um	₹13,125
C190	456	WORKZONE	456 High accuracy Coating Thickness gauge with separated probe, measuring range 3000um, Resolution 0-99.9um	₹23,200
C191		WORKZONE	Coat-PRO Coating Thickness gauge, range 0-1500um/0-60mil, Resolution 0.1um/, temp -10°C to 50°C	₹20,445

 **YSP- SCIENTIFIC SOLUTION**
COATING THICKNESS GAUGE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C192		GENERAL	Coating Thickness gauge AA C1 Non Ferrous	₹26,100

COATMETER / COATING THICKNESS GAUGE - METRAVI



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C193	CTG-01	METRAVI	Model CTG-01 Measures coating thickness on both ferrous and non-ferrous materials	₹21,025
C194	CTG-02	METRAVI	Model CTG-02 Measures coating thickness on both ferrous and non-ferrous materials with extended flexible cable with sensor	₹39,150

COUNTING CHAMBER



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C195		GENERAL	Roheem Silverlite	₹1,925
C196		GENERAL	China Chamber plain without silver coating	₹744

COD MEASUREMENT SYSTEM - LOVIBOND



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C197	2892602	LOVIBOND	COD Measurement system COD set up MD 200 (2892502)COD VARIO COD photometer, COD reactor RD 125Code (2418940), COD tube 0 - 150 mg/l and 0- 1500 mg/l cuvette rack, manual with Combined code 2892602	₹315,375
C198	COD Reagent	LOVIBOND	Lovibond COD Vario tubes 0-150 mg/l or 0- 1500 mg/l /15000 mg/l (Reagents) 2420710 / 2420711/2420712 each pack of 25 each	₹14,350

COD DIGESTER - YSP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
C199	YSPCOD15	YSP	YSP-COD-15 COD Digester with digital temp controller & timer along with glass parts Indian make 15 Samples	₹50,750
C200	YSPCOD6	YSP	YSP-COD-6 COD Digester with digital temp controller & timer along with glass parts Indian make 6 Samples	₹39,875

COMBUSTION BOAT & TUBE



C201		GENERAL	Kumar Make Boat Size 10 X 12 X 88 (Pkt of 150) Price per pcs	₹16
C202		GENERAL	Tube Kumar Make Size 19 X 24 X 600	₹963
C203		GENERAL	Tube Kumar Make Size 17 X 22 X 600	₹875
C204		GENERAL	Tube Kumar Make Size 22 X 29 X 600	₹1,313

D.O METER PEN TYPE - LUTRON



D1	PDO519	LUTRON	Model PDO 519 Lutron DO meter Pen type 0-20 ppm, Taiwan with built in probe	₹20,663
D2	PDO520	LUTRON	Model PDO 520 Lutron DO meter Pen type 0-20 ppm, Taiwan with Detachable probe	₹23,925

D.O METER PORTABLE MODEL & SPARES



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
D3	5509	LUTRON	Model 5509 Lutron DO Meter, 0 to 20mg/L Calibration	₹18,125
D4	5510	LUTRON	Model 5510 Lutron DO Meter 0 to 20mg/L measure oxygen in air	₹28,275
D5	5519	LUTRON	Model 5519 Lutron DO meter 0 - 20 mg/l / Oxygen in Air - 0 - 100 %	₹27,550
D6	OXEL03	LUTRON	Spare Electrolyte OXEL-03 for Lutron DO meters 50 ml bottle	₹2,144
D7	DIAPHRAMM EMBRANE	LUTRON	Do Spare Membrane LUTRON per pcs.	₹1,575

DO METER PORTABLE - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
D8	LH28	LABOHOLIC	LH-28 Portable Digital Dissolved Oxygen Meter 3 1/3 digit LCD display, battery operated, ABS Plastic body. (Against order)	₹8,575

D.O METER BENCH MODEL - EUTECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
D9	ECDO70042S	EUTECH	Model ECDO70042S DO/Temp RANGE 0 to 30 mg/l, NO Memory	POR

D10	ECDO270042	EUTECH	Model ECDO270042 DO/BOD/Temp. RANGE o to 50 mg/l Memory 500 Output RS232	POR
D11	DO2700AN270	EUTECH	ANTECH AN 270 Printer Suitable for ECDO270042S model .	₹36,000

 **YSP- SCIENTIFIC SOLUTION**
DISSOLVED OXYGEN METER
PORTABLE MODEL - EUTECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
D12	DO602K	EUTECH	Model DO602K DO Range 0 to 19.9 mg/l;0 TO199.9% IP Rating IP 54	₹109,620

 **YSP- SCIENTIFIC SOLUTION**

DISSOLVED OXYGEN METER TABLE MODEL -
LABOHOLIC (Against order - No Guarantee)



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
D13	LH52	LABOHOLIC	LH-52 Microprocessor D.O. Meter with 16x 2 alpha numeric LCD display, simultaneous display of O.O. & temp. Data Storage facility upto 1000 samples, computer interlace through USB, Supplied with O.O. probe, temp. probe, electrode stand, dust cover & instructions manual.	₹19,459
D14	LH18	LABOHOLIC	LH-18 Digital Dissolved Oxygen Meter 3.5 Digit LED display with Gold/Silver O.O.Probe	₹10,675
D15	LH27	LABOHOLIC	LH-27 Digital Dissolved Oxygen Analyser cum temperature Meter along with Gold/Silver O.O. probe and RTD temp.probe	₹10,675

DIFFRACTION GRATING HILGER



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
D16		GENERAL	HILGER 15000 Lines	₹3,500
D17		GENERAL	Diffraction Grating Student	₹525

DRY BATH - DLAB



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
D18	HB60-S	DLAB	Cat.no.5032141212 HB60-S LED digital dry bath, with fixed block 5mlx30 +10mlx 8 pcs	₹38,280

DRY BATH INCUBATOR - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
D19	LH42	LABOHOLIC	LH-42 Digital Dry Bath Incubator 24 Holes for 12 mm tubes, +5°C above ambient to70°C. Temperature sensor PT100, Maximum heating power 125W, Temperature Controller PID, Timer 1 min to 99 minutes, LCD Display for time & temperature, Anodized Aluminum block. 220VAC, Available with Buzzer alarm & Overheat protection.	₹10,675



DEEP FREEZER(VERTICAL) REMI PLUS SERIES (Against Order with Advance only) - REMI

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
D20	RQV200PLUSLED	REMI	RQV-200 Plus Capacity 200 liters, Temperature up to -20°C, Microprocessor (VS 02) LED	₹266,945
D21	RQV200PLUSLCD	REMI	RQV-200 Plus Capacity 200 liters, Temperature up to -20°C, Microprocessor (VS 02) LCD	₹307,255
D22	RQV300PLUSLED	REMI	RQV-300 Plus Capacity 300 liters, Temperature up to -20°C, Microprocessor (VS 02) LED	₹307,255
D23	RQV300PLUSLCD	REMI	RQV-300 Plus Capacity 300 liters, Temperature up to -20°C, Microprocessor (VS 02) LCD	₹344,448
D24	RQV400PLUSLED	REMI	RQV-400 Plus Capacity 400 liters, Temperature up to -20°C, Microprocessor (VS 03) LED	₹322,770
D25	RQV400PLUSLCD	REMI	RQV-400 Plus Capacity 400 liters, Temperature up to -20°C, Microprocessor (VS 03) LCD	₹363,080
D26	VS03	REMI	VS-03 Servo controlled Voltage Stabilizer	₹23,345
D27	VS02	REMI	VS-02 Servo controlled Voltage Stabilizer	₹20,155



**DEEP FREEZER(VERTICAL) REMI PLUS SERIES
(Against Order with Advance only) - REMI
YSP- SCIENTIFIC SOLUTION**


SL	CODE	BRAND NAME	PRODUCT NAME	RATE
D28	RFV245	REMI	RFV-245, 245 Liters capacity, Plastic inner chamber & outer body of G.I, temperature range -16°C to -24°C	₹130,355
D29	RFV340	REMI	RFV-340, 340 Liters capacity, Plastic inner chamber & outer body of G.I, temperature range -16°C to -24°C	₹207,930



**ULTRA DEEP FREEZERS (VERTICAL)
(Against Order with Advance) - REMI**


SL	CODE	BRAND NAME	PRODUCT NAME	RATE
D30	ULT90	REMI	ULT-90 Capacity 90 liters, Temp low upto -86c, microprocessor control lcd version(VS 05)	D30
D31	ULT185	REMI	ULT-185 Capacity 185 liters, Temp low upto -86c, microprocessor control lcd version(VS 06)	D31
D32	ULT360	REMI	ULT-360 Capacity 360 liters, Temp low upto -86c, microprocessor control lcd version(VS 07)	₹899,870
D33	VS05	REMI	VS-05 Servo controlled Voltage Stabilizer	₹57,420
D34	VS06	REMI	VS-06 Servo controlled Voltage Stabilizer	₹86,420
D35	VS07	REMI	VS-07 Servo controlled Voltage Stabilizer	₹100,848

ESRITE TUBE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
E1	ESR TUBE	HBG	ESR Tube HBG Make	₹219
E2		TOP TECH	ESR Tube Top Tech	₹105
E3		GENERAL	ESR Tube Superior Make	₹131
E4		GENERAL	ESR Stand 6 Tube powder coated	₹350

ESRITE KIT - TOP TECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
E5	ESRITE KIT E-5	TOP TECH	ESRITE KIT E-5	₹1,313
E6	ESR PIPETTE	TOP TECH	ESR Pipette Plain for ESRITE Kit	₹105
E7	T24	TOP TECH	ESR Pipette Graduated for ESRITE Kit	₹131

ELECTRODES - EUTECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
E8	ECFC7252101B	EUTECH	ECFC7252101B General Purpose Plastic-body Single Junction Gel-Filled pH Combination 12x90mm, BNC 1m for PH700/Tutor etc	₹10,938
E9	ECFC72521R01B	EUTECH	ECFC72521R01B General Plastic-body Single Junction Refillable pH Combination Electrode 12x90mm, BNC, 1m & 10ml refilling electrolyte	₹16,965
E10	ECFC7252201B	EUTECH	ECFC7252201B General Purpose Plastic-body Double Junction Gel-Filled pH Combination 12x90mm,BNC,1m	₹17,325
E11	ECFG7350401B	EUTECH	ECFG7350401B General Purpose Glass-body Single Junction Refillable pH Combination Electrode, annular ceramic reference junction with protective sensor guard, 12x110mm, BNC connector, 1m cable length & 10ml refilling electrolyte	₹17,483
E12	ECFC7352901B	EUTECH	ECFC7352901B General Purpose Plastic-body "3-in-1" pH/Temp. Combination Electrode, 12x90mm, BNC , 1m for Tutor DS	₹16,965
E13	ECFG7370101B	EUTECH	ECFG7370101B Glass body double junction Ag/AgCl refillable pH electrode, 12 x 110 mm, BNC connector, 1m cable (for pH 2700, PC 2700, Ion 2700)	₹43,065
E14	ECFG4390401B	EUTECH	ECFG4390401B Glass-body pH Combination Micro-Electrode, reference Calomel Refillable ,4x55 mm, BNC connector, 1m cable length. Ideal for Micro sample, such as semi-micro cuvettes and NMR tube, Research measurements, fruit juices, beer, milk and yogurt, Tris buffers, clinical, and biological media containing proteins, creams, fats and cosmetics	₹48,285

E15	ECFG7252001B	EUTECH	ECFG7252001B General Purpose, Flat-surface, Glass-body Single Junction Refillable pH Combination Electrode, 12x110mm, BNC connector, 1m cable length. Used for surface measurements such as paper, skin, textile, leather, agar	₹48,285
E16	ECPH5TEMB01P	EUTECH	ECPH5TEMB01P ATC Probe (for Tutor, Ion 510, pH 510, pH 700, Ion 700, pH 1100, pH 2100, pH 2700, & Ion 2700), 1m cable length	₹14,963



YSP- SCIENTIFIC SOLUTION

ELECTRODES - EUTECH

YSP- SCIENTIFIC SOLUTION



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
E17	Lutron PE-03	LUTRON	Lutron PE-03 pH combined electrode Gel filled BNC jack with protection cover	₹3,500
E18	PC22R1	YSP	pH glass combined electrode KCL filling type BNC jack PC22R1	₹875
E19	PC22R2	YSP	pH glass combined electrode BNC with 4 mm banana pin PC22R2	₹919
E20	GCEE1	YSP	pH unbreakable Gel filled combined electrode EPOXY body GCEE1	₹1,260
E21	GCEE2	YSP	pH unbreakable Gel filled combined electrode EPOXY body with 4mm Banana pin GCEE2	₹1,313
E22	CD06R1	YSP	TDS/COND Cell at Cell constant K=1.00(appx.) with BNC Connector CD06R1	₹2,494
E23	CD06R2	YSP	TDS/COND Cell at Cell constant K=1.00(appx.) with 4mm Banana Pin CD06R2	₹2,293
E24	CD06U1	YSP	TDS/COND Cell at Cell constant K=1.00(appx.) Unbreakable with BNC Connector CD06U1	₹2,844

E25	CD06U2	YSP	TDS/COND Cell at Cell constant K=1.00(appx.) Unbreakable with 4mm Banana pin CD06U2	₹2,713
E26	CD04R1	YSP	TDS/COND Cell at Cell constant K=0.50(appx.) with BNC connector CD04R1	₹4,095
E27	CD04R2	YSP	TDS/COND Cell at Cell constant K=0.50(appx.) with 4mm Banana pin CD04R2	₹3,850
E28	CD08R1	YSP	TDS/COND Cell at Cell constant K=0.10(appx.) with BNC connector CD08R1	₹6,388
E29	CD08R2	YSP	TDS/COND Cell at Cell constant K=0.10(appx.) with 4mm Banana pin CD08R2	₹6,300
E30	PE77R1	YSP	Calamel (reference) Electrode with BNC connector for Potentiometer PE77R1	₹805
E31	PE77R2	YSP	Calamel (reference) Electrode with 4mm Banana pin for Potentiometer PE77R2	₹613
E32	PG11R1	YSP	Glass Electrode with BNC Connector for Potentiometer PG11R1	₹840
E33	PM55R1	YSP	Platinum Metallic Electrode with BNC connector for Potentiometer PM55R1	₹1,313
E34	PM55R2	YSP	Platinum Metallic Electrode with 4mm Banana pin for Potentiometer PM55R2	₹1,225
E35	SM66R1	YSP	Silver Metallic Electrode with BNC connector for Potentiometer SM66R1	₹919
E36	SM66R2	YSP	Silver Metallic Electrode with 4mm Banana pin for Potentiometer SM66R2	₹963
E37	ECUR2	YSP	Copper Metallic Electrode with 4mm Banana pin ECUR2	₹700
E38	ECDR2	YSP	Cadmium Metallic Electrode with 4mm Banana pin ECDR2	₹875
E39	EALR2	YSP	Aluminium Metallic Electrode with 4mm Banana pin EALR2	₹569

E40	EZNR2	YSP	Zinc Metallic Electrode with 4mm Banana pin EZNR2	₹744
E41	ECAR2	YSP	Carbon Electrode with 4mm Banana pin ECAR2	₹569
E42	PE88R2	YSP	Reference Electrode (Ag/ AgCL) with 4mm Banana pin PE88R2	₹1,313

 **YSP- SCIENTIFIC SOLUTION**

ELECTRODES - TOSHNIWAL



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
E43	CA11B	GENERAL	Toshniwal Refillable pH combination Electrode Glass Body CA-11 with Banana Pin	₹2,625

 **YSP- SCIENTIFIC SOLUTION**

ELECTRODE STAND - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
E44	SWINGARMST AND	LABOHOLIC	LABOHOLIC SWING ARM STAND WITH GOOD BASE AT BEST RATE	₹1,313



+91 96262 97366
+91 73387 21875



yspscientific@gmail.com
sales.yspscientific@gmail.com



B: 25 G1, Ground Floor, #7/10
Annai Avanthika Vengaivasal,
Chennai - 600100

ELISA READER MICROPLATE - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
E45	LH1260	LABOHOLIC	"LH-1260 Microplate Elisa Reader Designed with long life LED light technology, 7Touch TFT LCD Screen, Test results simultaneous displayed in screen, Bichromatic, 8-channel optical system, 96-well plate, Light source Monochromatic LED, 3 shaking mode with variable time, Large storage capacity up to 200 test protocols and 100,000 results, Self-Check & malfunction alarm, Windows operating system. (Imported)"	₹276,406

ELISA WASHER MICROPLATE - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
E46	LH1261	LABOHOLIC	LH-1261 Microplate Elisa Washer Large LCD Display, Intelligent software, easy to operate, Strong washing function, bottom washing and two positions aspiration make completed washing effect., Plate - 48/96 well and trips (U,V, Flat bottom), Residual volume $\leq 2\mu\text{l}$, Dispensing volume : Adjustable From 50 to 900 μl , 1-16 optional wash cycles, Soak time upto 990 seconds adjustable, Shake time upto 990 sec, 3 channel for washing liquid, 1 channel for waste liquid, Self-calibration and trouble alarm function, Large memory to store up to 100 program, RS232 interface NEW ADDITION	₹276,406

FLAME PHOTOMETER - LABOHOLIC (ALL FLAME WITH 2 FILTER SODIUM (Na) & POTASSIUM (K))



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F1	LH6710	LABOHOLIC	LH-6710 Microprocessor Flame Photometer Touch Screen, Auto Gas Cut-off Facility, Auto Flame Failure Detection Facility, μ P Based with Printer Interface, Five Element Measurement in single Aspiration of the sample, Curve Calibration Programmability using Maximum of 25 Standards Data Processing Through CurveFitting Techniques - Least Squares, Quadratic, Multiple Calibration Curves Saving Facility, Samples Data Storage Facility, Online Help Key, Printer Attachment Facility, Measurement in ppm & meq, Useful for both Medical & Soil Testing Applications, Setup Storage facility upto 10 setups files (with 2 Filters Na & K)NEW ADDITION	₹1,14,985
F2	LH672	LABOHOLIC	LH-672 Microprocessor flame photometer (Graphical Display) Four Element measurement in a single aspiration (Na & Kas standard, Li and Ca optional at extra cost). Calibration curve fit routines first order and second order. Data displayed in concentration units i.e. ppm & meq, automatic filter selection; auto-calibration feature. Printer attachment, PC connecting compatibility, Sample data storage more than 500 results. Alphanumerical keyboard for easy data entry. Graphical LCD display thus easy to read, Curve calibration using up to 5 standard, Multiple calibration curve saving facility, Stainless steel burner. Auto Gas Cut off in case of power failure. (with 2 Filters Na & K)	₹92,873

F3	LH671	LABOHOLIC	<p>LH-671 Microprocessor flame photometer (Graphical Display) Four Element measurement in a single aspiration (Na & K as standard, Li and Ca optional at extra cost). Calibration curve fit routines first order and second order. Data displayed in concentration units i.e. ppm & meq, automatic filter selection; auto-calibration feature. Printer attachment, PC connecting compatibility, Sample data storage more than 500 results. Alphanumeric keyboard for easy data entry. Graphical LCD display thus easy to read, Curve calibration using upto 5 standard, Multiple calibration curve saving facility, Stainless steel burner. Auto Gas Cut off in case of power failure. PROVIDED WITH 2 FILTERS & COMPRESSOR. (with 2 Filters Na & K)</p>	₹77,394
F4	LH66	LABOHOLIC	<p>LH-66 Digital Clinical Flame Photometer Digital Flame Photometer (Dual Channel) single aspirations, with filters (Na & K). Final result in ppm, 2.5 digit dual LED display for simultaneous result along with compressor & accessories, autoignition. Ca & Li filter can be provided on request at extra cost. (with 2 Filters Na & K)</p>	₹43,138
F5	LH65	LABOHOLIC	<p>LH-66 Digital Clinical Flame Photometer Digital Flame Photometer (Dual Channel) single aspirations, with filters (Na & K). Final result in ppm, 2.5 digit dual LED display for simultaneous result along with compressor & accessories, autoignition. Ca & Li filter can be provided on request at extra cost. (with 2 Filters Na & K)</p>	₹37,591
F6	LHCAI	LABOHOLIC	<p>Spare Filter (Calcium or Lithium) for Flame Photometer Against order can be fitted any of the above models.</p>	₹9,625
F7	LHBA	LABOHOLIC	<p>Spare Filter (Barium) for Flame Photometer Against order can be fitted any of the above models</p>	₹17,500

PHOTO FLUOROMETER - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F8	LH330	LABOHOLIC	LH-330 Digital Photo Fluorometer measures fluorophores like vitamins, quinine, steroids and metal. Supplied with 2 primary and 3 secondary filters mounted on holders. Results are displayed on 3 digit 7 segment LED. Sensitivity in 4 ranges	₹71,050

FLOW CUP VISCOMETER



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F9		GENERAL	B-2 Cup with stand and Level	₹5,688
F10		GENERAL	B-3 Cup with stand and Level	₹5,688
F11		GENERAL	B-4 Cup with stand and Level	₹5,688
F12		GENERAL	B-5 Cup with stand and Level	₹5,688
F13		GENERAL	B-6 Cup with stand and Level	₹5,688
F14		GENERAL	B-8 Cup with stand and Level	₹6,563
F15		GENERAL	B-10 Cup with stand and Level	₹5,688
F16		GENERAL	B-2 Cup only	₹4,375
F17		GENERAL	B-3 Cup only	₹4,375

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F18		GENERAL	B-4 Cup only	
F19		GENERAL	B-5 Cup only	
F20		GENERAL	B-6 Cup only	
F21		GENERAL	B-8 Cup only	
F22		GENERAL	B-10 Cup only	
F23		GENERAL	Flow cup as per IS 101 with stand B4 cup only	

 **YSP- SCIENTIFIC SOLUTION**

ZAHN CUP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F24		GENERAL	ZAHN CUP WITH HANDLE NO. 1,2,3,4,5,OR 6 BOTTOM ROUND SHAPE	₹7,000

 **YSP- SCIENTIFIC SOLUTION**

FLOW CUP VISCOMETER



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
SPECIFIC GRAVITY CUP				
F25		GENERAL	S.S. WT./ LITR.CUP 50 CC	₹2,625
F26		GENERAL	S.S.WT./ LITR. CUP 100 CC	₹3,238

FILTER PAPER QUALITATIVE - WHATMAN



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
GRADE 1				
F27	1001042	CYTIVA WHATMAN	CAT NO. 1001 042 GRADE 1 DIAMETER 42.5 PKT OF 100 NO. WHATMAN	₹2,765
F28	1001047	CYTIVA WHATMAN	CAT NO. 1001 047 GRADE 1 DIAMETER 47 PK OF 100 NO. WHATMAN	₹2,905
F29	1001055	CYTIVA WHATMAN	CAT NO. 1001 055 GRADE 1 DIAMETER 55 PK OF 100 NO. WHATMAN	₹3,063
F30	1001070	CYTIVA WHATMAN	CAT NO. 1001 070 GRADE 1 DIAMETER 70 PK OF 100 NO. WHATMAN	₹3,045
F31	1001090	CYTIVA WHATMAN	CAT NO. 1001 090 GRADE 1 DIAMETER 90 PK OF 100 NO. WHATMAN	₹2,555
F32	1001110	CYTIVA WHATMAN	CAT NO. 1001 110 GRADE 1 DIAMETER 110 PK OF 100 NO. WHATMAN	₹2,783
F33	1001125	CYTIVA WHATMAN	CAT NO. 1001125 GRADE 1 DIAMETER 125 PK OF 100 NO. WHATMAN	₹3,185
F34	1001150	CYTIVA WHATMAN	CAT NO. 1001 150 GRADE 1 DIAMETER 150 PK OF 100 NO. WHATMAN	₹3,728
F35	1001185	CYTIVA WHATMAN	Cat No. 1001185 Grade 1 Diameter 185 pk of 100 NO. Whatman	₹5,548
F36	1001240	CYTIVA WHATMAN	Cat No. 1001 240 Grade 1 Diameter 240 pk of 100 NO. Whatman	₹11,988
F37	1001270	CYTIVA WHATMAN	Cat No. 1001 270 Grade 1 Diameter 270 pk of 100 NO. Whatman	₹13,825

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F38	1001320	CYTIVA WHATMAN	Cat No. 1001 320 Grade 1 Diameter 320 pk of 100 NO. Whatman	₹14,910
F39	1001917	CYTIVA WHATMAN	Cat No. 1001 917 Grade 1 Diameter 460x570 Sheets pk of 100 NO. Whatman	₹25,375
	1001917LS	CYTIVA WHATMAN	Cat No. 1001 917 Grade 1 Diameter 460x570 Sheets pk of 100 NO. Whatman LOOSE SHEET - PER SHEET	
F40	1001918	CYTIVA WHATMAN	Cat No. 1001 918 Grade 1 Diameter 460X570 Sheets pk of 500 NO. Whatman	₹92,003

Grade 2

F41	1002055	CYTIVA WHATMAN	Cat No. 1002 055 Grade 2 Diameter 55 pk of 100 NO. Whatman	₹2,940
F42	1002070	CYTIVA WHATMAN	Cat No. 1002 070 Grade 2 Diameter 70 pk of 100 NO. Whatman	₹3,185
F43	1002090	CYTIVA WHATMAN	Cat No. 1002 090 Grade 2 Diameter 90 pk of 100 NO. Whatman	₹3,098
F44	1002110	CYTIVA WHATMAN	Cat No. 1002 110 Grade 2 Diameter 110 pk of 100 NO. Whatman	₹3,483
F45	1002125	CYTIVA WHATMAN	Cat No. 1002 125 Grade 2 Diameter 125 pk of 100 NO. Whatman	₹3,763
F46	1002150	CYTIVA WHATMAN	Cat No. 1002 150 Grade 2 Diameter 150 pk of 100 NO. Whatman	₹5,705
F47	1002185	CYTIVA WHATMAN	Cat No. 1002 185 Grade 2 Diameter 185 pk of 100 NO. Whatman	₹7,823
F48	1002917	CYTIVA WHATMAN	Cat No. 1002 917 Grade 2 Diameter 460x570 Sheets pk of 100 NO. Whatman	₹43,268

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

Grade 3

F49	1003070	CYTIVA WHATMAN	Cat No. 1003 070 Grade 3 Diameter 70 pk of 100 NO. Whatman	₹4,848
F50	1003090	CYTIVA WHATMAN	Cat No.1003 090 Grade 3 Diameter 90 pk of 100 NO. Whatman	₹5,338
F51	1003110	CYTIVA WHATMAN	Cat No. 1003 110 Grade 3 Diameter 110 pk of 100 NO. Whatman	₹4,918
F52	1003125	CYTIVA WHATMAN	Cat No. 1003 125 Grade 3 Diameter 125 pk of 100 NO. Whatman	₹6,930
F53	1003150	CYTIVA WHATMAN	Cat No. 1003 150 Grade 3 Diameter 150 pk of 100 NO. Whatman	₹7,543
F54	1003185	CYTIVA WHATMAN	Cat No. 1003 185 Grade 3 Diameter 185 pk of 100 NO. Whatman	₹14,508
F55	1003917	CYTIVA WHATMAN	Cat No. 1003 917 Grade 3 Diameter 460x570Sheets pk of 100 NO. Whatman	₹77,778

Grade 4

F56	1004042	CYTIVA WHATMAN	Cat No. 1004 042 Grade 4 Diameter 42.5 pk of 100 NO. Whatman	₹2,643
F57	1004047	CYTIVA WHATMAN	Cat No. 1004 047 Grade 4 Diameter 47 pk of 100 NO. Whatman	₹2,643
F58	1004055	CYTIVA WHATMAN	Cat No. 1004 055 Grade 4 Diameter 55 pk of 100 NO.Whatman	₹2,765
F59	1004070	CYTIVA WHATMAN	Cat No. 1004 070 Grade 4 Diameter 70 pk of 100 NO.Whatman	₹2,678

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F60	1004090	CYTIVA WHATMAN	Cat No. 1004 090 Grade 4 Diameter 90 pk of 100 NO. Whatman	₹3,028
F61	1004110	CYTIVA WHATMAN	Cat No. 1004 110 Grade 4 Diameter 110 pk of 100 NO. Whatman	₹3,605
F62	1004125	CYTIVA WHATMAN	Cat No. 1004 125 Grade 4 Diameter 125 pk of 100 NO. Whatman	₹4,270
F63	1004150	CYTIVA WHATMAN	Cat No. 1004 150 Grade 4 Diameter 150 pk of 100 NO. Whatman	₹5,898
F64	1004185	CYTIVA WHATMAN	Cat No. 1004 185 Grade 4 Diameter 185 pk of 100 NO. Whatman	₹8,698
F65	1004240	CYTIVA WHATMAN	Cat No. 1004 240 Grade 4 Diameter 240 pk of 100 NO. Whatman	₹15,365
F66	1004917	CYTIVA WHATMAN	Cat No. 1004 917 Grade 4 Diameter 460x570 Sheets pk of 100 NO. Whatman	₹52,447

Grade 5

F67	1005070	CYTIVA WHATMAN	Cat No. 1005 070 Grade 5 Diameter 70 pk of 100 NO. Whatman	₹3,483
F68	1005090	CYTIVA WHATMAN	Cat No. 1005 090 Grade 5 Diameter 90 pk of 100 NO. Whatman	₹3,763
F69	1005110	CYTIVA WHATMAN	Cat No. 1005 110 Grade 5 Diameter 110 pk of 100 NO. Whatman	₹4,060
F70	1005125	CYTIVA WHATMAN	Cat No. 1005 125 Grade 5 Diameter 125 pk of 100 NO. Whatman	₹5,058
F71	1005150	CYTIVA WHATMAN	Cat No. 1005 150 Grade 5 Diameter 150 pk of 100 NO. Whatman	₹7,998
F72	1005185	CYTIVA WHATMAN	Cat No. 1005 185 Grade 5 Diameter 185 pk of 100 NO. Whatman	₹12,758

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Grade 40				
F73	1440012	CYTIVA WHATMAN	Cat No.1440 012 Grade 40 Diameter 12.7 / 1.27 cm pk of 400 NO.Whatman	₹5,355
F74	1440042	CYTIVA WHATMAN	Cat No.1440 042 Grade 40 Diameter 42.5 k of 100 NO. Whatman	₹6,668
F75	1440047	CYTIVA WHATMAN	Cat No. 1440 047 Grade 40 Diameter 47 pk of 100 NO. Whatman	₹6,965
F76	1440055	CYTIVA WHATMAN	Cat No. 1440 055 Grade 40 Diameter 55 pk of 100 NO. Whatman	₹4,935
F77	1440070	CYTIVA WHATMAN	Cat No. 1440 070 Grade 40 Diameter 70 pk of 100 NO. Whatman	₹5,338
F78	1440090	CYTIVA WHATMAN	Cat No. 1440 090 Grade 40 Diameter 90 pk of 100 NO. Whatman	₹6,755
F79	1440110	CYTIVA WHATMAN	Cat No. 1440 110 Grade 40 Diameter 110 pk of 100 NO. Whatman	₹7,875
F80	1440125	CYTIVA WHATMAN	Cat No. 1440 125 Grade 40 Diameter 125 pk of 100 NO. Whatman	₹9,765
F81	1440150	CYTIVA WHATMAN	Cat No. 1440 150 Grade 40 Diameter 150 pk of 100 NO. Whatman	₹12,740
F82	1440185	CYTIVA WHATMAN	Cat No. 1440 185 Grade 40 Diameter 185 pk of 100 NO. Whatman	₹17,748
F83	1440240	CYTIVA WHATMAN	Cat No. 1440 240 Grade 40 Diameter 240 pk of 100 NO. Whatman	₹36,294
F84	1440917	CYTIVA WHATMAN	Cat No. 1440 917 Grade 40 Diameter 460x570 Sheets 100 NO. Whatman	₹1,13,318

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

Grade 41

F85	1441070	CYTIVA WHATMAN	Cat No. 1441 070 Grade 41 Diameter 70 pk of 100 NO Whatman	₹5,425
F86	1441090	CYTIVA WHATMAN	Cat No. 1441 090 Grade 41 Diameter 90 pk of 100 NO. Whatman	₹6,825
F87	1441110	CYTIVA WHATMAN	Cat No. 1441 110 Grade 41 Diameter 110 pk of 100 NO. Whatman	₹9,083
F88	1441125	CYTIVA WHATMAN	Cat No. 1441 125 Grade 41 Diameter 125 pk of 100 NO. Whatman	₹10,395
F89	1441150	CYTIVA WHATMAN	Cat No. 1441 150 Grade 41 Diameter 150 pk of 100 NO. Whatman	₹13,265
F90	1441185	CYTIVA WHATMAN	Cat No. 1441 185 Grade 41 Diameter 185 pk of 100 NO. Whatman	₹17,299
F91	1441240	CYTIVA WHATMAN	CatNo. 1441 240 Grade 41 Diameter 240 pk of 100 NO. Whatman	₹34,394
F92	1441917	CYTIVA WHATMAN	Cat No. 1441 917 Grade 41 Diameter 460x570 Sheets pk of 100 NO. Whatman	₹1,04,806

Grade 42

F93	1442042	CYTIVA WHATMAN	Cat No. 1442 042 Grade 42 Diameter 42.5 pk of 100 NO. Whatman	₹4,848
F94	1442055	CYTIVA WHATMAN	Cat No. 1442 055 Grade 42 Diameter 55 pk of 100 NO. Whatman	₹4,428
F95	1442070	CYTIVA WHATMAN	Cat No. 1442 070 Grade 42 Diameter 70 pk of 100 NO. Whatman	₹5,425
F96	1442090	CYTIVA WHATMAN	Cat No. 1442 090 Grade 42 Diameter 90 pk of 100 NO. Whatman	₹6,650

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F97	1442110	CYTIVA WHATMAN	Cat No. 1442 110 Grade 42 Diameter 110 pk of 100 NO. Whatman	₹4,640
F98	1442125	CYTIVA WHATMAN	Ca No. 1442 125 Grade 42 Diameter 125 pk of 100 NO. Whatman	₹5,580
F99	1442150	CYTIVA WHATMAN	Cat No. 1442 150 Grade 42 Diameter 150 pk of 100 NO. Whatman	₹7,330
F100	1442185	CYTIVA WHATMAN	Cat No. 1442 185 Grade 42 Diameter 185 pk of 100 NO. Whatman	₹11,820
F101	1442240	CYTIVA WHATMAN	CatNo. 1442 240 Grade 42 Diameter 240 100 NO. Whatman	₹19,010
F102	1442917	CYTIVA WHATMAN	Cat No. 1442 917 Grade 42 Diameter 460x570 Sheets pk of 100 NO. Whatman	₹68,130

Grade 44

F103	1444070	CYTIVA WHATMAN	Cat No. 1444 070 Grade 44 Diameter 70 pk of 100 NO. Whatman	₹3,810
F104	1444090	CYTIVA WHATMAN	Cat No. 1444 090 Grade 44 Diameter 90 pk of 100 NO. Whatman	₹4,480
F105	1444110	CYTIVA WHATMAN	Cat No. 1444 110 Grade 44 Diameter 110 pk of 100 NO. Whatman	₹5,300
F106	1444125	CYTIVA WHATMAN	Cat No. 1444 125 Grade 44 Diameter 125 pk of 100 NO. Whatman	₹5,700
F107	1444150	CYTIVA WHATMAN	Cat No. 1444 150 Grade 44 Diameter 150 pk of 100 NO. Whatman	₹10,090
F108	1444185	CYTIVA WHATMAN	Cat No. 1444 185 Grade 44 Diameter 185 pk of 100 NO. Whatman	₹14,590

WHATMAN - HARDENED FILTERS



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Grade 50				
F109	1450070	CYTIVA WHATMAN	Cat No. 145 070 Grade 50 Diameter 70 pk of 100 NO. Whatman	₹6,668
F110	1450090	CYTIVA WHATMAN	Cat No. 1450 090 Grade 50 Diameter 90 pk of 100 NO. Whatman	₹7,928
F111	1450110	CYTIVA WHATMAN	Cat No. 1450 110 Grade 50 Diameter 110 pk of 100 NO. Whatman	₹12,180
F112	1450125	CYTIVA WHATMAN	CatNo. 1450 125 Grade 50 Diameter 125 pk of 100 NO. Whatman	₹12,845
F113	1450150	CYTIVA WHATMAN	Cat No. 1450 150 Grade 50 Diameter 150 pk of 100 NO. Whatman	₹15,617
F114	1450185	CYTIVA WHATMAN	Cat No. 1450 185 Grade 50 Diameter 185 pk of 100 NO. Whatman	₹23,171
F115	1450917	CYTIVA WHATMAN	Cat No. 1450 917 Grade 50 Diameter 460x570 Sheets pk of 100 NO. Whatman	₹1,23,265
Grade 52				
F116	1452090	CYTIVA WHATMAN	Cat No. 1452 090 Grade 52 Diameter 90 pk of 100 NO. Whatman	₹10,728
F117	1452110	CYTIVA WHATMAN	Cat No. 1452 110 Grade 52 Diameter 110 pk of 100 NO. Whatman	₹13,405
F118	1452125	CYTIVA WHATMAN	Cat No. 1452 125 Grade 52 Diameter 125 pk of 100 NO. Whatman	₹15,313
F119	1452150	CYTIVA WHATMAN	Cat No. 1452 150 Grade 52 Diameter 150 pk of 100 NO. Whatman	₹16,095

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

Grade 54

F120	1454070	CYTIVA WHATMAN	Cat No. 1454 070 Grade 54 Diameter 70 pk of 100 NO. Whatman	₹6,475
F121	1454090	CYTIVA WHATMAN	Cat No. 1454 090 Grade 54 Diameter 90 pk of 100 NO. Whatman	₹8,698
F122	1454110	CYTIVA WHATMAN	Cat No. 1454 110 Grade 54 Diameter 110 pk of 100 NO. Whatman	₹11,078
F123	1454125	CYTIVA WHATMAN	Cat No. 1454 125 Grade 54 Diameter 125 pk of 100 NO. Whatman	₹13,790
F124	1454150	CYTIVA WHATMAN	Cat No. 1454 150 Grade 54 Diameter 150 pk of 100 NO. Whatman	₹16,588
F125	1454185	CYTIVA WHATMAN	Cat No. 1454 185 Grade 54 Diameter 185 pk of 100 NO. Whatman	₹23,795

Grad 540

F126	1540090	CYTIVA WHATMAN	Cat No. 1540 090 Grade 540 Diameter 90 pk of 100 NO. Whatman	₹8,698
F127	1540110	CYTIVA WHATMAN	Cat No. 1540 110 Grade 540 Diameter 110 pk of 100 NO. Whatman	₹11,900
F128	1540125	CYTIVA WHATMAN	Cat No. 1540 125 Grade 540 Diameter 125 pk of 100 NO. Whatman	₹12,933
F129	1540150	CYTIVA WHATMAN	Cat No. 1540 150 Grade 540 Diameter 150 pk of 100 NO. Whatman	₹16,110
F130	1540185	CYTIVA WHATMAN	Cat No. 1540 185 Grade 540 Diameter 185 pk of 100 NO. Whatman	₹27,652

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

Grade 541

F131	1541070	CYTIVA WHATMAN	Cat No. 1541 070 Grade 541 Diameter 70 pk of 100 NO. Whatman	₹6,668
F132	1541090	CYTIVA WHATMAN	Cat No. 1541 090 Grade 541 Diameter 90 pk of 100 NO. Whatman	₹8,120
F133	1541110	CYTIVA WHATMAN	Cat No. 1541 110 Grade 541 Diameter 110 pk of 100 NO. Whatman	₹11,603
F134	1541125	CYTIVA WHATMAN	Cat No. 1541 125 Grade 541 Diameter 125 pk of 100 NO. Whatman	₹11,358
F135	1541150	CYTIVA WHATMAN	Cat No. 1541 150 Grade 541 Diameter 150 pk of 100 NO. Whatman	₹16,197
F136	1541185	CYTIVA WHATMAN	Cat No. 1541 185 Grade 541 Diameter 185 pk of 100 NO. Whatman	₹26,202

Grade 542

F137	1542070	CYTIVA WHATMAN	Cat No. 1542 070 Grade 542 Diameter 70 pk of 100 NO. Whatman	₹6,668
F138	1542090	CYTIVA WHATMAN	Cat No. 1542 090 Grade 542 Diameter 90 pk of 100 NO. Whatman	₹9,573
F139	1542110	CYTIVA WHATMAN	Cat No. 1542 110 Grade 542 Diameter 110 pk of 100 NO. Whatman	₹13,055
F140	1542125	CYTIVA WHATMAN	Cat No. 1542 125 Grade 542 Diameter 125 pk of 100 NO. Whatman	₹13,790
F141	1542150	CYTIVA WHATMAN	Cat No. 1542 150 Grade 542 Diameter 150 pk of 100 NO. Whatman	₹17,545
F142	1542185	CYTIVA WHATMAN	Cat No. 1542 185 Grade 542 Diameter 185 pk of 100 NO. Whatman	₹26,202

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

Grade 2V - Folded Prepleated Grades

F143	1202125	CYTIVA WHATMAN	Cat NO. 1202-125 Grade GR 2 FF Diameter 125mm pk of 100 NO Whatman	₹14,123
------	---------	-------------------	-----------------------------------------------------------------------	---------

EPM 2000

F144	1882047	CYTIVA WHATMAN	Cat No. 1882 047 Grade EPM 2000 Diameter 47 pk of 100 NO. Whatman	₹15,855
------	---------	-------------------	----------------------------------------------------------------------	---------

F145	1882866	CYTIVA WHATMAN	Cat No. 1882 866 Grade EPM 2000 Diameter 203 X 254* Sheets pk of 100 NO. Whatman	₹1,04,864
------	---------	-------------------	----------------------------------------------------------------------------------------	-----------

934 - AH

F146	1827047	CYTIVA WHATMAN	Cat No.1827 047 Grade 934-AH Diameter 47 pk of 100 NO. Whatman	₹16,765
------	---------	-------------------	-------------------------------------------------------------------	---------

F147	1827055	CYTIVA WHATMAN	Cat No. 1827 055 Grade 934-AH Diameter 55 pk of 100 NO. Whatman	₹16,023
------	---------	-------------------	--------------------------------------------------------------------	---------

F148	1827070	CYTIVA WHATMAN	Cat No. 1827 070 Grade 934-AH Diameter 70 pk of 100 NO. Whatman	₹20,141
------	---------	-------------------	--------------------------------------------------------------------	---------

F149	1827090	CYTIVA WHATMAN	Cat No. 1827 090 Grade 934-AH Diameter 90 pk of 100 NO. Whatman	₹30,102
------	---------	-------------------	--------------------------------------------------------------------	---------

F150	1827110	CYTIVA WHATMAN	Cat No. 1827 110 Grade 934-AH Diameter 110 pk of 100 NO. Whatman	₹36,787
------	---------	-------------------	---------------------------------------------------------------------	---------



Grade 2V - Folded Prepleated Filters offer consistent performance and easy handling. Ideal for precise filtration needs, they're a dependable choice for any lab setup.



WHATMAN - CELLULOSE THIMBLES



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
CELLULOSE EXTRACTION THIMBLES SINGLE				
F151	2800199	CYTIVA WHATMAN	Cat No. 2800199 CELLULOSE EXTRACTION THIMBLES Diameter 19X90 pk of 25 NO. Whatman	₹20,677
F152	2800228	CYTIVA WHATMAN	Cat No. 2800228 CELLULOSE EXTRACTION THIMBLES Diameter 22X80 pk of 25 NO. Whatman	₹20,909
F153	2800258	CYTIVA WHATMAN	Cat No. 2800258 CELLULOSE EXTRACTION THIMBLES Diameter 25X80 pk of 25 NO. Whatman	₹20,822
F154	2800250	CYTIVA WHATMAN	Cat No. 2800250 CELLULOSE EXTRACTION THIMBLES Diameter 25X100 pk of 25 NO. Whatman	₹20,677
F155	2800288	CYTIVA WHATMAN	Cat No. 2800288 CELLULOSE EXTRACTION THIMBLES Diameter 28X80 pk of 25 NO. Whatman	₹21,011
F156	2800308	CYTIVA WHATMAN	Cat No. 2800308 CELLULOSE EXTRACTION THIMBLES Diameter 30X80 pk of 25 NO. Whatman	₹21,011
F157	2800330	CYTIVA WHATMAN	Cat No. 2800330 CELLULOSE EXTRACTION THIMBLES Diameter 33X100 pk of 25 NO. Whatman	₹25,114
F158	2800432	CYTIVA WHATMAN	Cat No. 2800432 CELLULOSE EXTRACTION THIMBLES Diameter 43X12z pk of 25 NO. Whatman	₹31,277
F159	2800280	CYTIVA WHATMAN	Cat No. 2800280 CELLULOSE EXTRACTION THIMBLES Diameter 28X100 pk of 25 NO. Whatman	₹20,677

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F160	2800300	CYTIVA WHATMAN	CAT NO. 2800300 CELLULOSE EXTRACTION THIMBLES DIAMETER 30X100 PK OF 25 NO. WHATMAN	₹21,141

ASHLESS CLIPPINGS

F161	1703050	CYTIVA WHATMAN	Cat No.1703 050 ASHLESS CLIPPINGS 500g	₹13,808
------	---------	-------------------	-------------------------------------------	----------------

PARAFILM 'M' USA MAKE SPECIAL PRICE

F162		GENERAL	PARAFILM 4X 125 feet	₹4,069
F163		GENERAL	PARAFILM 2X 250 feet	₹4,069

NO RETURN / REPLACEMENT for Whatman. Kindly Order with the Correct order codes.

 **YSP- SCIENTIFIC SOLUTION**

WHATMAN - GLASS MICROFIBER FILTERS



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
GF / A				
F164	1820021	CYTIVA WHATMAN	Cat No. 1820 021 Grade GF/A Diameter 21 pk of 100 NO. Whatman	₹6,143
F165	1820024	CYTIVA WHATMAN	Cat No. 1820 024 Grade GF/A Diameter 24 pk of 100 NO. Whatman	₹6,738
F166	1820025	CYTIVA WHATMAN	Cat No. 1820 025 Grade GF/A Diameter 25 pk of 100 NO. Whatman	₹6,528

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F167	1820037	CYTIVA WHATMAN	CAT NO. 1820 037 GRADE GF/A DIAMETER 37 PK OF 100 NO. WHATMAN	₹8,960
F168	1820042	CYTIVA WHATMAN	CAT NO.1820 042 GRADE GF/A DIAMETER 42.5 PK OF 100 NO.WHATMAN	₹9,678
F169	1820047	CYTIVA WHATMAN	CAT NO.1820 047 GRADE GF/A DIAMETER 47 PK OF 100 NO.WHATMAN	₹9,503
F170	1820055	CYTIVA WHATMAN	CAT NO. 1820 055 GRADE GF/A DIAMETER 55 PK OF 100 NO.WHATMAN	₹13,178
F171	1820070	CYTIVA WHATMAN	CAT NO. 1820 070 GRADE GF/A DIAMETER 70 PK OF 100 NO.WHATMAN	₹14,123
F172	1820090	CYTIVA WHATMAN	CAT NO. 1820 090 GRADE GF/A DIAMETER 90 PK OF 100 NO.WHATMAN	₹15,182
F173	1820110	CYTIVA WHATMAN	CAT NO. 1820 110 GRADE GF/A DIAMETER 110 PK OF 100 NO.WHATMAN	₹18,198
F174	1820125	CYTIVA WHATMAN	CAT NO. 1820 125 GRADE GF/A DIAMETER 125 PK OF 100 NO. WHATMAN	₹25,317
F175	1820150	CYTIVA WHATMAN	CAT NO. 1820 150 GRADE GF/A DIAMETER 150 PK OF 100 NO.WHATMAN	₹34,235
F176	1820866	CYTIVA WHATMAN	CAT NO. 1820 866 GRADE GF/A DIAMETER 203 X 254* SHEETS 100 NO.PK SIZE IN MM WHATMAN	₹52,577

GF/ B

F177	1821025	CYTIVA WHATMAN	CAT NO. 1821 025 GRADE GF/B DIAMETER 25 PK OF 100 NO.WHATMAN	₹8,418
------	---------	-------------------	--------------------------------------------------------------------	--------

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F178	1821042	CYTIVA WHATMAN	CAT NO. 1821 042 GRADE GF/B DIAMETER 42.5 PK OF 100 NO.WHATMAN	₹17,226
F179	1821047	CYTIVA WHATMAN	CAT NO. 1821 047 GRADE GF/B DIAMETER 47 PK OF 100 NO. WHATMAN	₹15,278
F180	1821055	CYTIVA WHATMAN	CAT NO. 1821 055 GRADE GF/B DIAMETER 55 PK OF 100 NO.WHATMAN	₹19,662
F181	1821070	CYTIVA WHATMAN	CAT NO. 1821 070 GRADE GF/B DIAMETER 70 PK OF 100 NO.WHATMAN	₹23,055
F182	1821090	CYTIVA WHATMAN	CAT NO. 1821 090 GRADE GF/B DIAMETER 90 PK OF 100 NO. WHATMAN	₹13,370
F183	1821110	CYTIVA WHATMAN	CAT NO. 1821 110 GRADE GF/B DIAMETER 110 PK OF 100 NO.WHATMAN	₹16,643
F184	1821125	CYTIVA WHATMAN	CAT NO. 1821 125 GRADE GF/B DIAMETER 125 PK OF 100 NO.WHATMAN	₹18,198
F185	1821150	CYTIVA WHATMAN	CAT NO. 1821 150 GRADE GF/B DIAMETER 150 PK OF 100 NO.WHATMAN	₹23,403

GF / C

F186	1822025	CYTIVA WHATMAN	CAT NO. 1822 025 GRADE GF/C DIAMETER 25 PK OF 100 NO. WHATMAN	₹6,440
F187	1822037	CYTIVA WHATMAN	CAT NO. 1822 037 GRADE GF/C DIAMETER 37 PK OF 100 NO. WHATMAN	₹9,870
F188	1822042	CYTIVA WHATMAN	CATNO. 1822 042 GRADE GF/C DIAMETER 42.5 PK OF 100 NO.WHATMAN	₹10,395

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F189	1822047	CYTIVA WHATMAN	CAT NO. 1822 047 GRADE GF/C DIAMETER 47 PK OF 100 NO. WHATMAN	₹9,520
F190	1822055	CYTIVA WHATMAN	CAT NO. 1822 055 GRADE GF/C DIAMETER 55 PK OF 100 NO.WHATMAN	₹15,418
F191	1822070	CYTIVA WHATMAN	CAT NO. 1822 070 GRADE GF/C DIAMETER 70 PK OF 100 NO.WHATMAN	₹16,835
F192	1822090	CYTIVA WHATMAN	CAT NO. 1822 090 GRADE GF/C DIAMETER 90 PK OF 100 NO. WHATMAN	₹17,357
F193	1822110	CYTIVA WHATMAN	CAT NO. 1822 110 GRADE GF/C DIAMETER 110 PK OF 100 NO. WHATMAN	₹23,302
F194	1822125	CYTIVA WHATMAN	CAT NO. 1822 125 GRADE GF/C DIAMETER 125 PK OF 100 NO. WHATMAN	₹28,478
F195	1822150	CYTIVA WHATMAN	CAT NO. 1822 150 GRADE GF/C DIAMETER 150 PK OF 100 NO. WHATMAN	₹39,397

GF/ D

F196	1823025	CYTIVA WHATMAN	CAT NO. 1823-025 GRADE GF/D DIAMETER 25 PK OF 100 NO HATMAN	₹8,418
F197	1823047	CYTIVA WHATMAN	CAT NO. 1823-047 GRADE GF/D DIAMETER 47 PK OF 100 NO WHATMAN	₹16,590

GF/F

F198	1825025	CYTIVA WHATMAN	CAT NO1825-025 GF/F GLASS CIRCLES, 25MM 100/PK WHATMAN	₹8,418
F199	1825047	CYTIVA WHATMAN	CAT NO 1825-047 GF/F GLASS CIRCLES, 47MM 100/PK WHATMAN	₹16,590

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
GLASS MICROFIBRE THIMBLES				
F200	2814199	CYTIVA WHATMAN	Cat No.2814199 Diameter 19X90 pk of 25 NO. Whatman	₹36,874

 **YSP- SCIENTIFIC SOLUTION**

WHATMAN - CHROMATOGRAGHY



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Pure Cellulose Sheets & Reels				
1 Chr (0.18mm) thickness				
F201	3001845	CYTIVA WHATMAN	Cat No. 3001845 Size(mm) 100X300 Micro granular 100 Sheets Whatman	₹13,983
F202	3001861	CYTIVA WHATMAN	Cat No. 3001861 Size(mm) 200X200 Micro granular 100 Sheets Whatman	₹17,238
F203	3001878	CYTIVA WHATMAN	Cat No. 3001878 Size(mm) 250x250 Micro granular 100 Sheets Whatman	₹22,374
F204	3001917	CYTIVA WHATMAN	Cat No. 3001917 Size(mm) 460x570 Micro granular 100 Sheets Whatman	₹48,735
F205	3001931	CYTIVA WHATMAN	Cat No. 3001931 Size(mm) 580x680 Micro granular 100 Sheets Whatman	₹75,139
F206	3001604	CYTIVA WHATMAN	Cat No. 3001604 Size(mm) 10X100m Micro granular Reels Whatman	₹6,335
F207	3001640	CYTIVA WHATMAN	Cat No. 3001640 Size(mm) 30X100m Micro granular Reels Whatman	₹11,498

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

2 Chr (0.18 mm Thickness)

F208	3002917	CYTIVA WHATMAN	Cat No.3002917 Size(mm) 460x570 Micro granular 100 Sheets Whatman	₹70,238
------	---------	-------------------	----------------------------------------------------------------------	---------

3 Chr (0.34 mm Thickness)

F209	3003917	CYTIVA WHATMAN	Cat No.3003917 Size(mm) 460x570 Micro granular 100 Sheets Whatman	₹1,18,625
------	---------	-------------------	----------------------------------------------------------------------	-----------

 **YSP- SCIENTIFIC SOLUTION**
WHATMAN - NYLON
MEMBRANE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F210	7402001	CYTIVA WHATMAN	Cat No 7402001 NYLON MEMBRANE 13mm 0.2um pk 100 Whatman	₹20,228
F211	7402002	CYTIVA WHATMAN	Cat No 7402002 NYLON MEMBRANE 25 mm 0.2um pk 100 Whatman	₹20,286
F212	7402004	CYTIVA WHATMAN	Cat No 7402004 NYLON MEMBRANE 47 mm 0.2um pk 100 Whatman	₹15,182
F213	7404001	CYTIVA WHATMAN	Cat No 7404001 NYLON MEMBRANE 13 mm 0.45um pk 100 Whatman	₹20,228
F214	7404002	CYTIVA WHATMAN	Cat No 7404002 NYLON MEMBRANE 25 mm 0.45um pk 100 Whatman	₹22,852
F215	7404004	CYTIVA WHATMAN	Cat No 7404004 NYLON MEMBRANE 47 mm 0.45um pk 100 Whatman	₹19,271

WHATMAN - CELLULOSE NITRATE MEMBRANE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F216	7182001	CYTIVA WHATMAN	Cat No 7182001 CELLULOSE NITRATE MEMBRANE 13mm 0.2um pk 100 Whatman	₹21,663
F217	7182002	CYTIVA WHATMAN	Cat No 7182002 CELLULOSE NITRATE MEMBRANE 25 mm 0.2um pk 100 Whatman	₹15,488
F218	7182004	CYTIVA WHATMAN	Cat No 7182004 CELLULOSE NITRATE MEMBRANE 47 mm 0.2um pk 100 Whatman	₹11,183
F219	7184001	CYTIVA WHATMAN	Cat No 7184001 CELLULOSE NITRATE MEMBRANE 13 mm 0.45um pk 100 Whatman	₹18,082
F220	7184002	CYTIVA WHATMAN	Cat No 7184002 CELLULOSE NITRATE MEMBRANE 25 mm 0.45um pk 100 Whatman	₹18,082
F221	7184004	CYTIVA WHATMAN	Cat No 7184004 CELLULOSE NITRATE MEMBRANE 47 mm 0.45um pk 100 Whatman	₹10,640

We're your go-to source for all laboratory products.



Probe Sonicator



ice maker



ph meter

Upgrade your lab with top-quality equipment! Get the Probe Sonicator for efficient sample preparation, an Ice Maker for reliable cooling, and a PH Meter for accurate readings. Everything you need for precise lab work!

WHATMAN - PM2.5 / LENS CLEANNING TISSUE / QM A / BENCHKOTE




SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F222	1851025	CYTIVA WHATMAN	Cat No 1851-025 QM-A Quartz Circles,25mm 100/pk Whatman	₹31,030
F223	1851037	CYTIVA WHATMAN	Cat No 1851-037 QM-A Circles, 37mm 100/pk Whatman	₹39,165
F224	1851047	CYTIVA WHATMAN	Cat No 1851-047 QM-A Circles, 47mm 100/pk Whatman	₹40,948
F225	7592104	CYTIVA WHATMAN	Cat No 7592-104 PM2.5 Membrane, PTFE 46.2mm with support ring, sequentially numbered 50/pk Whatman	₹32,683
F226	2017006	CYTIVA WHATMAN	Cat No 2017-006 Antibiotic Assay Discs,6mm 1000/pk Whatman	₹19,575
F227	2105841	CYTIVA WHATMAN	Cat No 2105-841 105 Sheets, Lens Cleaning Tissue, 100 x 150mm, 25 wallets of 25 sheets 25/pk Whatman	₹16,695
F228	2105862	CYTIVA WHATMAN	Cat No 2105-862 105 Sheets, Lens Cleaning Tissue, 200 x 300mm 100/pk Whatman	₹5,565
F229	2300772	CYTIVA WHATMAN	Cat No 2300-772 Benchkote Surface Protector, Reel, 920mm x 50 metres 1/pk Whatman	₹66,120
F230	2300916	CYTIVA WHATMAN	Cat No 2300-916 Benchkote Surface Protector, Sheets, 460 x 570mm 50/pk Whatman	₹12,705

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F231	23016150	CYTIVA WHATMAN	Cat No 2301-6150 Benchkote Plus Surface Protector, Sheets, 500 x 600mm 50/pk Whatman	₹19,952
F232	23016160	CYTIVA WHATMAN	Cat No 2301-6160 Benchkote Plus Surface Protector, Reel, 600mm x 50 metres 1/pk Whatman	₹32,611
F233	2814199	CYTIVA WHATMAN	Cat No 2814-199 Glass Thimbles, HPGF, 19 x 90mm 25/pk Whatman	₹36,874
F234	10370394	CYTIVA WHATMAN	Cat No 10370394 GF10, 50mm x 100m Core25mm 1/pk Whatman	₹2,95,554

 **YSP- SCIENTIFIC SOLUTION**

WHATMAN - Specimen Collection Devices (903 DBS Cards)




SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F235	10534612	CYTIVA WHATMAN	Cat No 10534612 903 PROTEINSAVER US 100/PK Whatman	₹31,654
F236	10535097	CYTIVA WHATMAN	Cat No 10535097 903 Sheets 210 x 297mm 100/pk Whatman	₹68,672
F237	10538017	CYTIVA WHATMAN	Cat No 10538017 903 Sheets 46 x 57cm 100/pk Whatman	₹1,47,697
F238	10538018	CYTIVA WHATMAN	Cat No 10538018 903 Sheets 58 x 58cm 100/pk Whatman	₹1,19,872

WHATMAN - MEMBRANE FILTER



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F241	99091302	CYTIVA WHATMAN	Cat No 9909-1302 UNIFLO 13/0.2 PVDF 500/PK Whatman	₹24,259
F242	99091304	CYTIVA WHATMAN	Cat No 9909-1304 UNIFLO 13/0.45 PVDF500/PK Whatman	₹24,259
F243	99092502	CYTIVA WHATMAN	Cat No 9909-2502 UNIFLO 25/0.2 PVDF 500/PK Whatman	₹36,801
F244	99092504	CYTIVA WHATMAN	Cat No 9909-2504 UNIFLO 25/0.45 PVDF 500/PK Whatman	₹34,177
F245	99101302	CYTIVA WHATMAN	Cat No 9910-1302 UNIFLO 13/0.2 NYL 500/PK Whatman	₹23,679
F246	99101304	CYTIVA WHATMAN	Cat No 9910-1304 UNIFLO 13/0.45 NYL 500/PK Whatman	₹23,679
F247	99102502	CYTIVA WHATMAN	Cat No 9910-2502 UNIFLO 25/0.2 NYL 500/PK Whatman	₹26,28
F248	99102504	CYTIVA WHATMAN	Cat No 9910-2504 UNIFLO 25/0.45 NYL 500/PK Whatman	₹26,289
F249	99111302	CYTIVA WHATMAN	Cat No 9911-1302 UNIFLO 13/0.2 PTFE 500/PK Whatman	₹28,638
F250	99111304	CYTIVA WHATMAN	500/PK Whatman	₹28,638
F251	99112502	CYTIVA WHATMAN	Cat No 9911-2502 UNIFLO 25/0.2 PTFE 500/PK Whatman	₹39,005
F252	99112504	CYTIVA WHATMAN	Cat No 9911-2504 UNIFLO 25/0.45 PTFE 500/PK Whatman	₹38,382

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F253	99121302	CYTIVA WHATMAN	Cat No 9912-1302 UNIFLO 13/0.2 PES 500/PK Whatman	₹35,337
F254	99121304	CYTIVA WHATMAN	Cat No 9912-1304 UNIFLO 13/0.45 PES 500/PK Whatman	₹35,337
F255	99122502	CYTIVA WHATMAN	Cat No 9912-2502 UNIFLO 25/0.2 PES 500/PK Whatman	₹35,931
F256	99122504	CYTIVA WHATMAN	Cat No 9912-2504 UNIFLO 25/0.45 PES 500/PK Whatman	₹35,931
F257	9917-2504	CYTIVA WHATMAN	Cat no-9917-2504 UNIFLO 25 MM 0.45 PVDF NON STERILE PACK OF 100	₹9,135
F258	9918-2504	CYTIVA WHATMAN	Cat no-9918-2504 UNIFLO 25 MM 0.45 NYLON NON STERILE PACK OF 100	₹7,543

 **YSP- SCIENTIFIC SOLUTION**
**WHATMAN SYRINGE
 FILTERS STERILE**



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F259	99132502	CYTIVA WHATMAN	Cat No 99132502 UNIFLO 25 mm 0.22 PVDF Strtile Pak of 45 Whatman	₹6,668
F260	99132504	CYTIVA WHATMAN	CatNo 9913 2504 UNIFLO 25 mm 0.45 PVDF Strtile Pak of 45 Whatman	₹6,668
F261	99142502	CYTIVA WHATMAN	Cat No 9914 2502 UNIFLO 25 mm 0.22 PES Strtile Pak of 45 Whatman	₹6,563
F262	99142504	CYTIVA WHATMAN	CatNo 9914 2504 UNIFLO 25 mm 0.45 PES Strtile Pak of 45 Whatman	₹6,563
F263	99152502	CYTIVA WHATMAN	Ca No 9915 2502 UNIFLO 25 mm 0.22 PES Strtile Pak of 200 Whatman	₹25,752

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F264	99152504	CYTIVA WHATMAN	Cat No 9915 2504 UNIFLO 25 mm 0.45 PES Strtile Pak of 200Whatman	₹25,752
F265	99161302	CYTIVA WHATMAN	Cat No 9916 2502 UNIFLO 13 mm 0.22 PES Strtile Pak of 100Whatman	₹15,418
F266	99161304	CYTIVA WHATMAN	Cat No 9916 2504 UNIFLO 13 mm 0.45 PES Strtile Pak of 100Whatman	₹14,543

 **YSP- SCIENTIFIC SOLUTION**
FILTER HOLDER -
GLOLAB



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F267	GLFH47	GLOLAB	Stainless steel Funnel Type Filter Holder 47MM / 300ml capacity without Filter flask	₹4,900
F268	GLFH47GL4 5	GLOLAB	Glass With Filter Flask & Clamp 47MM GL- 45 with cap & adaptor 300ml, alluminium cap GLOLAB	₹6,125
F269	GLFH47B40	GLOLAB	Glass With Filter Flask & Clamp 47MM with B40 joint and filter flask GLOLAB	₹6,125
F270	GLFH47PLAS TIC	GLOLAB	Spare Glass filter holder 47MM with Plastic adaptor GLOLAB	₹3,938
F271	FH-3 HOLDERS	GLOLAB	Filter Holder 47 MM 3 IN 1 SET W/O S.S HOLDER (47 mm/300 ml)	₹21,750
F272		GLOLAB	Spare Glass filter holder 47MM FLASK B- 40 GLOLAB	₹2,625

SDI KIT



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F273	GLSDIKIT	GLOLAB	Filter Holder All Glass 47 mm, Vacuum Pump Oil Free TID15(1/20hp) ,Whatman Cat.no.10406870, 47mm dia 0.45um, 5 meter silicon tube 8/11	₹26,825

TSS KIT



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F274		GLOLAB	Filter Holder All Glass 47 mm ,Vacuum Pump Oil Free TID 15,Whatman GF/A (1.6)& GF/C 47MM (1.2um) & 934-AH (1.5) 5 meter silicon tube 8/11	POR

SILT DENSITY INDEX



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F275		GLOLAB	SDI Kit Mini with filtration unit, 47mm membranes, ball valve, gauge and regulator without vacuum pump	₹23,200

FUME HOOD



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F276		YSP	YSPF 0900 Fume Exhaust Hood G I sheet powder coated working size 3 x 2x 2	₹1,41,375

FOGTECH ULV FOGGERS / FUMIGATOR



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F277		GENERAL	FOGTECH NEO AB housing and HDPE Tank,Tank Capacity 8 Ltrs Mechanical Timer 1-60 Minutes	₹31,900
F278		GENERAL	FOGTECH NEO MAXX ABS housing and SS 304 Grade Tank,Tank Cap 6 Ltrs Mechanical Timer 1-60 Minutes	₹36,801
F279		GENERAL	FOGTECH PRO Complete SS 304 Grade Tank,Capacity 6 Ltrs Mechanical Timer 1-60 Minutes	₹40,238
F280		GENERAL	FOGTECH ADVANCE Complete Stainless Steel 304 Grade Tank,Tank Capacity 6 Ltrs Digital Timer 1-99 Minutes with Delay Timer for Operator Safety	₹42,050
F281	FUMIGATOR	GENERAL	Aerosol Fumigator 5 Liter capacity SS body	₹24,288

FRYING OIL MONITOR - ATAGO



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
F282	DOM24	ATAGO	Model DOM-24 used to periodically check and monitor the quality of oil, Total Polar Materials(TPM): 0.5 to 40.0% Acid Value(AV): 0.00 to 9.99,Temp.: +0 to+225°C (In continuous mode +0 to +180°C)	₹63,075

GLASS DESICCATORS CHINA MAKE



KNOBBLED COVER, PORCELAIN PLATE CHINA

KNOBBLED COVER, PORCELAIN PLATE CHINA				
G1		GENERAL	Desiccator China Make Glass Plain 100mm	₹1,925
G2		GENERAL	Desiccato China Make Glass Plain 120mm	₹2,100
G3		GENERAL	Desiccator China Make Glass Plain 150mm	₹2,450
G4		GENERAL	Desiccator China Make Glass Plain 210mm	₹4,375
G5		GENERAL	Desiccator China Make Glass Plain 240mm	₹5,688
G6		GENERAL	Desiccator China Make Glass Plain 300mm	₹8,313
G7		GENERAL	Desiccator China Make Glass Plain 400mm	₹17,400
G8		GENERAL	Vaccum Desiccator Glass Amber colour 150mm	₹7,000
G9		GENERAL	Vaccum Desiccator Glass Amber colour 210mm	₹12,250

GLASS VACCUM DESICCATOR



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
G10		GENERAL	Vaccum Desiccator Glass 120mm	₹2,800
G11		GENERAL	Vaccum Desiccator Glass 150mm	₹3,500
G12		GENERAL	Vaccum Desiccator Glass 210mm	₹6,125
G13		GENERAL	Vaccum Desiccator Glass 240mm	₹7,875
G14		GENERAL	Vaccum Desiccator Glass 300mm	₹11,813

GLASS BEADS (PER BOTTLE 500 GMS)



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
G15		GENERAL	Glass Beads (Per Bottle 500 gms) 1mm	₹613
G16		GENERAL	Glass Beads (Per Bottle 500 gms) 2mm	₹569
G17		GENERAL	Glass Beads (Per Bottle 500 gms) 3mm	₹569
G18		GENERAL	Glass Beads (Per Bottle 500 gms) 4mm	₹569
G19		GENERAL	Glass Beads (Per Bottle 500 gms) 5mm	₹569
G20		GENERAL	Glass Beads (Per Bottle 500 gms) 6mm	₹569
G21		GENERAL	Glass Beads (Per Bottle 500 gms) 8mm	₹788



GERBER CENTRIFUGE FOR MILK

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
G22		GENERAL	Gerber Centrifuge electrical for Milk Butymeter mini for 8 tests	₹16,625
G23		GENERAL	Gerber Centrifuge electrical for Milk Butymeter mini for 12 tests	₹15,225
G24		GENERAL	Gerber Centrifuge electrical for Milk Butymeter mini for 24 tests	₹16,675

R-86 for R-8k Centrifuge REMI with R86 Head SEE "C" Centrifuge REMI

GAS ANALYSERS - TESTO



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
G25	TESTO3173	TESTO	Testo 317-3 CO Monitor CO sensor, Adjustable alarm threshold, CO zeroing at site, No initialization phase, instrument ready to operate immediately	POR
G26	TESTO3162	TESTO	Testo 316-2 Gas Leak Detector Detects CH ₄ , C ₃ H ₈ , H ₂ , Optical and audible alarm, Trend display shows maximum leakage, Integrated pump, flexible measurement probe in accessible places	POR

HOT AIR OVEN - YSP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H1	YSPOVEN12	YSP	YSP-OVEN-12 Chamber size 300x300x300mm 12&x12&x12&, capacity 28 Ltr, 750 Rating Shelves Analog Regulator Temp.	₹16,675
H2	YSPOVEN12DIGI	YSP	YSP-OVEN-12 DIGI 12&x12&x12& with Digital Temperature Controller	₹25,375
H3	YSPOVEN12DIGIFAN	YSP	YSP-OVEN-12 DIGI & FAN 12&x12&x12& with Digital Temperature Controller and Fan for Circulating air	₹31,900
H4	YSPOVEN14	YSP	YSP-OVEN-14 Chamber Size 355x355x355mm 14&x14&x14&, capacity 45 Ltr, 1000 Rating 2 Shelves Analog Regulator Temp.	₹23,925
H5	YSPOVEN14DIGI	YSP	YSP-OVEN-14 DIGI 14&x14&x14& with Digital Temperature Controller	₹31,900
H6	YSPOVEN14DIGIFAN	YSP	YSP-OVEN-14 DIGI & FAN 14&x14&x14& with Digital Temperature Controller and Fan for Circulating air	₹39,875
H7	YSPOVEN18	YSP	YSP-OVEN-18 Chamber size 455x455x455mm 18&x18&x18& capacity 95 Ltr, 1500 Rating, 2 Shelves, Analog Regulator Temp.	₹30,450
H8	YSPOVEN18DIGI	YSP	YSP-OVEN-18 DIGI 18&x18&x18 with Digital Temperature Controller	₹38,425

H9	YSPOVEN18DI GI&FAN	YSP	YSP-OVEN-18 DIGI & FAN 18'x18"x18" with Digital Temperature Controller and Fan for Circulating air	₹46,400
H10	YSPOVEN24	YSP	YSP-OVEN- 24 Inner Chamber 605x605x605mm 24&x24&x24& capacity 224 Ltr, 2200 Rating, 3 Shelves Analog Regulator Temp.	₹47,850
H11	YSPOVEN24DI GI	YSP	YSP-OVEN-24 DIGI 24&x24&x24& with Digital Temperature Controller	₹55,825
H12	YSPOVEN24DI GIFAN	YSP	YSP-OVEN-24 DIGI & FAN 24&x24&x24& with Digital Temperature Controller and Fan for Circulating air	₹63,800

 **YSP- SCIENTIFIC SOLUTION**

HOT AIR OVEN - REMI



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H13	RDHO50	REMI	RDHO-50 Dry Hot Air Oven (50 liters) with temperature setting from RT + 10°C up to +300°C	₹1,08,605
H14	RDHO80	REMI	RDHO-80 Dry Hot Air Oven (80 Liters) with temperature setting from RT + 10°C up to +300°C	₹1,24,120

 **YSP- SCIENTIFIC SOLUTION**

HOT PLATE - YSP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H15	YSPHP8MS	YSP	YSP-HP-8 MS Hot Plate8& Dia (1000W) M.S. (Energy Regulator)	₹4,813
H16	YSPHP9MS	YSP	YSP-HP-9 MS Hot Plate9& Dia (1000W) M.S. (Energy Regulator)	₹5,775
H17	YSPHP12MS	YSP	YSP-HP-12 MS Hot Plate12& Dia (1500W) M.S. (Energy Regulator)	₹7,175
H18	YSPHP1012MS	YSP	YSP-HP(10x12) MS Hot Plate10x12 (1000W) M.S. (Energy Regulator)	₹8,313
H19	YSPHP1016MS	YSP	YSP-HP(10x16) MS Hot Plate10x16 (1500W) M.S. (Energy Regulator)	₹8,313
H20	YSPHP1218MS	YSP	YSP-HP(12x18) MS Hot Plate12x18 (2000W) M.S. (Energy Regulator)	₹11,375
H21	YSPHP1824MS	YSP	YSP-HP(18x24) MS Hot Plate18x24 (4000W) M.S. (Energy Regulator)	₹17,400
H22	YSPHP8MSDIGI	YSP	YSP-HP-8 MS DIGI Hot Plate 8" Digital Controller MS BODY	₹15,225
H23	YSPHP1218MSDIGI	YSP	YSP-HP(12x18) MS DIGI Hot Plate12X18" Digital Controller M.S Digital	₹18,125
H24	YSPHP7MS	YSP	YSP-HP-7 MS Hot Plate7" [ORDINARY] M.S Body Commercial	₹2,625
H25	YSPHP8SS	YSP	YSP-HP-8 SS Hot Plate 8& Dia (1000W) S.S. (Energy Regulator)	₹7,350
H26	YSPHP9SS	YSP	YSP-HP-9 SS Hot Plate9& Dia (1000W) S.S. (Energy Regulator)	₹9,625
H27	YSPHP12SS	YSP	YSP-HP-12 SS Hot Plate12& Dia (1500W) S.S. (Energy Regulator)	₹11,813

H28	YSPHP1012SS	YSP	YSP-HP(10x12) SS Hot Plate 10x12 (1000W) S.S. (Energy Regulator)	₹11,813
H29	YSPHP1016SS	YSP	YSP-HP(10x16) SS Hot Plate 10x16 (1500W) S.S. (Energy Regulator)	₹12,600
H30	YSPHP1218SS	YSP	YSP-HP(12x18) SS Hot Plate 12x18 (2000W) S.S. (Energy Regulator)	₹14,000
H31	YSPHP1224SS	YSP	YSP-HP(12x24) SS Hot Plate 12x24 (4000W) S.S. (Energy Regulator)	₹24,650
H32	YSPHP8SSDIGI	YSP	YSP-HP-8 SS DIGI Hot Plate 8" Digital Controller SS body	₹17,500
H33	YSPHP1012SS DIGI	YSP	YSP-HP(10x12) SS DIGI Hot Plate 10x12 Digital Controller SS Body	₹18,488
H34	YSPHP1218SS DIGI	YSP	YSP-HP(12x18) SS DIGI Hot Plate 12x18" SS Top Digital Controller	₹23,563

 **YSP- SCIENTIFIC SOLUTION**

HEATING MANTLE - YSP



H35	YSPHM50	YSP	YSP-HM-50 Heating Mantels Capacity 50ml	₹2,975
H36	YSPHM100	YSP	YSP-HM-100 Heating Mantels Capacity 100ml	₹3,150
H37	YSPHM250	YSP	YSP-HM-250 Heating Mantels Capacity 250ml	₹3,325
H38	YSPHM500	YSP	YSP-HM-500 Heating Mantels Capacity 500ml	₹3,500
H39	YSPHM1000	YSP	YSP-HM-1000 Heating Mantels Capacity 1000ml	₹3,850

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H40	YSPHM2000	YSP	YSP-HM-2000 Heating Mantels Capacity 2000ml	₹4,375
H41	YSPHM3000	YSP	YSP-HM-3000 Heating Mantels Capacity 3000ml	₹5,950
H42	YSPHM5000	YSP	YSP-HM-5000 Heating Mantels Capacity 5000ml	₹6,650
H43	YSPHM10000	YSP	YSP-HM-10000 Heating Mantels Capacity 10000ml	₹8,400
H44	YSPHM20000	YSP	YSP-HM-20000 Heating Mantels Capacity 20000ml	₹11,900
H45	YSPSM50	YSP	Spare Mantel 50ml	₹963
H46	YSPSM100	YSP	Spare Mantel 100ml	₹1,050
H47	YSPSM250	YSP	Spare Mantel 250ml	₹1,225
H48	YSPSM500	YSP	Spare Mantel 500ml	₹1,444
H49	YSPSM1000	YSP	Spare Mantel 1000ml	₹1,663
H50	YSPSM2000	YSP	Spare Mantel 2000ml	₹1,925
H51	YSPSM3000	YSP	Spare Mantel 3000ml	₹2,319
H52	YSPSM5000	YSP	Spare Mantel 5000ml	₹3,430
H53	YSPSM10000	YSP	Spare Mantel 10000ml	₹5,303
H54	YSPSM20000	YSP	Spare Mantel 20000ml	₹6,475

HYDROMETERS PERCENTAGE - GLOSIL



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H55	GHCAUSTIC1	GLOSIL HYDROMETER	CAUSTIC % in 30-35, 35-40, 40-45, 45-50 % any type	₹919
H56	GHCAUSTIC2	GLOSIL HYDROMETER	CAUSTIC percentage meter 30-50%	₹429
H57	GHSULPHURIC	GLOSIL HYDROMETER	SULPHURIC 60-80%, 80-98%, 65-75%	₹429
H58	GHHYDROCLORIC	GLOSIL HYDROMETER	HYDROCLORIC % Meter 25-35%	₹429
H59	GHNITRIC	GLOSIL HYDROMETER	NITRIC % 50-75%, 75-100%	₹429
H60	GHACETIC	GLOSIL HYDROMETER	ACETIC ACID % 90-100%	₹429
H61	GHHYDROGENPER OXIDE	GLOSIL HYDROMETER	HYDROGE PEROXIDE % 30-50%	₹613
H62	GHSALINOMETER	GLOSIL HYDROMETER	SALINOMETER (BRINOMETER)	₹788
H63	GHCARTIER	GLOSIL HYDROMETER	CARTIER 10-44	₹394
H64	GHLATEXOMETER	GLOSIL HYDROMETER	LATEXOMETER	₹569

HYDROMETER - GLOSIL



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
SPECIFIC GRAVITY,BAUME,BRIX,SOIL,SALINO,ALCOHOL,SIKES,TWADDEL INDIAN				
H65	GHSP1	GLOSIL HYDROMETER	Hydrometer Range SP,Gravity 1.000-2.000 (0-70) (For Heavy Liquid)	₹193
H66	GHSP2	GLOSIL HYDROMETER	Hydromete Range SP,Gravity 0.700 - 1.000 (10- 70) (For Heavy Light Liquid)	₹228
H67	GHSP3	GLOSIL HYDROMETER	Hydrometer Range SP,Gravity 1.000 - 1.500 (0- 70) (For Heavy Liquid bulb type)	₹333
H68	GHSP4	GLOSIL HYDROMETER	Hydrometer Range SP,Gravity 1.500 - 2.000 & BE	₹228
H69	GHSP5	GLOSIL HYDROMETER	Hydrometer Range SP,Gravity 1.000 - 1.300 & BE	₹263
H70	GHSP6	GLOSIL HYDROMETER	Hydrometer 0.002 Range SP,Gravity 0.600 - 0.800	₹306
H71	GHSP7	GLOSIL HYDROMETER	Hydrometer 0.002 Range SP,Gravity 0.800 - 1.000	₹228
H72	GHSP8	GLOSIL HYDROMETER	Hydrometer 0.002 Range SP,Gravity 1.000 - 1.200	₹228
H73	GHSP9	GLOSIL HYDROMETER	Hydrometer 0.002 Range SP,Gravity 1.200 - 1.400	₹228
H74	GHSP10	GLOSIL HYDROMETER	Hydrometer 0.002 Range SP,Gravity 1.400 - 1.600	₹228

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H76	GHSP12	GLOSIL HYDROMETER	Hydrometer 0.002 Range SP,Gravity 1.800 - 2.000	₹394
H77	GHSP13	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 0.600 - 0.700	₹394
H78	GHSP14	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 0.700 - 0.800	₹228
H79	GHSP15	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 0.800 - 0.900 / 0.900 - 1.000	₹228
H80	GHSP16	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 1.000 - 1.100	₹228
H81	GHSP17	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 1.100 - 1.200	₹228
H82	GHSP18	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 1.200 - 1.300	₹228
H83	GHSP19	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 1.300 - 1.400	₹228
H84	GHSP20	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 1.400-1.500	₹228
H85	GHSP21	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 1.500-1.600	₹228
H86	GHSP22	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 1.600 - 1.700	₹394
H87	GHSP23	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 1.700 - 1.800	₹394
H88	GHSP24	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 1.800 - 1.900	₹394
H89	GHSP25	GLOSIL HYDROMETER	Hydrometer 0.001 Range SP,Gravity 1.900 - 2.000	₹394

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H90	GHSP26	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 0.600 - 0.650	₹394
H91	GHSP27	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 0.650 - 0.700	₹394
H92	GHSP28	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 0.700 - 0.750	₹306
H93	GHSP29	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 0.750 - 0.800	₹306
H94	GHSP30	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 0.800 - 0.850	₹306
H95	GHSP31	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 0.850 - 0.900	₹306
H96	GHSP32	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 0.900 - 0.950	₹306
H97	GHSP33	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 0.950 - 1.000	₹306
H98	GHSP34	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 1.000 - 1.050	₹394
H99	GHSP35	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 1.050 - 1.100	₹394
H100	GHSP36	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 1.100 - 1.150	₹394
H101	GHSP37	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 1.150 - 1.200	₹394
H102	GHSP38	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 1.200 - 1.250	₹394
H103	GHSP39	GLOSIL HYDROMETER	Hydrometer 0.005 Range SP,Gravity 1.250 - 1.300	₹394

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
BAUME Hydrometer				
H104	GHBE1	GLOSI HYDROMETER	Hydrometer 0 - 10	₹228
H105	GHBE2	GLOSIL HYDROMETER	Hydrometer 10 - 20	₹228
H106	GHBE3	GLOSIL HYDROMETER	Hydrometer 20-30	₹228
H107	GHBE4	GLOSIL HYDROMETER	Hydrometer 30-40	₹228
H108	GHBE5	GLOSIL HYDROMETER	Hydrometer 40-50	₹228
H109	GHBE6	GLOSIL HYDROMETER	Hydrometer 50-60	₹228
H110	GHBE7	GLOSIL HYDROMETER	Hydrometer 60-70	₹228
H111	GHBE8	GLOSIL HYDROMETER	Hydrometer 0-20	₹228
H112	GHBE9	GLOSIL HYDROMETER	Hydrometer 20-40	₹228
H113	GHBE10	GLOSIL HYDROMETER	Hydrometer 40-60	₹228
H114	GHBE11	GLOSIL HYDROMETER	Hydrometer 0-30	₹228
H115	GHBE12	GLOSIL HYDROMETER	Hydrometer 30-60	₹228
H116	GHBE13	GLOSIL HYDROMETER	Hydrometer 0-40	₹228
H117	GHBE14	GLOSIL HYDROMETER	Hydrometer 0-50	₹228

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
BRIX Hydrometer				
H118	GHBX1	GLOSI HYDROMETER	BRIX Hydrometer 0-10	₹315
H119	GHBX2	GLOSIL HYDROMETER	BRIX Hydrometer 10-20	₹315
H120	GHBX3	GLOSIL HYDROMETER	BRIX Hydrometer 20-30	₹315
H121	GHBX4	GLOSIL HYDROMETER	BRIX Hydrometer 0-30	₹315
H122	GHBX5	GLOSIL HYDROMETER	BRIX Hydrometer 30-40	₹315
H123	GHBX6	GLOSIL HYDROMETER	BRIX Hydrometer 30 - 60	₹315
H124	GHBX7	GLOSIL HYDROMETER	BRIX Hydrometer 40 - 50	₹315
H125	GHBX8	GLOSIL HYDROMETER	BRIX Hydrometer 50 - 60	₹315
H126	GHBX9	GLOSIL HYDROMETER	BRIX Hydrometer 60 - 70	₹315
H127	GHBX10	GLOSIL HYDROMETER	BRIX Hydrometer 60 - 90	₹315

**We specialize in the sales of all
laboratory products.**

For inquiries or orders,



**96262 97366
73387 21875**



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
BRIX With Thermo				
H128	GHBXT1	GLOSIL HYDROMETER	BRIX With Thermo 0 - 10	₹963
H129	GHBXT2	GLOSIL HYDROMETER	BRIX With Thermo 10 - 20	₹963
H130	GHBXT3	GLOSIL HYDROMETER	BRIX With Thermo 20 - 30	₹963
H131	GHBXT4	GLOSIL HYDROMETER	BRIX With Thermo 0 - 30	₹963
H132	GHBXT5	GLOSIL HYDROMETER	BRIX With Thermo 30 - 60	₹963
H133	GHBXT6	GLOSIL HYDROMETER	BRIX With Thermo 60 - 90	₹963

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
SOIL Hydrometer				
H134	GSHS	GLOSIL HYDROMETER	SOIL Hydrometer 0.995 TO 1.030	₹744

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Alcohol Meter				
H135	GHAH	GLOSIL HYDROMETER	0-100% Alcoholometer Gay Lussac	₹175

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

LACTOMETER / URINOMETER GLASS HYDROMETERS

H136	GHIJ	GLOSIL HYDROMETER	1.000-1.040 Lactometer with jar	₹140
H137	GHL	GLOSIL HYDROMETER	1.000-1.040 Lactometer only	₹114
H138	GHU	GLOSIL HYDROMETER	1.000-1.060 Urinometer with jar	₹140

Sikes Hydrometers

H139	GHSI1	GLOSIL HYDROMETER	Sikes 0 -20A	₹718
H140	GHSI2	GLOSIL HYDROMETER	Sikes 0 - 20	₹718
H141	GHSI3	GLOSIL HYDROMETER	Sikes 20 - 40	₹718
H142	GHSI4	GLOSIL HYDROMETER	Sikes 40 - 60	₹718
H143	GHSI5	GLOSIL HYDROMETER	Sikes 60 - 80	₹718
H144	GHSI6	GLOSIL HYDROMETER	Sikes 80 - 100	₹718



The LACTOMETER/URINOMETER glass hydrometer exceeded our expectations with its precision and durability. It's an essential tool for our testing needs.



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

Barko Meter

H145	GHBM	GLOSIL HYDROMETER	Barko Meter 0 -120	₹744
------	------	----------------------	--------------------	-------------

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

Twaddle

H146	GHT1	GLOSIL HYDROMETER	Twaddle No.1 (0 - 26)	₹228
H147	GHT2	GLOSIL HYDROMETER	Twaddle No.2 (24 - 50)	₹228
H148	GHT3	GLOSIL HYDROMETER	Twaddle No.3 (48 - 74)	₹228
H149	GHT4	GLOSIL HYDROMETER	Twaddle No.4 (72 - 104)	₹280
H150	GHT5	GLOSIL HYDROMETER	Twaddle No.5 (104 - 138)	₹280
H151	GHT6	GLOSIL HYDROMETER	Twaddle No.6 (136 - 174)	₹280
H152		GLOSIL HYDROMETER	Battery Hydrometer Big Leimco-Range 1.100 to 1.300 specific Gravity	₹341
H153		GLOSIL HYDROMETER	Calibration Certificate by private for Hydrometer BAUME and Specific Gravity	₹438

Calibration Certificate charges for Brix hydro With Thermometer

Calibration Certificate charges for Brix hydrrometer

Calibration Certificate charges for Sikes hydrrometer

HYDROMETER - LEIMCO



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
SPECIFIC GRAVITY,BAUME,BRIX,SOIL,SALINO,ALCOHOL,SIKES,TWADDEL INDIAN (LEIMCO / LIEMCO BRANDS) Any Available One				
H154	LE1	LEIMCO	Hydrometer Range SP,Gravity 1.000-2.000 (0-70) LEIMCO	₹481
H155	LE2	LEIMCO	Hydrometer Range SP,Gravity .700 - 1.000 LEIMCO	₹481
H156	LE3	LEIMCO	Hydrometer Range SP,Gravity 1.000 - 1.500 LEIMCO	₹481
H157	LE4	LEIMCO	Hydrometer Range SP,Gravity 1.500 - 2.000 LEIMCO	₹481
H158	LE5	LEIMCO	Hydrometer 0.002 Range 0.600-.800 LEIMCO	₹735
H159	LE6	LEIMCO	Hydrometer 0.002 Range 0.800 - 1.000 LEIMCO	₹481
H160	LE7	LEIMCO	Hydrometer 0.002 Range 1.000 - 1.200 LEIMCO or 0.800 - 1.000	₹481
H161	LE8	LEIMCO	Hydrometer 0.002 Range SP,Gravity 1.200 - 1.400 LEIMCO	₹481
H162	LE9	LEIMCO	Hydrometer 0.002 Range SP,Gravity 1.400 - 1.600 LEIMCO	₹481
H163	LE10	LEIMCO	Hydrometer 0.002 Range SP,Gravity 1.600 - 1.800 LEIMCO	₹735
H164	LE11	LEIMCO	Hydrometer 0.002 Range SP,Gravity 1.800 - 2.000 LEIMCO	₹735

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H165	LE12	LEIMCO	Hydrometer 0.001 Range 0.600 -.700	₹735
H166	LE13	LEIMCO	Hydrometer 0.001 Range 0.700 - .800	₹481
H167	LE14	LEIMCO	Hydrometer 0.001 Range 0.800 - 0.900 / 0.900 - 1.000	₹481
H168	LE15	LEIMCO	Hydrometer 0.001 Range 1.000 - 1.100	₹481
H169	LE16	LEIMCO	Hydrometer 0.001 Range 1.100 - 1.200	₹481
H170	LE17	LEIMCO	Hydrometer 0.001 Range 1.200 - 1.300	₹481
H171	LE18	LEIMCO	Hydrometer 0.001 Range 1.300 - 1.400	₹481
H172	LE19	LEIMCO	Hydrometer 0.001 Range 1.400-1.500	₹648
H173	LE20	LEIMCO	Hydrometer 0.001 Range 1.500-1.600	₹648
H174	LE21	LEIMCO	Hydrometer 0.001 Range 1.600 - 1.700	₹648
H175	LE22	LEIMCO	Hydrometer 0.001 Range 1.700 - 1.800	₹648
H176	LE23	LEIMCO	Hydrometer 0.001 Range 1.800 - 1.900	₹735
H177	LE24	LEIMCO	Hydrometer 0.001 Range 1.900 - 2.000	₹735
H178	LE25	LEIMCO	Hydrometer 0.001 Range 0.600 - 0.650	₹1,050

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H179	LE26	LEIMCO	Hydrometer 0.001 Range 0.650 - 0.700	₹735
H180	LE27	LEIMCO	Hydrometer 0.001 Range 0.700 - 0.750	₹735
H181	LE28	LEIMCO	Hydrometer 0.001 Range 0.750 - 0.800	₹735
H182	LE29	LEIMCO	Hydrometer 0.001 Range 0.800 - 0.850	₹735
H183	LE30	LEIMCO	Hydrometer 0.001 Range 0.850 - 0.900	₹735
H184	LE31	LEIMCO	Hydrometer 0.001 Range 0.900 - 0.950	₹735
H185	LE32	LEIMCO	Hydrometer 0.001 Range 0.950 - 1.000	₹735
H186	LE33	LEIMCO	Hydrometer 0.001 Range 1.000 - 1.050	₹788
H187	LE34	LEIMCO	Hydrometer 0.001 Range 1.050 - 1.100	₹788
H188	LE35	LEIMCO	Hydrometer 0.001 Range 1.100 - 1.150	₹788
H189	LE36	LEIMCO	Hydrometer 0.001 Range 1.150 - 1.200	₹788
H190	LE37	LEIMCO	Hydrometer 0.001 Range 1.200 - 1.250	₹788
H191	LE38	LEIMCO	Hydrometer 0.001 Range 1.250 - 1.300	₹788

We're your go-to source for all laboratory products.



Probe Sonicator



ice maker



ph meter

Upgrade your lab with top-quality equipment! Get the Probe Sonicator for efficient sample preparation, an Ice Maker for reliable cooling, and a PH Meter for accurate readings. Everything you need for precise lab work!

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
BAUME Hydrometer				
H192	LE39	LEIMCO	Hydrometer 0 -10	₹481
H193	LE40	LEIMCO	Hydrometer 10 - 20	₹481
H194	LE41	LEIMCO	Hydrometer 20-30	₹481
H195	LE42	LEIMCO	Hydrometer 30-40	₹481
H196	LE43	LEIMCO	Hydrometer 40-50	₹656
H197	LE44	LEIMCO	Hydrometer 50-60	₹656
H198	LE45	LEIMCO	Hydrometer 60-70	₹656
H199	LE46	LEIMCO	Hydrometer 0-20	₹525
H200	LE47	LEIMCO	Hydrometer 20-40	₹525
H201	LE48	LEIMCO	Hydrometer 40-60	₹525
H202	LE49	LEIMCO	Hydrometer 0-30	₹481
H203	LE50	LEIMCO	Hydrometer 30-60	₹525
H204	LE51	LEIMCO	Hydrometer 0-40	₹525
H205	LE52	LEIMCO	Hydrometer 0-50	₹525
H206	LE53	LEIMCO	Hydrometer 0-70	₹525

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
BRIX Hydrometer				
H207	LE54	LEIMCO	BRIX Hydrometer 0-10	₹875
H208	LE55	LEIMCO	BRIX Hydrometer 10-20	₹875
H209	LE56	LEIMCO	BRIX Hydrometer 20-30	₹875

H210	LE57	LEIMCO	BRIX Hydrometer 0-30	₹770
H211	LE58	LEIMCO	BRIX Hydrometer 30-40	₹919
H212	LE59	LEIMCO	BRIX Hydrometer 30 - 60	₹770
H213	LE60	LEIMCO	BRIX Hydrometer 40 - 50	₹919
H214	LE61	LEIMCO	BRIX Hydrometer 50 - 60	₹1,006
H215	LE62	LEIMCO	BRIX Hydrometer 60 - 70	₹1,006
H216	LE63	LEIMCO	BRIX Hydrometer 60 - 90	₹936
H217	LE64	LEIMCO	BRIX Hydrometer 90-100	₹1,006

BRIX With Thermo

H218	LE65	LEIMCO	BRIX With Thermo 0 - 10	₹1,811
H219	LE66	LEIMCO	BRIX With Thermo 10 - 20	₹1,811
H220	LE67	LEIMCO	BRIX With Thermo 20 - 30	₹1,969
H221	LE68	LEIMCO	BRIX With Thermo 0 - 30	₹1,969
H222	LE69	LEIMCO	BRIX With Thermo 30 - 60	₹1,969
H223	LE70	LEIMCO	BRIX With Thermo 60 - 90	₹1,969

SOIL Hydrometer

H224	LE71	LEIMCO	SOIL Hydrometer 0.995 TO 1.030 LEIMCO	₹1,969
------	------	--------	------------------------------------------	--------



Shop our high-precision weights E1, E2, F1, F2, & M1

Perfect for calibration and lab accuracy. Ensure reliable measurements every time



96262 97366 - 73387 21875

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

Alcohol Meter

H225	LE72	LEIMCO	Alcohol Meter 0 - 100% (Gay Lucas)	₹481
H226	LE73	LEIMCO	Alcohol Meter 90 - 100%	₹1,698
H227	LE74	LEIMCO	Alcohol Meter 94 - 104	₹1,698
SL	CODE	BRAND NAME	PRODUCT NAME	RATE

Sikes

H228	LE75	LEIMCO	Sikes 0 -20A	₹1,138
H229	LE76	LEIMCO	Sikes 0 - 20	₹1,138
H230	LE77	LEIMCO	Sikes 20 - 40	₹1,208
H231	LE78	LEIMCO	Sikes 40 - 60	₹1,208
H232	LE79	LEIMCO	Sikes 60 - 80	₹1,208
H233	LE80	LEIMCO	Sikes 80 - 100	₹1,208
H234	LE81	LEIMCO	Sikes TABLE BOOK	₹831
SL	CODE	BRAND NAME	PRODUCT NAME	RATE

Barko Meter

H235	LE82	LEIMCO	Barko Meter 0 - 120	₹1,356
------	------	--------	---------------------	--------

We specialize in the sales of **all** laboratory products.

For inquiries or orders,



96262 97366


73387 21875



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Twaddle				
H236	LE83	LEIMCO	Twaddle No.1 (0 - 26)	₹525
H237	LE84	LEIMCO	Twaddle No.2 (24 - 50)	₹525
H238	LE85	LEIMCO	Twaddle No.3 (48 - 74)	₹525
H239	LE86	LEIMCO	Twaddle No.4 (72 - 104)	₹648
H240	LE87	LEIMCO	Twaddle No.5 (102 - 138)	₹648
H241	LE88	LEIMCO	Twaddle No.6 (136 - 174)	₹648
H242	LE89	LEIMCO	Twaddle No 7 (160-200)	₹823
H243	LE90	LEIMCO	Battery Hydrometer Big Leimco-Range 1.100 to 1.300 specific Gravity	₹858
Calibration Certificate by private for Hydrometer BAUME and Specific Gravity				
Calibration Certificate charges for Brix hydro With Thermometer				
Calibration Certificate charges for Brix hydrometer				
Calibration Certificate charges for Sikes hydrometer				

 **YSP- SCIENTIFIC SOLUTION**

GERMAN BRIX WITH THERMO



H244			GERMAN BRIX WITH THERMO 0-10% 27.5'C (ALCO BLUE)	₹6,738
H245			GERMAN BRIX WITH THERMO 10-20 27.5'C (ALCO BLUE)	₹6,738

HYDROMETER-ZEAL



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H246		ZEAL	Zeal Sikes Hydrometer B144/1 (0-20)	₹3,063
H247		ZEAL	Zeal Sikes Hydrometer B144/1A (0/20A) (0.2)	₹3,063
H248		ZEAL	Zeal Sikes Hydrometer B/144/2 (20-40)	₹3,063
H249		ZEAL	Zeal Sikes Hydrometer B/144/3 (40-60)	₹3,063
H250		ZEAL	Zeal Sikes Hydrometer B/144/4 (60-80)	₹3,063
H251		ZEAL	Zeal Sikes Hydrometer B/144/5 (80-100)	₹3,063

LACTOMETER / URINOMETER - LEIMCO



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H252	LE95	LEIMCO	Lactometer Liemco/ Indian make 1.000 - 1.040	₹481
H253	LE96	LEIMCO	Urinometer 1.000 -1.060 / 0.001	₹805

DN WITH LAB CERTIFICATE FOR 5 YEARS - LEIMCO



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H254	LE97	LEIMCO	DN M50 700-750 LEIMCO ,MUMBAI	₹1,006
H255	LE98	LEIMCO	DN M50 800-850 LEIMCO ,MUMBAI	₹1,006
H256	LE99	LEIMCO	Swastik Thermometer 50°C M (0.5)	₹788
H257	LE100	LEIMCO	OLEODIS OIL FINDING PASTE	₹245
H258	LE101	LEIMCO	AQUADIS WATER FINDING PASTE	₹245

Density Kit Small

H259	LE102	LEIMCO	2 Hydrometer (700 to 750 & 800 to 850), 1 no. Thermometer 50°C (0.5) with certificate and 1 measuring cylinder 500ml. Borosilicate glass and 1 each of water paste & oil pase	₹6,125
------	-------	--------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------



Durable and easy to clean, our SS 304 Table Top and Bench Scale are perfect for the pharmaceutical and chemistry industries. Reliable precision for your lab needs!



LEIMCO - DN 'M' SERIES HYDROMETER



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H260	LE103	LEIMCO	DN Range 0.600 to 0.650 LEIMCO M-50(0.001)	₹998
H261	LE104	LEIMCO	DN Range 0.650 to 0.700 LEIMCO M-50(0.001)	₹665
H262	LE105	LEIMCO	DN Range 0.700 to 0.750 LEIMCO M-50(0.001)	₹630
H263	LE106	LEIMCO	DN Range 0.750 to 0.800 LEIMCO M-50(0.001)	₹630
H264	LE107	LEIMCO	DN Range 0.800 to 0.850 LEIMCO M-50(0.001)	₹630
H265	LE108	LEIMCO	DN Range 0.850 to 0.900 LEIMCO M-50(0.001)	₹665
H266	LE109	LEIMCO	DN Range 0.900 to 0.950 LEIMCO M-50(0.001)	₹665
H267	LE110	LEIMCO	DN Range 0.950 to 1.000 LEIMCO M-50(0.001)	₹665
H268	LE111	LEIMCO	DN Range 1.000 to 1.050 LEIMCO M-50(0.001)	₹998
H269	LE112	LEIMCO	DN Range 1.0500 to 1.100 LEIMCO M-50(0.001)	₹998

LEIMCO - DN 'L' SERIES HYDROMETER



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H270	LE113	LEIMCO	DN Range 0.600 to 0.650 LEIMCO L-50(0.0005) (CERTIFICATE)	₹2,625
H271	LE114	LEIMCO	DN Range 0.650 to 0.700 LEIMCO L-50(0.0005) (CERTIFICATE)	₹2,013
H272	LE115	LEIMCO	DN Range 0.700 to 0.750 LEIMCO L-50(0.0005) (CERTIFICATE)	₹2,013
H273	LE116	LEIMCO	DN Range 0.750 to 0.800 LEIMCO L-50(0.0005) (CERTIFICATE)	₹2,013
H274	LE117	LEIMCO	DN Range 0.800 to 0.850 LEIMCO L-50(0.0005) (CERTIFICATE)	₹2,013
H275	LE118	LEIMCO	DN Range 0.850 to 0.900 LEIMCO L-50(0.0005) (CERTIFICATE)	₹2,013
H276	LE119	LEIMCO	DN Range 0.900 to 0.950 LEIMCO L-50(0.0005) (CERTIFICATE)	₹2,013
H277	LE120	LEIMCO	DN Range 0.950 to 1.000 LEIMCO L-50(0.0005) (CERTIFICATE)	₹2,013
H278	LE121	LEIMCO	DN Range 1.000 to 1.050 LEIMCO L-50(0.0005) (CERTIFICATE)	₹2,625
H279	LE122	LEIMCO	DN Range 1.0500 to 1.100 LEIMCO L-50(0.0005) (CERTIFICATE)	₹2,625

CERTIFICATION FOR HYDROMETER - LE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H280	LE122	LEIMCO	NABL CERTIFICATION FOR HYDROMETER	₹2,188
H281	LE123	LEIMCO	CERTIFICATE BY PRIVATE PARTY FOR NATIONAL STANDARD FOR HYDROMETER	₹613

HEAMOMETER ROUND - TOP TECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H282	HEAMOMETER	TOP TECH	Heamometer Round Top Tech full set	₹1,050
H283	Round Tube Top Tech	TOP TECH	Round Tube Top Tech	₹96



Durable and easy to clean, our SS 304 Table Top and Bench Scale are perfect for the pharmaceutical and chemistry industries. Reliable precision for your lab needs!



HEAMOMETER - HBG



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H284	HEAMOMETER	HBG	Heamometer SAHLIS HBG Make complete set with comparator, round tube 20cm pipette brush	₹1,400
H285	ROUND TUBE	HBG	Round Tube HBG Make	₹219
H286	20CM TUBE	HBG	20 cm Pipette HBG Make	₹219

HEAMOMETER AND HEAMOCYTOMETER ACCESSORIES - TOP TECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H287	T18	TOP TECH	RBC Pipette Top Tech	₹140
H288	T19	TOP TECH	WBC Pipette Top Tech	₹140
H289	T20	TOP TECH	20cm Pipette Top Tech	₹79
H290	T21	TOP TECH	Wintrobe Tube Top Tech	₹70

HEAMOCYTOMETER SET - HBG



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H291	HEAMOCYTO METER	HBG	Heamocytometer HBG with Counting Chamber, HBG WBC & RBC full set	₹2,713

HEAMOMETER SUPERIOR



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H292	HEAMOMETER	SUPERIOR	Heamometer Superior	₹2,538
H293		SUPERIOR	Round Tube HB	₹263
H294		SUPERIOR	20 cm Pipette Superior	₹394
H295		SUPERIOR	WBC Superior	₹744
H296		SUPERIOR	RBC Superior	₹744
H297		SUPERIOR	ESR Tube Superior	₹105

HEAMOCYTOMETER ACCESSORIES



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H298	COUNTING CHAMBER	HBG	Counting Chamber HBG Make	₹2,013
H299	H3	GENERAL	Counting Chamber Rohem Silverlite	₹2,013
H300	H4	GENERAL	Counting Chamber CHINA	₹744
H301	H5	GENERAL	Counting Chamber Accutek Brightline	₹1,575
H302	RBC AND WBC PIPETTE	HBG	RBC and WBC Pipette HBG Make	₹333
H303	H7	SUPERIOR	RBC and WBC Pipette SUPERIOR	₹656
	H8		RBC and WBC Pipette JAYNA	
H304	H9	GENERAL	20cm Pipette Top JAYNA	₹44
H305	WINDROP TUBE	HBG	Wintrobe Tube HBG Make	₹236
H306	H11	HBG	ESR Tube HBG Make	₹219

HEAMOCYTOMETER COMBINATION



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H307	H15	GENERAL	Accutek Chamber with HBG Pipette	₹3,675

H308	H16	GENERAL	Accutek Chamber with Indian Pipette	₹1,925
H309	H17	GENERAL	Accutek Chamber with Top Tech Pipette	₹2,013
H310	H18	GENERAL	China Chamber with HBG Pipette	₹1,356
H311	H19	GENERAL	China Chamber with Indian Pipette	₹1,094
H312	H20	GENERAL	China Chamber with Top Tech Pipette	₹1,225
H313	H21	GENERAL	HBG Chamber with HBG Pipette	₹2,625
H314	H22	GENERAL	HBG Chamber with Indian Pipette	₹2,188
H315	H23	GENERAL	HBG Chamber with Superior Pipette	₹3,456
H316	H24	GENERAL	HBG Chamber with Top Tech Pipette	₹2,363
H317	H25	GENERAL	Rohem/Top Tech Chamber with HBG Pipette	₹2,538
H318	H26	GENERAL	Rohem/Top Tech Chamber with Indian Pipette	₹2,231
H319	H27	GENERAL	Rohem/Top Tech Chamber with Superior Pipette	₹3,500
H320	H28	GENERAL	Rohem/Top Tech Chamber with Top Tech Pipette	₹2,450
H321	H29	GENERAL	Weber Chamber with HBG Pipette	₹4,900
H322	H30	GENERAL	Weber Chamber with Superior Pipette	₹5,906
H323	H31	GENERAL	Weber Chamber with Top Tech Pipette	₹4,813



+91 96262 97366
+91 73387 21875



yspscientific@gmail.com
sales.yspscientific@gmail.com



B: 25 G1, Ground Floor, #7/10
Annai Avanthika Vengaivasal,
Chennai - 600100

HEMOGLOBIN METER DIGITAL - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H324	LH113	LABOHOLIC	LH-113 Digital Haemoglobin Meter With finger touch zero setting, Memory for the storage of samples, 1 ml solution measurement.	₹10,675
H325	LH113P	LABOHOLIC	LH-113-P Portable Haemoglobin Meter Portable and high performance digital hemoglobin meter is fully automated process. This provides easy, fast & highly accurate monitoring of haemoglobin. In this the adjustment or calibration is not needed and ensures to provide results within 15seconds. 10 Testing strips, measuring range 5 - 25gm/dl, memory: 1000 test results, auto calibration.	₹9,625
H326	LHHB	LABOHOLIC	Spare HB Strips (Pack of 50)	₹4,375

We're your go-to source for all laboratory products.



Probe Sonicator



ice maker



ph meter

Upgrade your lab with top-quality equipment! Get the Probe Sonicator for efficient sample preparation, an Ice Maker for reliable cooling, and a PH Meter for accurate readings. Everything you need for precise lab work!

HUMIDITY METER - TESTO



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H327	TESTO622	TESTO	Testo 622 Large Display Hygrometer with Absolute Pressure Measures & displays % rH, °C/F, mBar,(abs), psychrometric parameters, Displays date and time, Reminder alarm for calibration, visual alarms for both the channels with Max./Min.values	POR
H328	TESTO623	TESTO	Testo 623 Large Display Hygrometer with History Function Measures & displays % rH, °C/F, mBar,(abs), psychrometric parameters, Displays date and time, Displays history of measurement for the last 90 days, visual alarms for both the channels with Max./Min.values	POR

 **YSP- SCIENTIFIC SOLUTION**
HIGH PRESSURE LIQUID CHROMATOGRAPH (HPLC)



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H329	LH3201	LABOHOLIC	LH-3201 High Pressure Liquid Chromatograph (HPLC) isocratic System with high pressure solvent delivery pump, Multi-wavelength UV Detector. Sample inject valve Rheodyne 7725i. LC Column (USA)with workstation software. Wavelength range 190-700 nm, Wavelength accuracy +1nm, spectral bandwidth 8nm, flow cell volume 8micro litre optical length-10nm. (Imported)	₹11,60,000

**volume 8micro litre optical length-
10nm. (Imported)**



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H330	GL6200	GLOLAB	<p>GL-6200 High-Pressure Liquid Chromatography system, Pump - Series dual- plunger pump, Pressure Accuracy $\pm 2\%$, Flow Rates 0.001 ~10.000 ml/min, Flow Accuracy $\leq \pm 1\%$ (2μL/min), Pressure 0 - 6000 psi, Pulsation: <15psi from 600 to 1,600 psi, Flow 1mL/min, Operating Temp.4~40°C, UV-Detector - Wavelength 190-700 nm,Spectral bandwidth 6 nm,Wavelength Accuracy $\leq \pm 1$ nm, Flow pool volume 10 μ L, Optical path 10 mm,Operating Temp. 4-40°C, Minimum detection concentration $\leq 1 \times 10^{-8}$g/mL(Naphthalene/methanol solution),Linear range $\geq 10^4$</p>	₹10,00,500
H331	GL6300	GLOLAB	<p>GL-6300 High-Pressure Liquid Chromatography system, Pump - Series dual- plunger pump, Pressure Accuracy $\pm 2\%$, Flow Rates 0.001 ~9.999 ml/min, Flow Accuracy $\leq \pm 0.5\%$, Pressure 0 - 9999 psi, Pulsation: <15psi from 600 to 1,600 psi, Flow 1mL/min, Operating Temp.4~40°C, UV-Detector - Wavelength 190-700 nm,Spectral bandwidth 6 nm,Wavelength Accuracy $\leq \pm 1$ nm, Flow pool volume 8 μ L, Optical path 10 mm,Operating Temp. 4-40°C</p>	₹15,05,000

HEMATOLOGY ANALYZER AUTOMATIC - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
H332	LH1290	LABOHOLIC	<p>LH-1290 Fully Automatic Hematology Analyzer (Double Chamber) for complete blood count Using Principle of electrical resistance for counting & SFT Method for HGB. Accurate & reliable results using digital filtering, blood cell overlap calibration dual way high pressure cleaning, to avoid aperture blocking. Electro hydraulic separation structure design alarm function, auto wiping/cleaning of sample needle. Large storage up to 30 thousand result with histogram. Easy operation using touch screen with 64Ç LCD Display. Automatic sample sipping, mixing testing & cleaning in Built thermal the printer and external printer Optional (Imported)</p>	₹5,30,700



Sterilize with confidence using the
**AUTOCLAVE VERTICAL LABORATORY
 RADIAL LOCKING - YSP.** Easy to use and
secure, perfect for any lab



INCUBATORS - YSP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
		YSP	Double walled cabinet provided with glass wool insulation with thermostat control. Temp. range 5'C above ambient to 80'C with a variation of 0.5'C complete with glass front door, expanded metal Shelves, indicator lamps rotary ON/OFF switch, to work on 220volts AC. Temp. range 5'C above room temp. to 80'C.	
11	YSPINC12	YSP	YSP-INC-12 Size Of Inner Chamber 300mm 12&x12&x12& Capacity 28 Ltr, 250 w Rating, 1 Shelves Analog Regulator Temp.	₹17,400
12	YSPINC12DIGI	YSP	YSP-INC-12 DIGI 12&x12&x12& with Digital Temperature Controller	₹26,100
13	YSPINC12DIGI FAN	YSP	YSP-INC-12 DIGI & FAN 12"x12"x12" with Digital Temperature Controller and Fan for Circulating air	₹33,350
14	YSPINC14	YSP	YSP-INC-14 Size Of Inner Chamber 355x355x355mm 14&x14&x14& Capacity 45 Ltr, 250 w Rating, 2 Shelves Analog Regulator Temp.	₹24,650
15	YSPINC14DIGI	YSP	YSP-INC-14 DIGI 14&x14&x14& with Digital Temperature Controller	₹32,625

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
16	YSPINC14DIGIFAN	YSP	YSP-INC-14 DIGI & FAN 14x14x14 with Digital Temperature Controller and Fan for Circulating air	₹40,600
17	YSPINC18	YSP	YSP-INC-18 Size Of Inner Chamber 455x455x455mm 18'x18'x18' Capacity 95 Ltr, 400 w Rating, 2 Shelves Analog Regulator Temp.	₹31,175
18	YSPINC18DIGI	YSP	YSP-INC-18 DIGI 18'x18'x18' with Digital Temperature Controller	₹39,150
19	YSPINC18DIGIFAN	YSP	YSP-INC-18 DIGI & FAN 18'x18'x18' with Digital Temperature Controller and Fan for Circulating air	₹46,400
110	YSPINC24	YSP	YSP-INC-24 Size Of Inner Chamber 605x605x605mm 24x24x24 Capacity 224 Ltr, 750 w Rating, 3 Shelves Analog Regulator Temp.	₹48,575
111	YSPINC24DIGI	YSP	YSP-INC-24 DIGI 24x24x24 with Digital Temperature Controller	₹56,550
112	YSPINC12DIGIFAN	YSP	YSP-INC-12 DIGI & FAN 24x24x24 with Digital Temperature Controller and Fan for Circulating air	₹65,250



+91 96262 97366
+91 73387 21875



yspscientific@gmail.com
sales.ypscientific@gmail.com



B: 25 G1, Ground Floor, #7/10
Annai Avanthika Vengaivasal,
Chennai - 600100

INCUBATORS - REMI



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
113	RHI50	REMI	RHI-50 Hot air Incubator with 50 Liters capacity & Temperature setting range from RT + 5°C up to +65°C	₹88,450
114	RHI80	REMI	RHI-80 Hot air Incubator with 80 Liters capacity & Temperature setting range from RT + 5°C up to +65°C	₹1,13,245

INCUBATORS - LOVIBOND



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
115	2438210	LOVIBOND	2438210 TC 140 G THERMOSTATICALLY CONTROLLED CABINET LOVIBOND	₹4,10,350

We're your go-to source for all laboratory products.



Probe Sonicator



ice maker



pH meter

Upgrade your lab with top-quality equipment! Get the Probe Sonicator for efficient sample preparation, an Ice Maker for reliable cooling, and a PH Meter for accurate readings. Everything you need for precise lab work!

ION METER & ION ELECTRODES - EUTECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
116	ECION70040S	EUTECH	ECION70040S Eutech Ion 700 Meter With Integral Electrode Holder & 100/240 VAC Adapter. pH and Ion Electrodes to be ordered separately	₹84,825
117	ECION270042GS	EUTECH	ECION270042GS Eutech ION 2700 Meter With pH Electrode (ECFG7370101B), ATC Probe (PH5TEMB01P), Integral Electrode Holder, 100/240 VAC Adapter, RS232 cable, & pH Electrode Refill Solution. Ion Electrodes to be ordered separately	POR
118	IONELECTRODES	EUTECH	ION Selective Electrodes ANTECH USA make with 500 ml Standard solution and Ion strength adjuster Price for One (FOR FLOURIDE, NITRATE, CHLORIDE, SODIUM ETC ANY ONE)	POR

INCUBATOR COOLING PLUS SERIES(Against Order With Advance Payment) - REMI



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
119	CI6PLUSLED	REMI	CI-6 Plus Capacity 200 liters, Temperature range 5°C to 60°C Microprocessor (VS-02) LED	₹2,76,225

120	CI6PLUSLCD	REMI	CI-6 Plus Capacity 200 liters, Temperature range 5°C to 60°C Microprocessor (VS-02) LCD	₹3,18,058
121	CI10PLUSLED	REMI	CI-10 Plus Capacity 280 liters, Temperature range 5°C to 60°C Microprocessor (VS-02) LED	₹2,88,623
122	CI10PLUSLCD	REMI	CI-10 Plus Capacity 280 liters, Temperature range 5°C to 60°C Microprocessor (VS-02) LCD	₹3,32,050
123	CI12PLUSLED	REMI	CI-12 Plus Capacity 370 liters, Temperature range 5°C to 60°C Microprocessor (VS-03) LED	₹3,04,138
124	CI12PLUSLCD	REMI	CI-12 Plus Capacity 370 liters, Temperature range 5°C to 60°C Microprocessor (VS-03) LCD	₹3,44,448
125	VS03	REMI	VS-03 Servo controlled Voltage Stabilizer	₹23,345
126	VS02	REMI	VS-02 Servo controlled Voltage Stabilizer	₹20,300

We're your go-to source for all laboratory products.



Probe Sonicator



ice maker



ph meter

Upgrade your lab with top-quality equipment! Get the Probe Sonicator for efficient sample preparation, an Ice Maker for reliable cooling, and a PH Meter for accurate readings. Everything you need for precise lab work!



JAR TEST APPARATUS - YSP

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
J1	YSPJAR2	YSP	YSP-JAR-2 Jar test apparatus with digital rpm & timer w/o glass beaker 2 test	₹22,475
J2	YSPJAR4	YSP	YSP-JAR-4 Jar test apparatus with digital rpm & timer w/o glass beaker 4 test	₹30,450
J3	YSPJAR6	YSP	YSP-JAR-6 Jar test apparatus with digital rpm & timer w/o glass beaker 6 test	₹34,800

JUNKER GAS CALORIMETER - YSP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
J4	YSPJUNKER	YSP	YSP-JUNKER For estimation of Calorific value of gaseous fuels with volume meter. Complete supply with all Spares & Gas flow meter.	₹203,000
J5	YSPJUNKERINSTALLATION	YSP	Installation charges for JUNKERS CALORIMETERS to be added separately with in chennai TamilNadu	₹10,500

KARL FISCHER TITRIMETER - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
K1	Kjeldhal Apparatus	YSP	Kjeldhal Apparatus Combine unit with Digestion and Distillation heating 6 Test & 6 Regulator With porcelain Coil type heater without glass parts	₹62,350
K2	Kjeldhal Apparatus glass parts 800 ml	YSP	Digestion Glass Part Full Set For 800 ml /Boss Head Clamp / Flask Clamp For Single Set	₹8,400
K3	Kjeldhal Apparatus glass parts 500ml	YSP	Digestion Glass Part Full Set For 500 ml /Boss Head Clamp / Flask Clamp For Single Set	₹7,263
K4	Kjeldhal Apparatus Distillation Glass Set of 6 with Flask for bottom side	YSP	Distillation Glass Set of 6 with Flask for bottom side	₹15,950
K5	Spare Silicon Tube 8 x11 mm 5 meter coil for Kjeldhal apparatus	YSP	Spare Silicon Tube 8 x11 mm 5 meter coil for Kjeldhal apparatus	₹1,313

We specialize in the sales of all laboratory products.

For inquiries or orders,



96262 97366
73387 21875



KARL FISCHER TITRIMETER - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
K6	LH6310	LABOHOLIC	LH-6310 Microprocessor K.F Moisture Titrator 7inch Coloured Touch screen, Microprocessor Based with Printer Interface, User Friendly and Menu Driven Software, Direct Results in ppm, % and mg H ₂ O on single display screen, Setup Data Storage Facility, Sample Data Storage Facility, 128x64 Graphical LCD, Batch-wise, Data-wise and Individual Result Data Printing, RS232/USB Interface (Optional), Measuring Range 10 ppm to 100% Moisture, Measuring method Karl fischer volumetric Titration method, Resolution of Titrant volume 0.01 ml, In built clock with date & time display, In built magnetic stirrer with adjustable speed, Dot Matrix Printer Interface. NEW ADDITION	₹117,813
K7	LH630	LABOHOLIC	LH-630 Microprocessor K.F Moisture Titrator (Graphical Display)with built-in magnetic stirrer, Auto computation of moisture in ppm, % age and mg H ₂ O, Simultaneous display of results in all three units, 1000 sample storage facility with batch no, Final result are flashed on 128*64 premium graphical LCD. Icon based Menu for easy navigation. K.F reagent dispensing by motor.	₹88,450
K8	LH63	LABOHOLIC	LH-63 Auto Karl Fishers Titrimeter Having array of LED's provided with inbuilt stirrer, auto zero burette,adjustable timer, dual platinum electrodes and other accessories.	₹28,746

LUX METER - LUTRON TESTO/ WORKZONE / HTC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
L1	LX101A	LUTRON	LX 101A LUTRON range 0 to 50000 LUX/One meter cable/Sensor with photodiode & Colour Correction Filters	₹5,950
L2	LX102	LUTRON	LX 102 LUTRON range 0 to 50000 LUX Fl. Candle. Data Hold	₹10,325
L3	LX1102	LUTRON	LX 1102 LUTRON 0 to 4 Lakh LUX Large LCD Display with RS 232 Port / cable and software to be ordered separately.	₹16,625
L4	LX1108	LUTRON	LX 1108 LUX METER with four light selection with RS 232 port / cable and software to be ordered separately.	₹16,675
L5	TESTO540	TESTO	Testo-540 Light Intensity Meter(Lux Meter) , Large backlit LCD display, Wide range upto 99,999 Lux,With Hold,Max./Min.function	₹17,255
L6	LX-101	WORKZ ONE	LX-101 Workzone LUX METER with built in sensor Range 1 - 50,000 Lux	₹3,500
L7	HTC LX-101A	HTC	HTC LX-101A 200,000 Lux Meter	₹3,500
L8	HTC LX-102A	HTC	HTC LX-102A 200,000 Lux Meter + FC	₹5,556
L9	HTC LX-103	HTC	HTC LX-103 200,000 Lux Meter + FC	₹4,506
L10	HTC LX-104	HTC	HTC LX-104 400,000 Lux Meter (Data Logging)	₹15,969
L11	HTC LX-107	HTC	HTC LX-107 Solar Power Meter	₹8,400

**LAMP ULTRA VIOLET - UV
GERMICIDAL TUBE - CLEAR**



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
L12		GENERAL	3 Feet 30 watts philips	₹2,188
L13		GENERAL	3 Feet 30 watts Osram Lamp	₹1,575
L14		GENERAL	3 Feet 30 watts Indian Fittings with & without reflector Philips	₹1,531
L15		GENERAL	1 1/2 Feet 15 watts Phillips	₹1,400
L16		GENERAL	1 1/2 Feet 15 watts Indian Fittings with & without reflector Philips	₹1,400

**LONGWAVE UV TUBES (CLEAR)
253 nm & SHORTWAVE (BLB) 360 nm**



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
L17		GENERAL	4 watts length 134.5mm (6") Phillips Clear Lamp	₹744
L18		GENERAL	4 watts length 134.5mm (6") Phillips / Osram Black Lamp	₹831
L19		GENERAL	6 watts length 210.5mm (9") Phillips / Osram Clear Lamp	₹656
L20		GENERAL	6 watts length 210.5mm (9") Phillips / Osram Black Lamp	₹770
L21		GENERAL	8 watts length 287 mm (12") Phillips Clear Lamp	₹525
L22		GENERAL	8 watts length 287 mm (12") Osram Clear Lamp	₹385

L23		GENERAL	8 watts length 287 mm (12") Phillips / Osram Black Lamp	₹700
L24		GENERAL	8 watts length 287 mm (12") OREVA FLOURESCENT TUBES	₹263
L25		GENERAL	18 watts length 2 feet Phillips / Osram Black Lamp	₹1,251
L26		GENERAL	18 watts length 2 feet Starlite	₹788
L27		GENERAL	36 watts length 4 feet Phillips / Osram Clear Lamp	₹2,100
L28		GENERAL	36 watts length 4 feet Starlite	₹1,050
L29		GENERAL	Fitting Indian 8 watts length 287 mm (12")	₹499

 **YSP- SCIENTIFIC SOLUTION**

LAMP - HALOGEN LAMP FOR MICROSCOPE/PROJECTORS/OHP & MV LAMP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
L30		GENERAL	Halogen Lamp 15 v 150 w Philips	₹735
L31		GENERAL	Halogen Lamp 15 v 150 w Osram Lamp	₹613
L32		GENERAL	Halogen Lamp 6 v 20 w Philips	₹499
L33		GENERAL	Halogen Lamp 6 v 20 w Osram Lamp	₹324
L34		GENERAL	Halogen Lamp 6 v 10 w Philips	₹359
L35		GENERAL	Halogen Lamp 6 v 10 w Osram Lamp	₹236
L36		GENERAL	Halogen Lamp 12 v 20 w Philips	₹2,363

L37		GENERAL	Halogen Lamp 12 v 20 w Osram Lamp	₹1,969
L38		GENERAL	Halogen Lamp 12 v 20 w Type 12345	₹2,363
L39		GENERAL	Halogen Lamp 24 v 250 w Philips	₹315
L40		GENERAL	Halogen Lamp 24 v 250 w Osram Lamp	₹210
L41		GENERAL	Halogen Lamp 24 v 100 w Osram Lamp	₹788
L42		GENERAL	Halogen Lamp 24 v 150 w Philips	₹315
L43		GENERAL	Halogen Lamp 24 v 150 w Osram Lamp	₹228
L44		GENERAL	Halogen Lamp 24 v 600 w	₹1,050
L45		GENERAL	Halogen Lamp 12 v 50 w Philips	₹333
L46		GENERAL	Halogen Lamp 12 v 50 w Osram Lamp	₹263
L47		GENERAL	Halogen Lamp 15 v 50 w	₹350
L48		GENERAL	Halogen Lamp 15 v 50 w with Reflector Philips	₹919
L49		GENERAL	Halogen Lamp 15 v 50 w with Reflector Osram Lamp	₹788
L50		GENERAL	Halogen Lamp 24 v 250 w with Reflector Philips	₹945
L51		GENERAL	Halogen Lamp 24 v 250 w with Reflector Osram Lamp	₹831

YSP- SCIENTIFIC SOLUTION
LAMP - MERCURY VAPOUR LAMP /
SODIUM VAPOUR LAMP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
L52		GENERAL	Mercury Vapour Lamp 125 watts	₹700
L53		GENERAL	Sodium Vapour Lamp Premium SOX-35 Watts	₹9,275

LABORATORY JACK



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
L54		GENERAL	Laboratory Jacks MS Top 15x15 cms	₹3,500
L55		GENERAL	Laboratory Jacks MS Top 12x14 cms	₹2,625
L56		GENERAL	Laboratory Jacks MS Top 20x20 cms	₹4,813
L57		GENERAL	Laboratory Jacks MS Top 25x25 cms	₹6,125
L58		GENERAL	Laboratory Jacks MS Top 30x30 cms	₹7,350
L59		GENERAL	Laboratory Jacks SS Top 12x14 cms	₹3,063
L60		GENERAL	Laboratory Jacks SS Top 20x20 cms	₹9,625
L61		GENERAL	Laboratory Jacks SS Top 25x25 cms	₹11,813
L62		GENERAL	Laboratory Jacks SS Top 30x30 cms	₹14,000
L63		GENERAL	Laboratory Jacks SS Top 15x15 cms	₹7,350

LABCOATS



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
L64		GENERAL	College school Lab coat made at Standard Quality (TC) Size S, M,L,XL any size Half Sleeve	₹613

LAMINAR AIR FLOW - HORIZONTAL



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
		YSP	Laminar Airflow system allows operation in sterile and particle free conditions because of the continuous flushing in the working area by a unidirectional and ultra filtered air flow, it assures a full product protection.	
L65	YSP 0600	YSP	YSP 0600 Horizontal Laminar Air Flow working size 2 x 2 x 2 GI sheet with powder coated 304 SS surface	₹99,325
L66	YSP 0900	YSP	YSP 0900 Horizontal Laminar Air Flow working size 3 x 2 x 2 GI sheet with powder coated 304 SS surface	₹116,725
L67	YSP 1200	YSP	YSP 1200 Horizontal Laminar Air Flow working size 4 x 2 x 2 GI sheet with powder coated 304 SS surface	₹137,460

LEAK TEST APPARATUS - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
L68	LH75	LABOHOLIC	LH-75 Digital Leak Test Apparatus (12DESICCATOR) with digit seven segment bright red LED display, programmable vacuum hold time up to 9999 secs. Extremely useful for packing industry and ideal for leak testing for strips, blisters & sachets. Please specify the desiccator size while ordering. Supplied in S.S. 304 Cabinet.	

LACTOMETER REFER "H"

MELTING POINT APPARATUS - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M1	LH108	LABOHOLIC	LH-108 Digital Melting Point Apparatus with resolution 0.1°C to measure melting point, boiling point of solid and liquid samples. Temperature range 0-300°C with 4.1" TFTScreen & Camera. NEW	₹71,050
M2	LH109	LABOHOLIC	LH-109 Automatic Melting Point Apparatus (Graphical Display) a research grade instrument with 128x64 Graphical LCD display, Icon based navigation menu for ease of operation, Programmable Ramp rates (slow, medium and Rapid). Storage up to 1000 results, Quick Cooling with Integrated fans. Audio and Visual Indications for System Status, Range ambient to 300°C, USB Interface to PC. Power 230V 210 % AC.	₹44,225
M3	LH115	LABOHOLIC	LH-115 Digital Melting Point Apparatus Within built Silicon Oil Bath and Stirrer Adjustable heating rate and stirrer speed. Range upto 275°C with readability of 1°C.	₹28,746

MELTING POINT APPARATUS - YSP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M4	YSPMPALDIE	YSP	Melting Point Apparatus	₹4,375
M5	YSPMPR	YSP	Melting Point Apparatus Superior quality	₹5,250
M6	YSPMPD	YSP	Melting Point Apparatus Digital	₹17,500

MOISTURE BALANCE IR - YSP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M7	YSPIR5	YSP	IR Moisture Balance 5 gms	₹25,375
M8	YSPIR10	YSP	IR Moisture Balance 10 gms	₹25,375
M9	YSPIR25	YSP	IR Moisture Balance 25 gms	₹25,375

MOISTURE BALANCE IR - YSP YSP- SCIENTIFIC SOLUTION



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M10		YSP	IR Moisture Tray Pkt of 80 (Spare Alluminium Foil Pan for Digital Balance)	₹8,750
M11		YSP	Spare disposable Alluminium pan for Digital IR Moisture Balance	₹123



+91 96262 97366
+91 73387 21875



yspscientific@gmail.com
sales.ypscientific@gmail.com



B: 25 G1, Ground Floor, #7/10
Annai Avanthika Vengaivasal,
Chennai - 600100

MOISTURE BALANCE - LABOHOLIC



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M12	LH76	LABOHOLIC	LH-76 Digital Moisture Balance with Weighing Range 50g, Readability 1mg, Repeatability (3g samples) 0.20%, Minimum Sample Amount 0.5g, Advised Sample Amount 3-10g, Heat up time-99 minutes, 1 minute Interval, Temperature Program Standard Terminal Control Timing, Automatic, Heating Temperature Range 50-180°C, Showing Content Moisture%, Solid%, Weight, Time, Data etc Heating Source Halogen Lamp Pan Size 100mm. Window matrix liquid crystal display (LCD) High precision heat insulation. Automatic peeling function. Halogen drying method. Highly accurate & stable. NEW	

MOISTURE BALANCE - SHIMADZU



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
UniBloc Electronic Moisture Balance				
M13	MOC63U	SHIMADZU	Model MOC63u External cali, Display 1mg/0.001g, capacity 62g, pan size 95 mm	₹246,500
M14	MOC120H	SHIMADZU	Model MOC120H External cali, Display 1mg/0.001, capacity 120g, pan size 130 mm	₹717,750



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M15	SES50MA	SAFFRON	Model SES50MA Capacity 50g, Resolution 0.001g, Pan size 90mmØ, Touch Screen, Auto external calibration	₹97,440
M16	SES110MA	SAFFRON	Model SES110MA Capacity 110g, Resolution 0.001g, Pan size 100mmØ, Touch Screen, Auto external calibration	₹90,625
M17	SESPRINTER	SAFFRON	Statistic Lab Data Printer	₹27,188
M18	DENSITYKIT SOLID	SAFFRON	Density Kit for Solid	₹18,125
M19	DENSITYKIT LIQUID	SAFFRON	Density Kit for Liquid	₹18,125
M20	WIFI21CFR	SAFFRON	WiFi Based 21 CFR with Ultimate Data saving	₹54,375



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M21	PMS713	LUTRON	Moisture Meter built in pins LUTRON Model PMS 713. Range 6 TO 40%, For Wood And 0-100% For cotton, Paper etc	₹16,625
M22	MS7000	LUTRON	MS-7000 Wood moisture meter, Range: 9% to 30% moisture content, 9 material species groups in memory, Max. min. Data hold, Self-calibration, RS-232 interface	₹17,400
M23	MS7003	LUTRON	MS-7003 Digital and LED Indicator, Range 0-40% for Wood And 0-100% Concrete and other material With Additional sensor	₹15,225

MOISTURE METER



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M24		GENERAL	Moisture Meter for Grains in bags with 12& inserting Probe China made GMM 7821 Model Range 8 - 20 %	₹7,000
M25		GENERAL	Moisture Meter for Grains in bags with 12& inserting Probe China make GMM 7002 (or) Equivalent Model Range 5 - 30 %	₹8,313
M26	WM-101	ACETEQ	Moisture Meter Digital 5 to 45% with built in pins handy China made for paper / wood etc.	₹1,750
M27	AR991	GENERAL	Grain Moisture Meter Smart Sensor Measuring range: 7.5- (rapeseed, sesame, soybean meal, cotton), 0-50 (soybean, peanut), -50 (barley, rice, sorghum, watermelon seed), -55 (paddy, wheat), .5-50 (maize, fodder).Resolution 0.1%	₹5,250

MOISTURE METER - TESTO



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M28	TESTO 606-1	TESTO	Testo 606-1 Wood/ Material Moisture Meter 8.8 to 54.8 % by weight beech, spruce, larch, birch, cherry, walnut, 7.0 to 47.9 % by weight oak, pine, maple, ash-tree, douglas fir, meranti 0.9 to 22.1 % by weight cement screed, concrete 0.0 to 11.0 % by weight anhydrite screed 0.7 to 8.6 % by weight cement mortar 0.6 to 9.9 % by weight lime mortar, plaster 0.1 to 16.5 % by weight bricks	M28

MOISTURE METER DIGITAL FOR GRAINS



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
			Digital Moisture Meter measure % moisture in the grains and immediately indicates the result on LCD display. Just pour the sample and read % moisture. Deliver 1-2 weeks after confirming the grains name to measured. Made to order for 1-8 combination of grains. Please Specify Grains like Maize, Groundnut, Paddy, Rice, Sunflower, Etc. while ordering.	₹13,650
M29		GENERAL	Moisture Meter Model DMA 1 for specific one grain	₹14,088
M30		GENERAL	Moisture Meter Model DMA 2 for specific two grains	₹14,525
M31		GENERAL	Moisture Meter Model DMA 3 for specific three grains	₹14,963
M32		GENERAL	Moisture Meter Model DMA 4 for specific four grains	₹15,400
M33		GENERAL	Moisture Meter Model DMA 5 for specific five grains	₹15,838
M34		GENERAL	Moisture Meter Model DMA 6 for specific six grains	₹16,275
M35		GENERAL	Moisture Meter Model DMA 7 for specific seven grains	₹16,713
M36		GENERAL	Moisture Meter Model DMA 8 for specific eight grains	₹13,650

MAGNETIC STIRRER WITH HOT PLATE - D LAB



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M37	MS-H280 PRO	DLAB	MS-H280-Pro Temperature up to 280°C - LED Digital Hotplate Magnetic Stirrer 5Circular Stainless steel with Ceramic Coated Hotplate	₹42,050

MAGNETIC STIRRER - YSP



M38	YSP1MLH	YSP	YSP-1MLH Magnetic Stirrer 1 litre Capacity	₹4,113
M39	YSP2MLH	YSP	YSP-2MLH Magnetic Stirrer 2 litre Capacity	₹5,250

MAGNETIC STIRRERS - REMI



		REMI	REMI ML - Digital Speed Indicator & Without Hot Plate REMI MLH -Digital Speed Indicator & With Hot Plate	
M40	MS500	REMI	Magnetic Stirrer REMI MS-500	₹6,650
M41	1ML	REMI	Magnetic Stirrer REMI 1ML(1 Litre)	₹12,775
M42	1MLH	REMI	Magnetic Stirrer REMI 1MLH	₹15,750
M43	2ML	REMI	Magnetic Stirrer REMI 2ML	₹13,125
M44	2MLH	REMI	Magnetic Stirrer REMI 2MLH(2 Litter)	₹17,150

M45	5ML	REMI	Magnetic Stirrer REMI 5ML(5 Liter)	₹16,800
M46	5MLH	REMI	Magnetic Stirrer REMI 5MLH	₹18,125

 **YSP- SCIENTIFIC SOLUTION**

MAGNETIC STIRRERS REMI PLUS SERIES



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M47	MS800PLUS	REMI	MS800 PLUS Magnetic Stirrer 800 ml,Capacity Speed 1500 RPM Remi	₹16,240
M48	5MLHPLUS	REMI	5 MLH PLUS Magnetic Stirrer with SS Hot Plate (Max.temp.340°C), LCD display of speed & temp with PT-1000 probe & Stand	₹40,600
M49	10MLHPLUS	REMI	10 MLH PLUS Magnetic Stirrer with CERAMIC Hot Plate (TEMP 550°C) LCD display of speed & temp with PT-1000 probe & stand	₹76,850
M50	20MLPLUS	REMI	20 ML Plus Magnetic stirrer 1500 RPM 20 Liter	₹49,300
M51	MPS1	REMI	Capacity	₹26,825

 **YSP- SCIENTIFIC SOLUTION**

MAGNETIC STIRRERS REMI PLUS SERIES



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M52	RQ121D	REMI	Model RQ121/D Shaft dia in mm 8X250 Speed 200 - 2000 Rpm, Stirring capacity 10 Ltr, Speed Regulator with digital RPM	₹16,240
M53	RQG121D	REMI	Model RQG121/D Shaft dia in mm 6X250, Speed 50 - 550 Rpm, Stirring capacity 5 Ltr, Speed Regulator Inbuilt with digital RPM	₹21,025
M54	RQM122R	REMI	Model RQM122/R Shaft dia in mm 6X250, Motor HP 1/20hp, Speed 800 - 4000 Rpm, Stirring capacity 7Ltr Mount On	₹15,750

M55	RQ122D	REMI	Model RQ122/D Shaft dia in mm 6X250, Motor HP 1/20hp Speed 800 - 4000 Rpm, Stirring capacity 7Ltr Remote	₹15,950
M56	RQ124AD	REMI	Model RQ 124A/D Shaft dia in mm 8X350, Motor HP 1/8HP Speed 300 - 4000 Rpm, Stirring capacity 20Ltr Remote	₹19,575
M57	RQT124AD	REMI	Model RQT124A/D Shaft dia in mm 8X350, Motor HP 1/8HP Speed 300 - 4000 Rpm, Stirring capacity 20Ltr Remote	₹26,825
M58	RQ126D	REMI	Model RQ126/D Shaft dia in mm 10X450mm, Speed 120 - 2000 Rpm, Stirring capacity 25 Ltr, Remote with digital speed	₹24,360
M59	RQ127D	REMI	Model RQ127/D Shaft dia in mm 8X240, Motor HP 1/8HP, Speed 300 - 6000 Rpm, Stirring capacity 8Ltr, Remote	₹22,910
M60	RQ127AD	REMI	Model RQ127A/D Shaft dia in mm 8X275, Motor HP 1/8HP, Speed 300 - 8000 Rpm, Stirring with 50 ml homogeniser	₹22,475
M61	RQT127AD	REMI	Model RQT127A/D Shaft dia in mm 8X275, Motor HP 1/8HP, Speed 300 - 8000 Stirring with 50 ml homogeniser	₹28,275
M62	HOMO510	REMI	5 ml or 10 ml capacity For homogeniser	₹2,625
M63	HOMO1525	REMI	15 ml or 25 ml capacity For homogeniser	₹2,800
M64	HOMO50100	REMI	50 ml or 100 ml capacity For homogeniser	₹3,325
M65	RQT127D	REMI	Model RQT127/D Medium duty high speed Emulsifier with 1/8 Motor, Digital Speed Indicator & speed regulator.	₹29,000
M66	RQ128D	REMI	Model RQ128/D Shaft dia in mm 8X310 Motor HP 1/20HP, Speed 50 - 150 Rpm Stirring capacity 3kgs, Remote	₹28,275
M67	RQ128DD	REMI	Model RQ128D/D Shaft dia in mm 12X450 Motor HP 1/2HP Speed 50 - 280 Rpm Stirring capacity 20kgs Remote	₹39,150

M68	RQG128CD	REMI	Model RQG128C/D Shaft dia in mm 10X450 Speed 50 - 500 Rpm Stirring capacity 20Lts Remote	₹28,275
M69	RQ129D	REMI	Model RQ129/D Shaft dia in mm 10X550 Speed 100 - 1400 Rpm, Stirring capacity 50Ltr Remote with digital speed	₹26,825
M70	RQG126D	REMI	Model RQG126/D Shaft dia in mm 8X450 Speed 50 - 750 Rpm Stirring capacity 15Ltr Remote with digital speed	₹29,000
M71	RQG129D	REMI	Model RQG129/D Shaft dia in mm 10X550, Speed 50 - 550 Rpm, Stirring capacity 25Ltr, Remote with digital speed	₹33,060
M72	RQ130BD	REMI	Model RQ130B/D Shaft dia in mm 12X475, Motor HP 2/3HP, Speed 300 - 5000 Rpm, Stirring capacity 35Lts, Remote	₹48,575
M73	RQ134LD	REMI	Model RQ134L/D Shaft dia in mm 12X450, Motor HP 1/2HP, Speed 900 - 5000 Rpm, Stirring capacity 60Lts, Remote	₹32,625
M74	RQ134 HD	REMI	Model RQ134 H/D Shaft dia in mm 12X450, Motor HP 2/3HP, Speed 500 - 5000 Rpm, Stirring capacity 80Lts, Remote	₹36,975
M75	RQT134HD	REMI	Model RQT134H/D Shaft dia in mm 12X450, Motor HP 2/3HP, Speed 500 - 5000 Rpm, Stirring capacity 80Lts, Remote	₹48,575
M76	RQ140D	REMI	Model RQ140/D Shaft dia in mm 10X730, Speed 100 - 1400 Rpm, Stirring capacity 40Ltr, Frequency Drive	₹84,100
M77	RQ140DFP	REMI	Model RQ140/DFP Shaft dia in mm 10X730, Speed 100 - 1400, Stirring capacity 40Ltr, Frequency Drive	₹130,500

MECHANICAL STIRRERS REMI PLUS SERIES



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M78	RQ5PLUS	REMI	RQ 5 PLUS Mini Stirrer With Burshless DC Motor, LCD Display of Speed (5 litre Capacity)	₹23,925
M79	RQ20PLUS	REMI	RQ-20 Plus Direct Drive Stirrer with brushless Motor, LCD Display of Speed & Torque (20 Letter Capacity), Speed 50 to 2200rpm	₹92,800
M80	RQ40PLUS	REMI	RQ-40 Plus Direct Drive Stirrer with brushless Motor, LCD Display of Speed & Torque (40 Letter Capacity) Speed 50 to 2200rpm	₹110,200

MICROSCOPE - LABOMED



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M81	VISION 2000 LED	LABOMED	Labomed (Old Body) VISION 2000 (Classic) Binocular Laboratory Microscope with LED Illumination	₹44,588
M82	(New Body) VISION 2000	LABOMED	Labomed (New Body) VISION 2000 Binocular Laboratory Microscope with Halogen Illumination	₹49,445
M83	LX-200 LED B	LABOMED	Labomed LX-200 LED Binocular Microscope with Battry back-up	₹51,620
M84	LX-200 LED T	LABOMED	Labomed LX-200 LED Trinocular Microscope with Battry back-up	₹58,943

M85	LX-300 LED B	LABOMED	Labomed LX-300 LED Binocular Microscope with Battery back-up with imported plan Achromatic obj.	₹64,453
M86	LX-300 LED T	LABOMED	Labomed LX-300 LED Trinocular Microscope with Battery back-up with imported plan Achromatic obj.	₹71,195
M87	LX-400 LED B	LABOMED	Labomed LX-400 LED Research Binocular Microscope with Battery back-up with infinity corrected plan Achromatic imported obj.	₹84,318
M88	LX-400 LED T	LABOMED	Labomed LX-400 LED Research Trinocular Microscope with Battery back-up with infinity corrected plan Achromatic imported obj.	₹97,440
M89	LX-500 LED B	LABOMED	Labomed LX-500 LED Research Binocular Microscope with Battery back-up with infinity corrected plan Achromatic imported obj.	₹110,345
M90	LX-500 LED T	LABOMED	Labomed LX-500 LED Research Trinocular Microscope with Battery back-up with infinity corrected plan Achromatic imported obj.	₹121,800

 **YSP- SCIENTIFIC SOLUTION**
DISSECTING MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M91	RYS-2	RAYSCOPE	RYS-2 Dissecting Monocular Microscope, heavy round cast iron base for greater stability, Eye piece pipe, knobs & lens holding arm, with black and white wide square glass stage size 85mm X 75mm having two stage clips. Magnification 10X and 20X, Illumination by Plano-concave reflector mirror fitted in fork for light reflection	₹2,144
M92	RYS-3	RAYSCOPE	RYS-3 Heavy round base providing convenience for lifting and movement. Smooth racks & pinion mechanism operated for easy working during dissection. Plano-concave mirror in fork provides illumination in day light or under artificial light.	₹3,413

STUDENT & MEDICAL MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
		RAYSCOPE	Monocular, Eye piece Huygenian 15x,10X, Objective JIS Achro 1, 4(SL), 120mmx120mm mechanical stage with clip, Separate coarse and fine movement, Condenser Abbe N A 1.25 with IRIS diaphragm, Plano concave mirror for light source, Optional LED sub stage	
M93	RYS-5	RAYSCOPE	RYS-5 Student Microscope Fixed condenser, fixed Abbe N A 1.25 with IRIS diaphragm	₹6,563
M94	RYS-7	RAYSCOPE	RYS-7 Student Microscope Movable condenser, movable Abbe N A 1.25 with IRIS diaphragm	₹7,000
M95	RYS-9	RAYSCOPE	RYS-9 Medical Microscope Movable condenser, movable Abbe N A 1.25 with IRIS diaphragm, with Blue filter & 100X	₹10,938

MONOCULAR MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M96	RYS-20M	RAYSCOPE	RYS-20 Monocular Microscope, Monocular Tube, 45° Inclined, 360° rotatable, Wide field Eye piece 10X, Objective JIS 4x, 1, 4(S) and 10 Oil(S), 140mmx140mm mechanical stage with cross travel 90(x) X 65mm (y) on ball bearing with coaxial controls, Separate coarse and fine movement, Movable abbe condenser N A 1.25 with aspheric lens, IRIS diaphragm, Illumination system LED or Halogen	₹15,750

M97	RYSINFYSTARM	RAYSCOPE	RYS-Infystar Monocular Reasearch Microscope, Monocular head 45' Inclined, 360' Rotatable, Wide field 1/18mm focusable Eyepiece with foldable eye guard, Objectives 4x,14(S), 10 Oil(S) Mechanical stage, single plate stage size 140mmx132mm X/Y Travel range 90mmx55mm, Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm, Illumination Halogen or LED	₹18,560
M98	RYS-INFINIX M LED	RAYSCOPE	RYS-Infinix Monocular LED Microscope , An Adaptable System for the Colleges & Laboratories. Sliding Monocular Head 45'Inclined, 360'rotatable, Eyepieces Wide field 1/18mm, lockable, antifungal, Semi -Plan Achromat DIN objectives 4x,1,4 (s),10(s,oil), Mechanical stage 140 x 132mm, X-Y travel range 90mm x 55mm with low-height controls, Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm with blue daylight filter, Intensity controlled LED (optional 6V-20W Halogen illumination)	₹20,119
M99	RYS-700M	RAYSCOPE	RYS-700 Research Monocular LED Microscope,360'rotatable 30'inclined Monocular head, Eyepieces High point paired WF 1/18mm, Semi-Plan Achromat DIN objectives 4x, 1, 4 (s), 10(s,oil), Mechanical stage doube plate stage size 190 x 140mm, X-Y travel range 90mm x 65mm, Double specimen SS clip, Sub stage Abbe condenser NA 1.25 with aspheric lens, Iris diaphragm with frosted white/blue filter, Illumination, halogen 6v 20w light source or LED 3w	₹21,750



★★★★★
BEST QUALITY

BINOCULAR MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M100	RYS-22B	RAYSCOPE	RYS-22 Binocular Microscope, Sliding Binocular head,45' Inclined,360' rotatable, Wide field Eye piece 10X / 20mm, Achro objectives 4x,1, 4(S) and 10 Oil(S), 190mm x140mm mechanical stage with XY travel range 90mmX65mm, Coaxial coarse and fine movement, Movable abbe condenser N A 1.25 with aspheric lens.IRIS diaphragm with Frosted white/Blue Filter, Illumination system with 6v 20W Halogen bulb or LED	₹14,863
M101	RYSINFYSTARB	RAYSCOPE	RYS-Infystar Binocular Reasearch Microscope, Binocular head 45' Inclined, 360' Rotatable, Wide field 1/18mm focusable Eyepiece with foldable eye guard, Objectives 4x,1,4(S),10 Oil(S) Mechanical stage, single plate stage size 140mmx132mm X/Y Travel range 90mmx55mm, Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm, Illumination Halogen or LED	₹19,575
M102	RYS-INFINIX-B	RAYSCOPE	RYS- Infinix Binocular LED Microscope , An Adaptable System for the Colleges & Laboratories. Sliding Binocular Head 45'Inclined, 360'rotatable, Eyepieces Wide field 1/18mm, lockable, antifungal, Semi -Plan Achromat DIN objectives 4x,1,4 (s),10(s,oil), Mechanical stage 140 x 132mm, X-Y travel range 90mm x 55mm with low-height controls, Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm with blue daylight filter, Intensity controlled LED (optional 6V-20W Halogen illumination)	₹21,750

M103	RYS-700-B	RAYSCOPE	RYS-700 Research Binocular LED Microscope, Single mould aluminium stand with hand rests for enhanced comfort and stability. Butterfly type Binocular head, 30° inclined, 360° Rotable, Eyepiece Wide field 1 / 20mm. Achromatic Objectives 4x, 1, 4(S), 10 Oil(S) Mechanical stage, double plate stage size 190mmx140mm X/Y Travel range 75mmx50mm, Double specimen holder clip, Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm with blue daylight filter, Illumination Halogen 6v-20w or 3w LED	₹26,825
M104	RYS-700 PREMIA-B	RAYSCOPE	RYS-700 Premia Binocular LED Microscope, 360° Rotable 30° sidetoff Binocular Head, Highpoint paired WF 1/28mm Eyepieces, Achromatic Objectives 4x, 1, 4(S), 10 Oil(S) Mechanical stage, double plate stage size 190x140mm X/Y Travel range 75mm x 50mm, Double specimen SS clip, Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm with frosted white/blue filter, Illumination Kohler Q-LED 3.5W 3V	₹33,350
M105	RYSICONIC55B	RAYSCOPE	RYS-ICONIC-55 Binocular Research Microscope Single mould aluminium stand with hand rests for enhanced comfort and stability, Sliding type Binocular head, 30° inclined, 360° Rotable, Eyepiece Wide field 1/18mm, DIN Semi Plan Achromatic Objectives 4x, 1, 4(S), 10 Oil(S) Mechanical stage, double plate stage size 190mmx140mm X/Y Travel range 75mmx50mm, Double specimen holder clip, Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm, Illumination Halogen or LED	₹34,438



Shop our high-precision weights E1, E2, F1, F2, & M1

Perfect for calibration and lab accuracy. Ensure reliable measurements every time



96262 97366 - 73387 21875

M106	RYS-ICONIC-77-B	RAYSCOPE	RYS-ICONIC-77 Binocular Research Microscope, Single mould aluminium stand with hand rests for enhanced comfort and stability, Butterfly type Binocular head, 30° inclined, 360° Rotable, Eyepiece Wide field 1/20mm, Achro Objectives 4x, 1,4(S), 10 Oil(S) Mechanical stage, double plate stage size 190mmx140mm X/Y Travel range 75mmx50mm, Double specimen holder clip, Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm, Illumination Halogen or LED	₹40,600
M107	RYS-ICONIC-99-B	RAYSCOPE	RYS-ICONIC-99 Binocular Research Microscope, Single mould aluminium stand with hand rests for enhanced comfort and stability, Butterfly type Binocular head, 30° inclined, 360° Rotatable, Eyepiece WF 1/20mm with anti-reflection & anti-fungal coatings. DIN Plan Achro Objectives 4x, 1,4(S), 10 Oil(S), Anti-fungal. Mechanical stage Rackless X Axis, double plate stage size 190mmx140mm X/Y Travel range 75mmx50mm, Double specimen holder clip. Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm with Blue daylight filter, Illumination Halogen 6v-20w or 3w LED	₹62,350
M108	RYS-NE-620-B	RAYSCOPE	RYS-NE-620 Binocular Research Microscope, 30° Inclined Binocular head, Eyepiece 10/22mm Infinity Plan Achro objectives 4, 10, 40(s), 100 Oil(s) or 10 water, Synchronous belt stage size 230 X 150mm, Moving range 78 X 54mm, Abbe condenser NA 1.25, Coaxial coarse & fine control, Koehler Q-LED 3.5W 3V illumination with automatic switching and sleep timer.	₹181,250



+91 96262 97366
+91 73387 21875



yspscientific@gmail.com
sales.yspscientific@gmail.com



B: 25 G1, Ground Floor, #7/10
Annai Avanthika Vengaivasal,
Chennai - 600100

TRINOCULAR MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M109	RYS-22-T	RAYSCOPE	RYS-22 Trinocular Microscope, Sliding Trinocular head,45' Inclined,360' rotatable, Wide field Eye piece 10X / 20mm, Achro objectives 4x,1, 4(S) and 10 Oil(S), 190mm x140mm mechanical stage with XY travel range 90mmX65mm.Coaxial coarse and fine movement, Movable abbe condenser N A 1.25 with aspheric lens.IRIS diaphragm with Frosted white/Blue Filter,Illumination system with 6v 20W Halogen bulb	₹15,950
M110	RYS-INFYSTAR-T	RAYSCOPE	RYS-Infystar Trinocular Research Microscope, Trinocular head 45' Inclined, 360' Rotatable. Wide field 1/18mm focusable Eyepiece with foldable eye guard.Objectives 4x,1,4(S),10 Oil(S) Mechanical stage,single plate stage size 140mmx132mm X/Y Travel range 90mmx55mm. Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm, Illumination Halogen or LED	₹21,750
M111	RYS-INFINIX-T	RAYSCOPE	RYS- Infinix Trinocular LED Microscope , An Adaptable System for the Colleges & Laboratories. Sliding Trinocular Head 45'Inclined, 360'rotatable. Eyepieces Wide field 1/18mm, lockable, antifungal. Semi -Plan Achromat DIN objectives 4x,1,4 (s),10(s,oil), Mechanical stage 140 x 132mm, X-Y travel range 90mm x 55mm with low-height controls, Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm with blue daylight filter, Intensity controlled LED (optional 6V-20W Halogen illumination)	₹23,925

M112	RYS-700-T	RAYSCOPE	RYS-700 Research Trinocular LED Microscope, Single mould aluminium stand with hand rests for enhanced comfort and stability, Butterfly type Trinocular head, 30'inclined,360' Rotable, Eyepiece Wide field 1 / 20mm.Achro Objectives 4x,1,4(S),10 Oil(S) Mechanical stage, double plate stage size 190mmx140mm X/Y Travel range 75mmx50mm, Double specimen holder clip.Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm with blue daylight filter, Illumination Halogen 6v-20w or 3w LED	₹28,275
M113	RYS-700 PREMIA-T	RAYSCOPE	RYS-700 Premia Trinocular LED Microscope,360' Rotable 30' sidentoff Trinocular Head, Highpoint paired WF 1/28mm Eyepieces, Achro Objectives 4x,1,4(S),10 Oil(S) Mechanical stage, double plate stage size 190X140mm X/Y Travel range 75mm x 50mm, Double specimen SS clip, Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm with frosted white/blue filter, Illumination Kohler Q-LED 3.5W 3V	₹34,800
M114	RYSICONIC55T	RAYSCOPE	RYS-ICONIC-55 Trinocular Research Microscope Single mould aluminium stand with hand rests for enhanced comfort and stability, Sliding type Trinocular head, 30'inclined, 360' Rotable, Eyepiece Wide field 1/18mm, DIN Semi Plan Achromatic Objectives 4x,1,4(S),10 Oil(S) Mechanical stage, double plate stage size 190mmx140mm X/Y Travel range 75mmx50mm, Double specimen holder clip, Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm, Illumination Halogen or LED	₹36,250
M115	RYS-ICONIC-77-T	RAYSCOPE	RYS-ICONIC-77 Trinocular Research Microscope, Single mould aluminium stand with hand rests for enhanced comfort and stability, Butterfly type Trinocular head, 30'inclined, 360' Rotable, Eyepiece Wide field 1/20mm, Achro Objectives 4x,1,4(S),10 Oil(S) Mechanical stage,	

			double plate stage size 190mmx140mm X/Y Travel range 75mmx50mm, Double specimen holder clip, Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm, Illumination Halogen or LED	₹43,500
M116	RYS-ICONIC-99-T	RAYSCOPE	RYS-ICONIC-99 Trinocular Research Microscope, Single mould aluminium stand with hand rests for enhanced comfort and stability, Butterfly type Trinocular head, 30' inclined, 360' Rotatable, Eyepiece WF 1/20mm with anti-reflection & anti-fungal coatings, DIN Plan Achromatic Objectives 4x, 10x, 40x(S), 100x Oil(S), Anti-fungal. Mechanical stage Rackless X Axis, double plate stage size 190mmx140mm X/Y Travel range 75mmx50mm, Double specimen holder clip, Sub stage Abbe condenser NA 1.25 with aspheric lens, IRIS diaphragm with Blue daylight filter, Illumination Halogen 6v-20w or 3w LED	₹65,250
M117	RYS-7 DELTA-T	RAYSCOPE	RYS-7 Delta Trinocular Research LED Microscope, Viewing head-Compensation Free Trinocular Head, Inclined at 30 ° 48mm-75mm, Eyepiece WF10x/22mm with rubber eye-cup, Infinite PLAN objectives 4x, 10x, 40x(s), 100x Oil, Double layers mechanical stage Double specimen holder clip, Stage Size 190mm x 160 mm, Moving Range 75mmx50mm, N.A.1.25 Abbe condenser with iris diaphragm & filter holder, Coaxial coarse & fine focusing adjustment, light Source Koehler Q-LED 3.5W, 3V illumination with variable control.	₹91,350
M118	RYS7ELITET	RAYSCOPE	RYS-7 Elite Trinocular Research Microscope, Trinocular Head Inclined at 30 °, Eyepiece WF 10x/22mm High contrast Infinity corrected PLAN Achromatic objectives 4x, 10x, 40x(s), 100x Oil, Double layers mechanical stage, Double specimen holder clip, Stage Size 216mmx140mm, X/Y Travel range 90mm X 65mm, Abbe condenser NA 1.25 with Aspheric lens, IRIS diaphragm with frosted white filter, Coaxial coarse & fine focusing, Koehler Q-LED 3.5W, 3V illumination with variable control.	₹91,350

MULTI VIEWING MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M119	RYS-NEXUS-3MV	RAYSCOPE	Nexus-3 Multi Viewing Microscope, Compensation free Trinocular Head, Inclined 30'(50mm-70mm), 2 Binocular head, Interpupillary distance adjustable 50mm-75mm, Eyepiece Plan WF 10X/22mm, Infinity Plan Objectives 4X,10X,40X,100X, X axis double layer mechanical stage size 190mm X 160mm, moving range 80mm X 50mm, Coaxial & fine focusing, Light source halogen bulb 12V/100w.	₹435,000
M120	RYS-NEXUS-5MV	RAYSCOPE	Nexus-5 Multi Viewing Microscope, Compensation free Trinocular Head. Inclined 30'(50mm-70mm), 4 Binocular head, Interpupillary distance adjustable 50mm-75mm, Eyepiece Plan WF 10X/22mm, Infinity Plan Objectives 4X,10X,40X,100X, X axis double layer mechanical stage size 190mm X 160mm, moving range 80mm X 50mm, Coaxial & fine focusing, Light source halogen bulb 12V/100w.	₹797,500
M121	RYS-NEXUS-10MV	RAYSCOPE	Nexus-10 Multi Viewing Microscope, Compensation free Trinocular Head. Inclined 30'(50mm-70mm), 9 Binocular head, Interpupillary distance adjustable 50mm-75mm, Eyepiece Plan WF 10X/22mm, Infinity Plan Objectives 4X,10X,40X,100X, X axis double layer mechanical stage size 190mm X 160mm, moving range 80mm X 50mm, Coaxial & fine focusing, Light source halogen bulb 12V/100w.	₹1,533,000

M122	RYS-PH-777MV	RAYSCOPE	<p>RYS-PH-777 Multi Viewing Microscope, Observation head one number high transmission, Trinocular head 45' inclined, 360' Rotatable, Four number compensation free Binocular head inclined 45' for multi viewing, Antifungal coated, 5 pairs of WF Eyepieces 10X/18mm with diopter adjustment, DIN Plan Achromatic Objectives 4X, 10X, 40X, 100X anti fungal coated, Double plate stage size 140 X 120, X/Y Travel range 100mm X 60mm, Abbe condenser NA 1.25, IRIS diaphragm with blue daylight filter, Coaxial & fine focusing, True Kohler Illumination, Halogen lamp 24V-100W/LED 3W with brightness controller.</p>	₹413,250
M123	RYS-PH-999MV	RAYSCOPE	<p>RYS-PH-999 Multi Viewing Microscope, Observation head one number high transmission, Trinocular head 30' inclined, 360' Rotatable. Four number compensation free Binocular head inclined 45' for multi viewing. Antifungal coated, 5 pairs of WF Eyepieces 10X/18mm with diopter adjustment, DIN Plan Achromatic Objectives 4X, 10X, 40X, 100X anti fungal coated, Double plate stage size 140 X 120, X/Y Travel range 100mm X 60mm, Abbe condenser NA 1.25, IRIS diaphragm with blue daylight filter, Coaxial & fine focusing, True Kohler Illumination, Halogen lamp 24V-100W/LED 3W with brightness controller.</p>	₹1,087,500

 **YSP- SCIENTIFIC SOLUTION**

STEREO ZOOM MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
	RYS-44S	RAYSCOPE	RYS-44 series Stereo Microscope has two Magnification levels(1/2 or 1/3 or 2/4). It is an Ideal Microscope for student learning and classroom teaching, also a perfect device for viewing Jewel, PCB, Insects and so on. WF Eyepiece 1/20mm, Binocular head 45' Inclined, 360'Rotatable, Interpupillary distance adjustment 55-75mm. 12V /10W Incandescent lamp, 12V/10W Halogen Lamp	
M124	RYS -44	RAYSCOPE	RYS -44 Ordering Information with Eyepiece 1 / 2	₹21,750
M125	RYS-44	RAYSCOPE	RYS -44 Ordering Information with Eyepiece 1/ 3	₹21,750
M126	RYS -44	RAYSCOPE	RYS -44 Ordering Information with Eyepiece 2 / 4	₹21,750
M127	RYS-45ES	RAYSCOPE	RYS-45e series Stereo Microscope can provide crystal clear sharp images with high resolution and good depth within a board field of view. Magnification 7X-45X,3.5X-270X,WF Eyepiece 1/20mm, Objective 0.7X-4.5X stereo zoom,Binocular head 45' Inclined, 360'Rotatable,Indicator LED cursor brightness adjustable.	₹50,025
M128	RYS-6745S	RAYSCOPE	RYS-6745 series Stereo Microscope is fully coated optical system with super high resolution. Ideally suitable for Industrial as well as for biological research where high precision, Maximum working distance and large field coverage is required. Eyepiece WF 15X/FN 15, WF 20X/FN 10, Binocular head 45' Inclined, 360' Rotatable, can be locked in any position desired. Interpupillary distance adjustment 54-75mm, working distance 30-165mm, Illumination transmitted & reflected 5v/3w LED, with independent switches and adjustable intensity.	₹67,425

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M129	RYS-45S	RAYSCOPE	RYS-45s series Stereo Microscope is fully coated optical system with super high resolution.Eyepiece WF 10X/22mm, Binocular head 45' Inclined, 360'Rotatable,can be locked in any position desired. Interpupillary distance adjustment 52-75mm, Objective 0.8X5. Stereo zoom, zoom ratio 6.3/Aux lens 0.5X,2.0X Illumination transmitted & reflected 5v/3w LED, with independent switches and adjustable intensity.	₹72,500
M130	RYS-DMSZ20S	RAYSCOPE	RYS-DMSZ20 series Digital Stereo Zoom Microscope integrated zoom optical module with high definition camera. It can be conveniently used in a narrow and small space. According to ergonomic principles, its internal optical corner design and angle adjustment of the display screen, allow operator to look at the display in a flat view. Through the HDMI HD interface, it can connect to another display simultaneously. Magnification 17x/1, Objective 0.7x-4.5X, Output HDMI,USB,SD card, Illumination comes with 56LED ring light.	₹152,250
M131	RYS-48A -BS	RAYSCOPE	RYS-48A High Performance parallel Optical Stereo Zoom Microscope - Binocular Head (Ask for Specifications)	₹333,500
M132	RYS-48A -TS	RAYSCOPE	RYS-48A High Performance parallel Optical Stereo Zoom Microscope - Trinocular Head (Ask for Specifications)	₹435,000
M133	RYS-SZX12S	RAYSCOPE	RYS-SZX12 Research Stereo Zoom Microscope (Ask for Specifications)	₹1,015,000
M134	RYS2D3DS	RAYSCOPE	RYS-2D/3D Inspection Microscope (Ask for Specifications)	₹176,900



YSP- SCIENTIFIC SOLUTION

RESEARCH INVERTED SYSTEM MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M135	RYS60RISM	RAYSCOPE	RYS-60 Research Inverted Sytem Microscope, viewing head 20-45 tilting binocular head, inverted image, interpupillary distance 50-76mm, High eye point plan eyepiece PL10X/22mm, diopter adjustable,reticle attachable	POR



YSP- SCIENTIFIC SOLUTION

MICROTOME - RAYSCOPE



M136	RYS-ARM-10B	RAYSCOPE	RYS-ARM-10B Erma Japan Type Rotary Microtome Regular Quality	₹44,950
M137	RYS-ARM-10A	RAYSCOPE	RYS-ARM-10A Spencer Type Rotary Microtome High Precision Quality	₹68,875



YSP- SCIENTIFIC SOLUTION

BIOLOGICAL RESEARCH MICROSCOPE - RAYSCOPE



M138	RYS-ICX-40	RAYSCOPE	RYS-ICX-40 Biological Research Microscope(Ask for Specifications)	POR
M139	RYS-IRX-50BIO	RAYSCOPE	RYS-IRX-50 Biological Research Microscope(Ask for Specifications)	POR

FLUORESCENT RESEARCH MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
		RAYSCOPE	RYS-Delta-7FLB LED Fluorescent Research Microscope (Ask for Specifications)	
M140	RYS7FLBBLUE	RAYSCOPE	Single Filter (Blue)	₹181,250
M141	RYS7FLBG	RAYSCOPE	Dual Filters(Blue & Green)	₹326,250
		RAYSCOPE	RYS-102 Fluorescent Research Microscope with Two Filter (Ask for Specifications)	
M142	RYS-102BG	RAYSCOPE	2 Filters (Blue & Green)	P.O.R
M143	RYS102BGUV	RAYSCOPE	4 Filters (Blue & Green , Ultra Violet and Violet)	P.O.R
M144	RYSSTAR9FL	RAYSCOPE	RYS-Star-9FL Research Biological Fluorescent Microscope (5+1 Filters) (Blue & Green , Ultra Violet and Violet , Blue Violet + One Empty Block) (Ask for Specifications)	₹1,135,050
		RAYSCOPE	RYS-106-FL Inverted Fluorescent Research Microscope (Ask for Specifications)	
M145	RYS-106-FLBG	RAYSCOPE	2 Filters (Blue & Green)	P.O.R
M146	RYS-106-FLFLBGUV	RAYSCOPE	4 Filters (Blue & Green , Ultra Violet and Violet)	P.O.R
M147	RYS-600-FL	RAYSCOPE	RYS-600-FL Research Biological Fluorescent Microscope (Ask for Specifications)	POR
M148	RYS-60-FL	RAYSCOPE	RYS-60-FL Inverted Fluorescent Research Microscope (Ask for Specifications)	₹978,750



YSP- SCIENTIFIC SOLUTION

BIOLOGICAL RESEARCH MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M149	RYS-105STISSUE	RAYSCOPE	RYS-105s Research Tissue Culture Microscope (Ask for specifications)	₹337,125
M150	RYS-41TISSUE	RAYSCOPE	RYS-41 Cell Culture Research Microscope (Ask for specifications)	POR



YSP- SCIENTIFIC SOLUTION

POLARIZING MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M151	RYS-PR-1	RAYSCOPE	RYS-PR-1 Monocular Polarizing Microscope (Ask for specifications)	₹43,500
M152	RYS-PR-2	RAYSCOPE	RYS-PR-2 Binocular Polarizing Microscope (Ask for specifications)	₹46,400
M153	RYS-PR-3	RAYSCOPE	RYS-PR-3 Trinocular Polarizing Microscope (Ask for specifications)	₹47,850
M154	RYS-PR-5	RAYSCOPE	RYS-PR-5 Research Polarizing Microscope (Ask for specifications)	₹60,900
M155	RYS-SMARTPOL	RAYSCOPE	RYS-Smart POL Research Polarizing Microscope (Ask for specifications)	₹224,025



YSP- SCIENTIFIC SOLUTION

METALLURGICAL MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M156	RYSITMM	RAYSCOPE	RYS-INFINIX Trinocular Metallurgical Microscope	₹79,750
M157	RYS-M40	RAYSCOPE	RYS-M40 Inverted Metal Analysis Microscope	₹152,250
M158	RYS-M50	RAYSCOPE	RYS-M50 Inverted Metallurgical Microscope	₹217,500
M159	RYS-M200	RAYSCOPE	RYS-M200 Inverted Metallurgical Microscope	₹195,750
M160	RYS-M500	RAYSCOPE	RYS-M500 Upright Metallurgical Microscope	₹253,750

 **YSP- SCIENTIFIC SOLUTION**

INVERTED DARK & BRIGHT FIELD MICROSCOPE - RAYSCOPE



RYS-20/20BD Inverted Dark and Bright Field Microscope

M161	RYS-20BF	RAYSCOPE	Bright Field	₹275,500
M162	RYS-20DF	RAYSCOPE	BF/DF	₹693,100
M163	RYS-20DIC	RAYSCOPE	RYS-20-DIC Inverted Dark and Bright Field Microscope	POR

 **YSP- SCIENTIFIC SOLUTION**

TOOL MAKERS MICROSCOPE - RAYSCOPE



M164	RYS-TM-200	RAYSCOPE	RYS-TM-200 Tool Makers Microscope	₹63,800
------	------------	----------	-----------------------------------	---------

TRAVELLING MICROSCOPE - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M165	RYS-TM2	RAYSCOPE	RYS-TM2 Two Motion Travelling Microscope	₹11,375
M166	RYS-TM4	RAYSCOPE	RYS-TM4 Three Motion Travelling Microscope	₹14,000
M167	RYS-TM6	RAYSCOPE	RYS-TM6 Digital Travelling Microscope	POR

GEMOLOGICAL MICROSCOPE - RAYSCOPE



M168	RYS-GM50	RAYSCOPE	RYS-GM50 Gemological Binocular Microscope	₹478,500
M169	RYS-GM52	RAYSCOPE	RYS-GM52 Gemological Trinocular Microscope	₹478,500

SURGICAL MICROSCOPE - RAYSCOPE



M170	RYSOPM400SURGICAL	RAYSCOPE	RYS-OPM-400 Surgical Operating Microscope	₹478,500
M171	RYS-OPM-400ENT	RAYSCOPE	RYS-OPM-400 Surgical ENT Microscope	₹478,500
M172	RYS-OPM-400DENTAL	RAYSCOPE	RYS-OPM-400 Dental Operating Microscope	₹478,500

UA SERIES CAMERA - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M173	C51MP	RAYSCOPE	UA Series 5.1 MP Camera	₹35,525
M174	C10MP	RAYSCOPE	UA Series 10.0 MP Camera	₹50,750
M175	C63MP	RAYSCOPE	1080P HDMI Camera 6.3 MP	₹100,050
M176	C8MP4K	RAYSCOPE	8MP 4K HDMI + USB CAMERA	₹137,750
M177	C5MPTOUCH	RAYSCOPE	RYS Touch Digital touch Screen with built in android tab and 5.0 MP Camera	₹104,219

MICROSCOPE CAMERA - LUCIDA



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M178		GENERAL	Microscope Spare Camera Lucida Mirror type	₹1,400
M179		GENERAL	Microscope Spare Camera Lucida Prism type	₹1,225

MD SERIES CAMERA - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M180	MD130	RAYSCOPE	MD 130 , 1.3 MP, Max. Resolution 1280*1024, Pixel size 2.2µm*2.2µm, Spectral range 400-650nm	₹18,488
M181	MD300	RAYSCOPE	MD 300 , 3.0 MP, Max. Resolution 2048*1536, Pixel size 2.2µm*2.2µm, Spectral range 400-650nm	₹20,663
M182	MD500	RAYSCOPE	MD 500 , 5.0 MP, Max. Resolution 2592*1944, Pixel size 2.2µm*2.2µm, Spectral range 400-650nm	₹23,925

MICROSCOPE SPARES - RAYSCOPE (AGAINST ORDER ONLY 7 DAYS TIME)



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
MECHANICAL STAGE				
M183	MSCSP	RAYSCOPE	Coaxial Single Plate Mechanical Stage	₹2,100
M184	MSCXM	RAYSCOPE	Coaxial X-axis rackless Mechanical Stage	₹2,625
HEAD				
M185	HEAD45B	RAYSCOPE	45 degree Bino Head Sliding Type	₹2,888
M186	HEAD30SB	RAYSCOPE	30 degree Sliding Head Bino	₹4,375

M187	HEAD30T	RAYSCOPE	30 degree Trino Head	₹5,425
M188	HEAD30BB	RAYSCOPE	30 Degree Siedentopf Binocular Head butterfly type	₹9,800
M189	HEAD30TB	RAYSCOPE	30 Degree Siedentopf Binocular Head butterfly type	₹12,950

LED Circuit

M190	LEDVC	RAYSCOPE	LED Circuit with Volume Control	₹2,450
M191	LEDBB	RAYSCOPE	LED Circuit With Volume Control & Battery Backup	₹4,813

Phase Contrast attachment

M192	PCAST	RAYSCOPE	Sliding Type phase contrast attachment with 10x , 40x phase objectives and centering Eyepiece and Green Filter	₹30,000
M193	PCATT	RAYSCOPE	Turret Type phase contrast attachment with 10x 20x, 40x, 100xphase objectives andcentering Eyepiece and Green Filter	₹35,000
M194	PCAIC	RAYSCOPE	Infinity Corrected Turret Type phase contrast attachment with 10x 20x,40x , 100xphase objectives andcentering Eyepiece and Green Filter	₹106,000

M195		RAYSCOPE	Spare Objective Superior 4X	₹1,000
M196		RAYSCOPE	Spare Objective Superior 10X	₹1,500
M197		RAYSCOPE	Spare Objective Superior 40X	₹2,200
M198		RAYSCOPE	Spare Objective Superior 100X (S) Oil	₹3,000

Objective Rayscope Quality Suitable for Our All Regular Bino / Trino Ranges

M199		RAYSCOPE	Spare Objective Rayscope 4X	₹875
M200		RAYSCOPE	Spare Objective Rayscope 10X	₹1,094
M201		RAYSCOPE	Spare Objective Rayscope 40X	₹1,400
M202		RAYSCOPE	Spare Objective Rayscope 100X (S) Oil	₹2,625

Objective Economic Quality Suitable for Student / Medical Microscope Rayscope

M203		RAYSCOPE	Spare Objective Student / Medical 4X	₹525
M204		RAYSCOPE	Spare Objective Student / Medicalc 10X	₹788
M205		RAYSCOPE	Spare Objective Student / Medical 40X	₹963
M206		RAYSCOPE	Spare Objective Student / Medical 100X (S) Oil	₹1,488

Eye Piece Rayscope

M207		RAYSCOPE	Eyepiece Superior Quality 10X / 18 mm WF	₹1,094
M208		RAYSCOPE	Eyepiece Ecomonic Quality 10X / 18 mm WF	₹875
M209		RAYSCOPE	Eyepiece Student Model 5X HYG	₹613
M210		RAYSCOPE	Eyepiece Student Model 10X HYG	₹525
M211		RAYSCOPE	Eyepiece Student Model 15X HYG	₹525

Dark Field Objective

M212	DFO100	RAYSCOPE	100x (S) Oil iris diaphragm dark field obejective with condenser	₹18,125
------	--------	----------	---------------------------------------------------------------------	---------

 **YSP- SCIENTIFIC SOLUTION**

**MICROSCOPE MAGNUS MFG BY OLYMPUS
(Other Model Against Request) - OLYMPUS**



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M213	HM100ECO	OLYMPUS	MagMaster Monocular HM 100 Eco Teaching Microscope All Metal rugged Stand with Individual coarse and fine Eco LED Illumination including battery backup 110 x120 mm stand 4x,1,4 Acromatic Objects 360' rotatable 45' head, with 1	₹24,962
M214	SM100ECO	OLYMPUS	MagMaster Monocular SM 100 Eco Teaching Microscope All Metal rugged Stand with Individual coarse and fine Eco LED Illumination including battery backup 115 x125 mm stand Mechanical stage 4x,1,4d 100 x Oil Imm. Acromatic Objects 360' rotatable 45' head, with WF 1	₹33,749
M215	MLXMPLUSPLE D	OLYMPUS	Model MLX-M PLUS (SP) co-axial coarse and fine focusing controls,high resolution long barrel Semi-plan achromatic objectives 4x,1,4 and 10 with an inclined Monocular 360 degree fitted with fungus resistant, anti-reflection coated prism,widefield paired eyepiece WF1 co-axial low drive, mechanical stage, LED Light	₹45,682
M216	MLXBPLUSSP	OLYMPUS	Model MLX-B PLUS (SP) co-axial coarse and fine focusing controls, high resolution long barrel Semi-plan achromatic objectives 4x,1,4 and 10 with an inclined binocular 360 degree fitted with fungus resistant, anti-reflection coated prism, widefield paired eyepiece WF1 co-axial low drive, mechanical stage, Halogen light	₹45,305
M217	MLXPLUSPLE D	OLYMPUS	Model MLX-B PLUS (SP) co-axial coarse and fine focusing control, high resolution long barrel Semi-plan achromatic obj. LED	₹49,249
M218	MLXTRPLUSSP	OLYMPUS	Model MLX-TR PLUS (SP) co-axial coarse and fine focusing controls, high resolution long barrel Semi-plan achromatic objectives 4x,1,4 and 10 with an inclined Trinocular 360 degree fitted with fungus resistant, anti-reflection coated prism, widefield paired eyepiece WF1 co-axial low drive, mechanical stage Halogen lamp	₹49,017

SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M219	LXTRPLUSSPLED	OLYMPUS	SAME AS ABOVE BUT WITH LED LIGHT MLX-TR PLUS (SP) LED	₹65,707
M220	CH20I	OLYMPUS	Model CH20i (BINOCULAR VERSION) in-built 6V 20W halogen light illuminator, with Anti-Fungus Achromatic Objectives, 1,4 (spring) & 10 (oil, spring), paired eyepieces wide field 1 HALOGEN LIGHT	₹67,338
M221	CH20ILED	OLYMPUS	Model CH20i (BINOCULAR VERSION) , with Anti-Fungus Achromatic Objectives, 1,4 (spring) & 10 (oil, spring), paired eyepieces wide field 1 LED LIGHT	₹69,593
M222	CH20ITR	OLYMPUS	MODLE CH20i TR (Trinocular Version) in bilt 6V 20W halogen light illuminator, with Anti-Fungus Plan Achromatic Objectives 4x, 1,4 (spring) & 10 (oil spring), paired eyepieces wide field 1 HALOGEN LIGHT	₹71,239
M223	CH20ITRLED	OLYMPUS	MODLE CH20i TR (Trinocular Version), with Anti-Fungus Plan Achromatic Objectives 4x, 1,4 (spring) & 10 (oil spring), paired eyepieces wide field 1 LED LIGHT	₹91,894
M224	MX21IB	OLYMPUS	Model MX21i (BINOCULAR VERSION) in-built 6V 20W halogen light illuminator, with Anti-Fungus Plan Achromatic Objectives 4x, 1,4 (spring) & 10 (oil spring), paired eyepieces wide field 1 Halogen lamp	₹94,170
M225	MX21IBLED	OLYMPUS	Model MX21i (BINOCULAR VERSION) LED Display illuminator with Anti-Fungus Plan Achromatic Objectives 4x, 1,4 (spring) & 10(oil spring) paired eyepieces wide field 1 LED lamp Illuminated	₹96,063
M226	MX21T	OLYMPUS	Model MX21i (Trinocular VERSION) in-built 6V 20W halogen light illuminator, with Anti-Fungus Plan Achromatic Objectives 4x, 1,4 (spring) & 10 (oil spring), paired eyepieces wide field 1 Halogen lamp	₹96,063

M227

MX21TLED

OLYMPUS

Model MX21i (Trinocular VERSION), with Anti-Fungus Plan Achromatic Objectives 4x, 1, 4 (spring) & 10 (oil spring), paired eyepieces wide field 1 LED lamp

₹98,245

 **YSP- SCIENTIFIC SOLUTION**

MICROSCOPE STEREO MAGNUS - OLYMPUS



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M228	MSZBI	OLYMPUS	Stereo Zoom MAGNUS Model MSZ BI , 45' degree Inclined Ratio 1:7, Obj Zoom Range 0.65x - 4.5x Hologen Lamp	₹99,238
M229	MSZBILED	OLYMPUS	Stereo Zoom MAGNUS Model MSZ BI , 45' degree Inclined Ratio 1:7, Obj Zoom Range 0.65x - 4.5x LED Light	₹100,130
M230	MSZTR	OLYMPUS	Stereo Zoom MAGNUS Model MSZ TR , 45' degree Inclined Ratio 1:7, Obj Zoom Range 0.65x - 4.5x Hologen Lamp	₹141,977
M231	MSZTRLED	OLYMPUS	Stereo Zoom MAGNUS Model MSZ TR , 45' degree Inclined Ratio 1:7, Obj Zoom Range 0.65x - 4.5x LED Lamp	₹150,619

 **YSP- SCIENTIFIC SOLUTION**

MICROSCOPE STEREO MAGNUS - OLYMPUS YSP- SCIENTIFIC SOLUTION




SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M232	MLXBPLUS	OLYMPUS	Phase Contrast attachment only for Olympus Model MLXB Plus or CH 20i	POR

MANOMETER (SUPPLIED WITHOUT MERCURY (OR) LIQUID) - LUTR



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M233	PM9107	LUTRON	Manometer digital PM9107 Lutron, Taiwan Range 7 mbar	POR
M234	TESTO510	TESTO	Testo-510 Differential Pressure & Flow Meter (Digital Manometer) , 8 pressure units selectable e.g. mbar, Pa, mm, H ₂ O, etc., Temperature compensated, can be connected with pitot tube for velocity measurement	₹19,104

MAGNIFIER - HAND TYPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M235		GENERAL	Torch Magnifier 75mm (3")	₹656
M236		GENERAL	Torch Magnifier 100mm (4")	₹744
M237		GENERAL	Torch Magnifier 138mm (5.5")	₹1,050

MAGNIFIER WITH STAND



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M238	L104 ST	GENERAL	TABLE STAND Model L104 ST with extra magnification attachment for 10X Lensel make With Flexible Arm Model STD	₹10,500
M239	L104 ST (LED)	GENERAL	TABLE STAND Model L104 ST (LED) with extra magnification attachment for 10X Lensel make With Flexible Arm Model STD	₹11,813

YSP- SCIENTIFIC SOLUTION
MICROPIPETTE FIXED VOLUME
(ECONOMY MODEL) - GLOLINE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M240	GLFV1	GLOLINE	Gloline Range 1 ul	₹700
M241	GLFV2	GLOLINE	Gloline Range 2 ul	₹770
M242	GLFV5	GLOLINE	Gloline Range 5 ul	₹770
M243	GLFV10	GLOLINE	Gloline Range 10 ul	₹770
M244	GLFV20	GLOLINE	Gloline Range 20 ul	₹770
M245	GLFV50	GLOLINE	Gloline Range 50 ul	₹700
M246	GLFV100	GLOLINE	Gloline Range 100 ul	₹700
M247	GLFV200	GLOLINE	Gloline Range 200 ul	₹700
M248	GLFV250	GLOLINE	Gloline Range 250ul	₹700
M249	GLFV500	GLOLINE	Gloline Range 500 ul	₹700
M250	GLFV1000	GLOLINE	Gloline Range 1000 ul	₹700
M251	GLFV2000	GLOLINE	Gloline Range 2000 ul	₹700
M252	GLFV5000	GLOLINE	Gloline Range 5000 ul	₹700
M253	GLFV10000	GLOLINE	Gloline Range 10000 ul	₹700

YSP- SCIENTIFIC SOLUTION
MICROPIPETTE VARIABLE VOLUME
(ECONOMY MODEL) - GLOLINE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M254	GLVV22	GLOLINE	Gloline Range 0.2 -2 ul	₹963
M255	GLVV510	GLOLINE	Gloline Range 0.5 -10 ul	₹963
M256	GLVV220	GLOLINE	Gloline Range 2 -20 ul	₹963
M257	GLVV550	GLOLINE	Gloline Range 5 -50 ul	₹963
M258	GLVV10100	GLOLINE	Gloline Range 10 -100 ul	₹963
M259	GLVV20200	GLOLINE	Gloline Range 20 -200 ul	₹963
M260	GLVV1001000	GLOLINE	Gloline Range 100 -1000 ul	₹963
M261	GLVV55	GLOLINE	Gloline Range 0.5 -5 ml	₹963
M262	GLVV110	GLOLINE	Gloline Range 1-10 ml	₹963

YSP- SCIENTIFIC SOLUTION
MICROPIPETTE PRIME FIXED VOLUME -
GLOLINE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M263	GLPFV1	GLOLINE	Gloline Prime Range 1 ul	₹1,313
M264	GLPFV2	GLOLINE	Gloline Prime Range 2 ul	₹1,313
M265	GLPFV5	GLOLINE	Gloline Prime Range 5 ul	₹1,313
M266	GLPFV10	GLOLINE	Gloline Prime Range 10 ul	₹1,313
M267	GLPFV20	GLOLINE	Gloline Prime Range 20 ul	₹1,313
M268	GLPFV25	GLOLINE	Gloline Prime Range 25 ul	₹1,313
M269	GLPFV50	GLOLINE	Gloline Prime Range 50 ul	₹1,313
M270	GLPFV100	GLOLINE	Gloline Prime Range 100 ul	₹1,313
M271	GLPFV200	GLOLINE	Gloline Prime Range 200 ul	₹1,313

M272	GLPFV250	GLOLINE	GloLine Prime Range 250ul	₹1,313
M273	GLPFV500	GLOLINE	GloLine Prime Range 500 ul	₹1,313
M274	GLPFV1000	GLOLINE	GloLine Prime Range 1000 ul	₹1,313

 **YSP- SCIENTIFIC SOLUTION**

MICROPIPETTE PRIME VARIABLE VOLUME - GLOLINE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M282	FAV2	MICROLIT	FAV- 2 Range 0.2 -2 ul Microlit	₹3,063
M283	FAV10	MICROLIT	FAV- 10 Range 0.5 -10 ul Microlit	₹3,063
M284	FAV20	MICROLIT	FAV- 20 Range 2 -20 ul Microlit	₹3,063
M285	FAV50	MICROLIT	FAV- 50 Range 5 -50 ul Microlit	₹3,063
M286	FAV100	MICROLIT	FAV- 100 Range 10 -100 ul Microlit	₹3,063
M287	FAV200	MICROLIT	FAV- 200 Range 20 -200 ul Microlit	₹3,063
	FAV1000	MICROLIT	FAV- 1000 Range 100 -1000 ul Microlit (NO STOCK - PRODUCT DISCONTINUED)	
M288	FAV5000	MICROLIT	FAV- 5000 Range 0.5 -5 ml Microlit	₹3,063
M289	FAV10000	MICROLIT	FAV- 10000 Range 1-10 ml Microlit	₹3,063

 **YSP- SCIENTIFIC SOLUTION**

MICROPIPETTE FIXED VOLUME- FAF SERIES (ECONOMY MODEL) (TILL STOCK LAST)- MICROLIT



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M290	FAF1	MICROLIT	FAF- 1 Range 1 ul Microlit	₹1,969
M291	FAF2	MICROLIT	FAF- 2 Range 2 ul Microlit	₹1,969
M292	FAF5	MICROLIT	FAF- 5 Range 5 ul Microlit	₹1,969
M293	FAF10	MICROLIT	FAF- 10 Range 10 ul Microlit	₹1,969
M294	FAF20	MICROLIT	FAF- 20 Range 20 ul Microlit	₹1,969
M295	FAF25	MICROLIT	FAF- 25 Range 25 ul Microlit	₹1,969
M296	FAF50	MICROLIT	FAF- 50 Range 50 ul Microlit	₹1,969
M297	FAF100	MICROLIT	FAF- 100 Range 100 ul Microlit	₹1,969
M298	FAF200	MICROLIT	FAF- 200 Range 200 ul Microlit	₹1,969
M299	FAF250	MICROLIT	FAF- 250 Range 250ul Microlit	₹1,969
M300	FAF500	MICROLIT	FAF- 500 Range 500 ul Microlit	₹1,969
M301	FAF1000	MICROLIT	FAF- 1000 Range 1000 ul Microlit	₹1,969
M302	FAF2000	MICROLIT	FAF- 2000 Range 2000 ul Microlit	₹1,969
M303	FAF5000	MICROLIT	FAF- 5000 Range 5000 ul Microlit	₹1,969
M304	FAF10000	MICROLIT	FAF- 10000 Range 10000 ul Microlit	₹1,969

 **YSP- SCIENTIFIC SOLUTION**

MICROPIPETTES MICROLIT FULLY AUTOCLAVABLE VARIABLE VOLUME- SLP SERIES



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M305	SLP-2	MICROLIT	SLP-2 Range 0.2 -2 ul Microlit	₹6,125
M306	SLP- 10	MICROLIT	SLP- 10 Range 0.5 -10 ul Microlit	₹6,125
M307	SLP-20	MICROLIT	SLP- 20 Range 2 -20 ul Microlit	₹6,125
M308	SLP-50	MICROLIT	SLP- 50 Range 5 -50 ul Microlit	₹6,125
M309	SLP-100	MICROLIT	SLP- 100 Range 10 -100 ul Microlit	₹6,125
M310	SLP-200	MICROLIT	SLP- 200 Range 20 -200 ul Microlit	₹6,125
M311	SLP-1000	MICROLIT	SLP- 1000 Range 100 -1000 ul Microlit	₹6,125
M312	SLP-5000	MICROLIT	SLP- 5000 Range 0.5 -5 ml Microlit	₹6,125
M313	SLP-10000	MICROLIT	SLP- 10000 Range 1-10 ml Microlit	₹6,125

YSP- SCIENTIFIC SOLUTION

MICROPIPETTES MICROLIT FULLY AUTOCLAVABLE FIXED VOLUME- SLP SERIES



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M314	SLP-F-1	MICROLIT	SLP-F-1 Range 1 ul Microlit	₹3,850
M315	SLP-F-2	MICROLIT	SLP-F-2 Range 2 ul Microlit	₹3,850
M316	SLP-F-5	MICROLIT	SLP-F- 5 Range 5 ul Microlit 	₹3,850
M317	SLP-F-10	MICROLIT	SLP-F-10 Range 10 ul Microlit	₹3,850
M318	SLP-F-20	MICROLIT	SLP-F-20 Range 20 ul Microlit	₹3,850
M319	SLP-F-25	MICROLIT	SLP-F-25 Range 25 ul Microlit	₹3,850
M320	SLP-F-50	MICROLIT	SLP-F-50 Range 50 ul Microlit	₹3,850
M321	SLP-F-100	MICROLIT	SLP-F-100 Range 100 ul Microlit	₹3,850
M322	SLP-F-200	MICROLIT	SLP-F- 200 Range 200 ul Microlit	₹3,850

M323	SLP-F-250	MICROLIT	SLP-F-250 Range 250ul Microlit	₹3,850
M324	SLP-F-500	MICROLIT	SLP-F-500 Range 500 ul Microlit	₹3,850
M325	SLP-F-1000	MICROLIT	SLP-F-1000 Range 1000 ul Microlit	₹3,850
M326	SLP-F-2000	MICROLIT	SLP-F- 2000 Range 2000 ul Microlit	₹3,850
M327	SLP-F-5000	MICROLIT	SLP-F-5000 Range 5000 ul Microlit	₹3,850
M328	SLP-F-10000	MICROLIT	SLP-F-10000 Range 10000 ul Microlit	₹3,850


YSP- SCIENTIFIC SOLUTION
MICROPIPETTES MICROLIT FULLY
AUTOCLAVABLE FIXED VOLUME- SLP SERIES
YSP- SCIENTIFIC SOLUTION



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M329	RBO2	MICROLIT	RBO - 2 Range 0.2-2 ul Microlit Fully Autoclavable	₹9,013
M330	RBO10	MICROLIT	RBO - 10 Range 0.5-10 ul Microlit Fully Autoclavable	₹9,013
M331	RBO20	MICROLIT	RBO - 20 Range 2-20 ul Microlit Fully Autoclavable	₹9,013
M332	RBO50	MICROLIT	RBO - 50 Range 5-50 ul Microlit Fully Autoclavable	₹9,013
M333	RBO100	MICROLIT	RBO - 100 Range 10-100 ul Microlit Fully Autoclavable	₹9,013
M334	RBO200	MICROLIT	RBO - 200 Range 20-200 ul Microlit Fully Autoclavable	₹9,013
M335	RBO1000	MICROLIT	RBO - 1000 Range 100-1000 ul Microlit Fully Autoclavable	₹9,013
M336	RBO5000	MICROLIT	RBO - 5000 Range 0.5-5 ml Microlit Fully Autoclavable	₹9,013
M337	RBO10000	MICROLIT	RBO - 10000 Range 1-10 mll Microlit Fully Autoclavable	₹9,013

YSP- SCIENTIFIC SOLUTION

**MICROPIPETTES MICROLIT FULLY
AUTOCLAVABLE FIXED VOLUME- SLP SERIES
YSP- SCIENTIFIC SOLUTION**



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M338	GLMC850	GLOLINE	MC-8/50 Range 5-50 ul Gline	₹16,459
M339	GLMC8300	GLOLINE	MC-8/300 Range 30-300 ul Gline	₹16,459

YSP- SCIENTIFIC SOLUTION

**MICROPIPETTES MICROLIT FULLY AUTOCLAVABLE
FIXED VOLUME- SLP SERIES YSP- SCIENTIFIC
SOLUTION**



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M340	MCA810	MICROLIT	MCA-8/10 Range 0.5-10 ul Microlit	₹16,501
M341	MCA820	MICROLIT	MCA-8/20 Range 2-20 ul Microlit	₹16,501
M342	MCA850	MICROLIT	MCA-8/50 Range 5-50 ul Microlit	₹16,501
M343	MCA8100	MICROLIT	MCA-8/100 Range 10-100 ul Microlit	₹16,501
M344	MCA8200	MICROLIT	MCA-8/200 Range 20-200 ul Microlit	₹16,501
M345	MCA8300	MICROLIT	MCA-8/300 Range 40-300 ul Microlit	₹16,501

YSP- SCIENTIFIC SOLUTION

**MOTORIZED BURETTE - e- BURETTE
Microlit (New Model)**



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M346	TI10	MICROLIT	TI-10 Range 10 ml Microlit	₹127,600
M347	TI25	MICROLIT	TI-25 Range 25 ml Microlit	₹127,600
M348	TI50	MICROLIT	TI-50 Range 50 ml Microlit	₹127,600
M349	TISWS	MICROLIT	MICROLIT E-BURETTE SOFTWARE 21 CFR Compliant (for Titration, Data Capturing & Calculations) Server set up & one client license (TI-SW-S)	₹79,750
M350	TISWC	MICROLIT	TI-SW-C Additional client license (required for each additional computer setup)	₹43,500

 YSP- SCIENTIFIC SOLUTION

TIPS FOR MICROPIPETTES - MICROLIT



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M351	FILTER TIPS WITH BOX OF 96 TIPS	MICROLIT	Microtips 10ML / 10000ul (PACK OF 100) MICROLIT 	₹2,625
M352		MICROLIT	Microtips 5ML / 5000ul(PACK OF 100) MICROLIT	₹1,750

 YSP- SCIENTIFIC SOLUTION

TIPS FOR MICROPIPETTES / MCP TUBE / STORAGE VIAL



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M353		GENERAL	Tips Vol.Range200μl Pack of 1000 Gloine / Sp / Vp colour Yellow or White	₹438
M354		GENERAL	Tips Vol.Range1000μl Pack of 500 Gloine / Sp / Vp colour Blue or White	₹438

M355		GENERAL	Tips Vol.Range5000μl Pack of 100 Gloine / Sp / Vp colour White	₹1,313
M356		GENERAL	Tips Vol.Range10μl Pack of 1000 Gloine / Sp / Vp colour White	₹875
M357		GENERAL	Tips Vol.Range200μl Pack of 1000 Gloine / Sp / Vp colour White	₹350

 YSP- SCIENTIFIC SOLUTION

MICROLIT - FILTER TIPS WITH BOX OF 96 TIPS



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M358	FM10	MICROLIT	Filtered Microtips 10 ul 96 no's with Tip Box	₹1,024
M359	FM20	MICROLIT	Filtered Microtips 20 ul 96 no's with Tip Box	₹1,024
M360	FM100	MICROLIT	Filtered Microtips 100 ul 96 no's with Tip Box	₹1,024
M361	FM200	MICROLIT	Filtered Microtips 200 ul 96 no's with Tip Box	₹1,024
M362	FM1000	MICROLIT	Filtered Microtips 1000 ul 96 no's with Tip Box	₹1,138

 YSP- SCIENTIFIC SOLUTION

MICROPIPETTE STAND - TOP TECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M363	MICROPIPETTE STAND 3	TOP TECH	Micropipette Stand 3 Holes	₹1,024
M364	MICROPIPETTE STAND 6	TOP TECH	Micropipette Stand 6 Holes	₹1,024

MICROPIPETTE STAND - TOP TECH



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M365	BS1	BLUESTAR	Cover Glass Size 18 x 18 Pkt of 10 gms Blue Star	₹166
M366	BS2	BLUESTAR	Cover Glass Size 22 x 22 Pkt of 10 gms Blue Star	₹166
M367	BS3	BLUESTAR	Cover Glass Size 22 x 30 Pkt of 10 gms Blue Star	₹170
M368	BS4	BLUESTAR	Cover Glass Size 22 x 40 Pkt of 10 gms Blue Star	₹170
M369	BS5	BLUESTAR	Cover Glass Size 22 x 50 Pkt of 10 gms Blue Star	₹170
M370	BS6	BLUESTAR	Cover Glass Size 22 x 60 Pkt of 10 gms Blue Star	₹180
M371	BS7	BLUESTAR	Cover Glass Size 18 mm Round	₹271
M372	BS14	BLUESTAR	Heamocytometer Cover Glass 20 X 25 mm (Pack of 20 no's)	₹473

COVER GLASS - ROHEM




SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M373		ROHEM	Cover Glass Size 18 x 18 Pkt of 10 gms Rohem	₹219
M374		ROHEM	Cover Glass Size 22 x 22 Pkt of 10 gms Rohem	₹219
M375		ROHEM	Cover Glass Size 22 x 40 Pkt of 10 gms Rohem	₹219
M376		ROHEM	Cover Glass Size 22 x 50 Pkt of 10 gms Rohem	₹219
M377		ROHEM	Cover Glass Size 22 x 60 Pkt of 10 gms Rohem	₹219
M378		ROHEM	Frosted Slide Rohem 1.35mm Pkt of 50 Slide	₹219

M379		ROHEM	ROHEM Rosenthal Counting Chamber bright line Indian	₹3,500
M380		ROHEM	ROHEM Sperm Counting Chamber bright line Indian	₹14,875

 YSP- SCIENTIFIC SOLUTION

MICROSLIDE - BLUESTAR



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M381	BS8	BLUESTAR	Blue Star PIC 1 (1.35 mm) Pkt of 50	₹131
M382	BS9	BLUESTAR	FROSTED SLIDE BLUE STAR 1.35MM PKT OF 50	₹175

 YSP- SCIENTIFIC SOLUTION

BLUE STAR - MICRO CONCAVITY SLIDE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
Size : 75mm x 25mm x 1.45mm thick. Edges are rounded to give better smear. Pk of 10 no's				
M383	BS10	BLUESTAR	One Concavity (per pkt of 10)	₹236
M384	BS11	BLUESTAR	Two Concavity (per pkt of 10)	₹245
M385	BS12	BLUESTAR	Three Concavity (per pkt of 10)	₹254

 YSP- SCIENTIFIC SOLUTION

BLUE STAR - V.D.R.L SLIDE (KLINE CONCAVITY SLIDES)



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

For Kline B.J.L., V.D.R.L. & other serodiagnostic tests. Size 75mm x 56mm x 3mm thick with spherical concavities, each 16mm in diameter x 1.75mm deep. Made from non-corrosive plate glass with bevelled & polished edges. Frosted surface is for notation.

M386	BS13	BLUESTAR	2, 4, 6, 8, 10 & 12 Concavities (each)	₹236
------	------	----------	----------------------------------------	------

 YSP- SCIENTIFIC SOLUTION

HEAMOCYTOMETER COVER GLASS



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

M387		BLUESTAR	Heamocytometer Cover Glass 20 X 25 mm (Pack of 20 no's)	₹473
------	--	----------	---------------------------------------------------------	------

 YSP- SCIENTIFIC SOLUTION

MICROTOME - RAYSCOPE



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
----	------	------------	--------------	------

M388	RYSARM10B	RAYSCOPE	RYS-ARM-10B Rayscope Microtome Regular Quality, Erma Japan Type Rotary Microtome Upgradable with Disposable Blades, Feed adjustment 1-25um, Minimum Thickness 1 um, Total Feed Excursions 28mm	₹44,950
------	-----------	----------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------

M389	RYSARM10A	RAYSCOPE	RYS-ARM-10A Rasyscope Microtome High Precision Quality , Spencer type Rotary Microtome Upgradable with Disposable Blades, Feed adjustment 1-50um, Minimum Thickness 1 um	₹68,875
------	-----------	----------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------

MULTIPARAMETER to check pH/Cond/TDS/DO etc. SEE PAGE " P"

MUFFLE FURNACE - YSP



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M390	YSPMF4	YSP	YSP-MF-4 Muffle Furnace Working size 4"x4"x9" Temperature up to 900c digital	₹42,050
M391	YSPMF4DIGI	YSP	YSP-MF-4 DIGI Muffle Furnace Working size 4"x4"x9" Temperature up to 1100c digital	₹56,550
M392	YSPMF5DIGI	YSP	YSP-MF-5 DIGI Muffle Furnace Working size 5"x5"x10" Temperature up to 900c digital	₹53,650
M393	YSPMF6DIGI	YSP	YSP-MF-6 DIGI Muffle Furnace Working size 6"x6"x12" Temperature up to 900c digital	₹71,775
M394	YSPMF6DIGI1100	YSP	YSP-MF-6 DIGI(1100) Muffle Furnace Working size 6"x6"x12" Temperature up to 1100c digital	₹94,250

ROTAMANTEL - REMI



SL	CODE	BRAND NAME	PRODUCT NAME	RATE
M395	1 RML	REMI	Heating Mantle With Magnetic Stirrer 1 RML 1 Litre Capacity	₹15,750
M396	2 RML	REMI	Heating Mantle With Magnetic Stirrer 2 RML 2 Litre Capacity	₹15,950
M397	3 RML	REMI	Heating Mantle With Magnetic Stirrer 3 RML 3 Litre Capacity	₹19,430



+91 96262 97366
+91 73387 21875



yspscientific@gmail.com
sales.yspscientific@gmail.com



B: 25 G1, Ground Floor, #7/10
Annai Avanthika Vengaivasal,
Chennai - 600100